



# Finnigan™ Kiel IV Carbonate Device

## Operating Manual

Revision A - 119 6820



**Finnigan™  
Kiel IV Carbonate  
Device**

**Operating Manual**

Revision A - 119 6820

© 2006 Thermo Electron Corporation. All rights reserved.

“Microsoft”, “Windows”, and “Excel” are registered trademarks of Microsoft Corporation. “Parafilm” is a registered trademark of American Can Company. “SNOOP” is a registered trademark of NUPRO Co. “Swagelok” is a registered trademark of the Crawford Fitting Company. “Teflon” is a registered trademark of E. I. du Pont de Nemours & Co. “Viton” is a registered trademark of Dupont Dow Elastomers. “Apiezon” denotes Apiezon Products, a business unit of M&I Materials Ltd.. “Pfeiffer” denotes Pfeiffer Vacuum, Inc.. “BOC Edwards” is a division of The BOC Group, Inc. “BOC Edwards” is a trade mark of The BOC Group.

All other trademarks are the property of Thermo Electron Corporation and its subsidiaries.

Technical information contained in this publication is for reference purposes only and is subject to change without notice. Every effort has been made to supply complete and accurate information; however, Thermo Electron Corporation assumes no responsibility and will not be liable for any errors, omissions, damage, or loss that might result from any use of this manual or the information contained therein (even if this information is properly followed and problems still arise).

This publication is not part of the Agreement of Sale between Thermo Electron Corporation and the purchaser of an Thermo Electron system. In the event of any conflict between the provisions of this document and those contained in Thermo Electron Corporation’s Terms and Conditions, the provisions of the Terms and Conditions shall govern.

System Configurations and specifications supersede all previous information and are subject to change without notice.

Printing History: Revision A printed in 8/16/06.

# Übereinstimmungserklärung gemäß EN 45014

Declaration of conformity according to EN 45014

Dichiarazione di conformità alla EN 45014

**Name des Herstellers:** Thermo Electron (Bremen) GmbH  
manufacturers name  
nome produttore

**Adresse des Herstellers:** Hanna-Kunath-Strasse 11  
manufacturers address  
indirizzo produttore  
28199 Bremen  
Germany

**erklärt, dass das Produkt**  
declares that the following product  
dichiara che il seguente prodotto

## **KIEL IV** Carbonate Device

**mit den folgenden Produktspezifikationen übereinstimmt:**

complies with the following product specifications

rispetta le seguenti specifiche del prodotto

**EMV (Störemissionen):** EN 50081-1; EN 55022 class B  
EMC (emissions)  
EMC (emissioni)

**EMV (Störfestigkeit):** EN 61000-3-2, -3; EN 61000-4-2, -3, -4, -5, -6, -11; EN 61000-6-2; EN 50204  
EMC (immunity)  
EMC (immunità)

**Elektrische Sicherheit:** EN 61010-1  
electrical safety  
sicurezza elettrica

**Ergänzende Informationen:**  
complementary information  
informazioni complementari

---

Dieses Produkt erfüllt die EMV-Richtlinie 89/336/EWG und Niederspannungsrichtlinie 73/23/EWG.

This product complies with EMC directive 89/336/EEC and Low Voltage Directive 73/23/EEC.

Questo prodotto rispetta la direttiva 89/336/EEC e la direttiva 73/23/EEC.

**Bremen, Germany, 23. März 2005**

**Technischer Leiter:**  
Operations Manager  
Direttore fabbricazione

**Thermo**  
ELECTRON CORPORATION

*ppa. G. Schröder*

**Notice on Lifting and Handling of  
Thermo Electron Bremen Instruments**

For your safety, and in compliance with international regulations, the physical handling of this Thermo Electron Bremen instrument *requires a team effort* for lifting and/or moving the instrument. This instrument is too heavy and/or bulky for one person alone to handle safely.

**Notice on the Proper Use of  
Thermo Electron Bremen Instruments**

In compliance with international regulations: If this instrument is used in a manner not specified by Thermo Electron Bremen, the protection provided by the instrument could be impaired.

## WEEE Compliance

This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive 2002/96/EC. It is marked with the following symbol:



Thermo Electron has contracted with one or more recycling/disposal companies in each EU Member State, and this product should be disposed of or recycled through them. Further information on Thermo Electron's compliance with these Directives, the recyclers in your country, and information on Thermo Electron products which may assist the detection of substances subject to the RoHS Directive are available at [www.thermo.com/WEEERoHS](http://www.thermo.com/WEEERoHS).

## WEEE Konformität

Dieses Produkt muss die EU Waste Electrical & Electronic Equipment (WEEE) Richtlinie 2002/96/EC erfüllen. Das Produkt ist durch folgendes Symbol gekennzeichnet:



Thermo Electron hat Vereinbarungen getroffen mit Verwertungs-/Entsorgungsanlagen in allen EU-Mitgliederstaaten und dieses Produkt muss durch diese Firmen wiederverwertet oder entsorgt werden. Mehr Informationen über die Einhaltung dieser Anweisungen durch Thermo Electron, die Verwerter und Hinweise die Ihnen nützlich sein können, die Thermo Electron Produkte zu identifizieren, die unter diese RoHS Anweisung fallen, finden Sie unter [www.thermo.com/WEEERoHS](http://www.thermo.com/WEEERoHS).

## Conformité DEEE

Ce produit doit être conforme à la directive européenne (2002/96/EC) des Déchets d'Equipements Electriques et Electroniques (DEEE). Il est marqué par le symbole suivant:



Thermo Electron s'est associé avec une ou plusieurs compagnies de recyclage dans chaque état membre de l'union européenne et ce produit devrait être collecté ou recyclé par celles-ci. Davantage d'informations sur la conformité de Thermo Electron à ces directives, les recycleurs dans votre pays et les informations sur les produits Thermo Electron qui peuvent aider la détection des substances sujettes à la directive RoHS sont disponibles sur [www.thermo.com/WEEERoHS](http://www.thermo.com/WEEERoHS).

# Read This First

Welcome to the Thermo Electron, Finnigan Kiel IV Carbonate Device system!

This *Finnigan Kiel IV Carbonate Device Operating Manual* describes how to setup and use your Kiel IV Carbonate Device. In addition, this manual describes principle hardware components.

It includes the following chapters:

- **Chapter 1: “Preinstallation Requirements”** summarizes requirements related to site, power, the reference sample and the various gases in use before operating your Finnigan Kiel IV Carbonate Device.
- **Chapter 2: “Hardware Components”** treats instrument layout, connections, vacuum system, valves and traps, oven and oven control, autosampler, Liquid Nitrogen Refill device, Reference Gas Refill, acid flow and pinch valve.
- **Chapter 3: “Isodat 2.5”** describes how to start Isodat 2.5 and subsequently how to create a Kiel IV Carbonate Device-related configuration.

This chapter further denotes how to create a new Kiel IV Carbonate Device method and a new Kiel IV Carbonate Device sequence in Isodat 2.5’s Acquisition Mode.

Finally, it outlines interpretation of result files, time slicing and interfering masses.

- **Chapter 4: “Basic Operations”** describes several test routines as leak check, bakeout, autosampler operation, capillary matching, cleaning the acid valve, adjusting the liquid nitrogen refill sensor, pinch valve operation, troubleshooting, elementary handling of the Finnigan Kiel IV Carbonate Device, vial test, phosphoric acid preparation and sample vial handling.

- **Chapter 5: “Measurement Procedures for Real Samples”** deals with the measurement principle, sample placement into a vial, preparation of carbonates and IRMS and the measurement procedure itself. Furthermore, it outlines how the quality of result data is checked.

The chapter comprises information about referencing vs. VPDB and about Reference Refill as well.

- **Chapter 6: “Technical Information”** outlines spare parts and consumables, the valve unit and valve replacement. It contains information about IAEA primary standards.

Furthermore, advice is given for internal leak checking, maintenance and programming. Finally, this chapter contains schematics of compressed air supply, vacuum system and circuit diagrams.

## Changes to the Manual

To suggest changes to this manual, please send your comments to:

Editors, Technical Documentation  
Thermo Electron (Bremen) GmbH  
Finnigan Advanced Mass Spectrometry  
Hanna-Kunath-Str. 11

D-28199 Bremen

Germany

documentation@thermo-bremen.com

You are encouraged to report errors or omissions in the text or index. Thank you.

## Typographical Conventions

Typographical conventions have been established for Thermo Electron manuals for the following:

- Data input
- Admonitions
- Topic headings

## Data Input

Throughout this manual, the following conventions indicate data input and output via the computer:

- Messages displayed on the screen are represented by capitalizing the initial letter of each word and by italicizing each word.
- Input that you enter by keyboard is identified by quotation marks: single quotes for single characters, double quotes for strings.
- For brevity, expressions such as “choose File > **Directories**” are used rather than “pull down the File menu and choose Directories.”
- Any command enclosed in angle brackets < > represents a single keystroke. For example, “press <**F1**>” means press the key labeled *F1*.
- Any command that requires pressing two or more keys simultaneously is shown with a plus sign connecting the keys. For example, “press <**Shift**> + <**F1**>” means press and hold the <Shift> key and then press the <F1> key.
- Any button that you click on the screen is represented in bold face letters. For example, “click on **Close**”.

## Admonitions

Admonitions contain information that is important, but not part of the main flow of text.

Admonitions can be of the following types:

- **Note** – information that can affect the quality of your data. In addition, notes often contain information that you might need if you are having trouble.
- **Caution** – information necessary to protect your instrument from damage.
- **Warning** – hazards to human beings. Each Warning is accompanied by a Warning symbol.

## Topic Headings

The following headings are used to show the organization of topics within a chapter:

# Chapter Name

The following headings appear in the left column of each page:

## Second Level Topics

### Third Level Topics

#### Fourth Level Topics

## Safety and EMC Information

In accordance with our commitment to customer service and safety, these instruments have satisfied the requirements for the European CE Mark including the Low Voltage Directive.

Designed, processor and tested in an ISO9001 registered facility, this instrument has been shipped to you from our manufacturing facility in a safe condition.

**Caution** This instrument must be used as described in this manual. Any use of this instrument in a manner other than described here may result in instrument damage and/or operator injury. ▲

## Identifying Safety Information

The *Finnigan Kiel IV Carbonate Device Operating Manual* contains precautionary statements that can prevent personal injury, instrument damage, and loss of data if properly followed. Warning symbols which alert the user to check for hazardous conditions. These appear throughout the manual, where applicable, and are defined in [Table i](#) on [page i-v](#).

Table i. Warning Symbols

Symbol	Description
	<b>General</b> This general symbol indicates that a hazard is present, which if not avoided, could result in injuries. The source of danger is described in the accompanying text. ▲
	<b>Cold Burns Hazard</b> Wear protective clothing. ▲
	<b>Electric Shock</b> High voltages capable of causing personal injury are used in the instrument. The instrument must be shut down and disconnected from line power before service or repair work is performed. ▲
	<b>Magnetic Field</b> Keep away from heart pacemakers, computers, credit cards, and any other magnetically sensitive device. ▲
	<b>Noxious</b> This symbol alerts to hazards resulting from noxious fumes. ▲
	<b>Hot Surface / Heat</b> Allow heated components to cool down before servicing them! ▲
	<b>Poisonous Gases</b> This symbols points to possible danger because of poisonous gases and vapors. ▲

### Instrument-Specific Hazards

Every instrument has specific hazards, so be sure to read and comply with the following precautions. They will help ensure the safe, long-term use of your system.

1. Before plugging in any of the instrument modules or turning on the power, always make sure that the voltage and fuses are set appropriately for your local line voltage.
2. Only use fuses of the type and current rating specified. Do not use repaired fuses and do not short-circuit the fuse holder.
3. The supplied power cord must be inserted into a power outlet with a protective earth contact (ground). When using an extension cord, make sure that the cord also has an earth contact.

4. Do not change the external or internal grounding connections. Tampering with or disconnecting these connections could endanger you and/or damage the system.

**Caution** The instrument is properly grounded in accordance with regulations when shipped. You don't need to make any changes to the electrical connections or to the instrument's chassis to ensure safe operation. ▲

5. Never run the system without the housing on. Permanent damage can occur.
6. Do not turn the instrument on if you suspect that it has incurred any kind of electrical damage. Instead, disconnect the power cord and contact a Service Representative for a product evaluation. Do not attempt to use the instrument until it has been evaluated. (Electrical damage may have occurred if the system shows visible signs of damage, or has been transported under severe stress.)
7. Damage can also result if the instrument is stored for prolonged periods under unfavorable conditions (e.g. subjected to heat, water, etc.).
8. Always disconnect the power cord before attempting any type of maintenance.
9. Capacitors inside the instrument may still be charged even if the instrument is turned off. The superconducting magnet is still charged even if the instrument is turned off.
10. Never try to repair or replace any component of the system that is not described in this manual without the assistance of your service representative.

# Contents

<b>Chapter 1</b>	<b>Preinstallation Requirements.....</b>	<b>1-1</b>
	Site Requirements .....	1-2
	Power Requirements .....	1-2
	Gas Requirements .....	1-2
	Further Requirements .....	1-3
	Reference Sample .....	1-3
<b>Chapter 2</b>	<b>Hardware Components .....</b>	<b>2-1</b>
	Layout.....	2-2
	Upper Section .....	2-2
	Lower Section .....	2-3
	Connections.....	2-12
	Connecting Kiel IV Carbonate Device to IRMS .....	2-12
	Vacuum System .....	2-16
	Valves and Traps.....	2-17
	Pneumatic Valves .....	2-18
	Valves of Dual Inlet System .....	2-18
	Trap 1 .....	2-21
	Microvolume (Trap 2) .....	2-23
	Autocool Unit.....	2-24
	Oven and Oven Control .....	2-27
	Programming - Step 1 .....	2-27
	Programming - Step 2.....	2-28
	Alternative - Automatic Programming.....	2-28
	Tune Scan after Programming .....	2-28
	Autosampler.....	2-29
	Installing and Removing Magazine .....	2-29
	Adjusting Magazine Position.....	2-30
	Adjusting Piston Speed .....	2-31
	Adjusting Piston Height.....	2-32
	Proximity Switch.....	2-32
	Turret Motor .....	2-33
	Liquid Nitrogen Refill Device.....	2-34
	Warning Note for Liquid Nitrogen Supply .....	2-36
	LN2 Transfer from Refill Device into Dewar .....	2-38
	Reference Gas Refill .....	2-39
	Filling Reference Refill Device from External Source ...	2-41
	Acid Flow.....	2-43
	Pinch Valve.....	2-44

<b>Chapter 3</b>	<b>Isodat 2.5</b> .....	<b>3-1</b>
	Creating a Kiel IV Carbonate Device Configuration .....	3-2
	Acquisition Mode .....	3-4
	Starting Acquisition Mode .....	3-4
	Activating Toolbars .....	3-4
	Accessories Bar and its Components.....	3-5
	Changing Visibility of its Components .....	3-8
	Kiel IV Carbonate Window .....	3-8
	Dual Inlet Window.....	3-8
	File Browser .....	3-13
	Info Window .....	3-16
	Creating a New Method.....	3-20
	Structure of Methods for Kiel IV Carbonate Device.....	3-21
	Instrument Tab.....	3-21
	Peripherals Tab .....	3-23
	Evaluation Tab.....	3-27
	Printout Tab .....	3-29
	Saving a Method .....	3-29
	Creating a New Sequence.....	3-31
	Saving a Sequence .....	3-33
	Starting a Sequence .....	3-35
	Standards .....	3-38
	Using Generic Editor .....	3-38
	Using Standard Editor .....	3-40
	Excel Export.....	3-42
	Dyn Externals .....	3-44
	Service Scripts .....	3-49
	Take Magazine.....	3-49
	Load Magazine.....	3-51
	Acid Drop Test .....	3-52
	Standby & Drop .....	3-53
	Standby and Pump.....	3-55
	Terminate .....	3-55
	Interpreting Results.....	3-56
	Raw Tab .....	3-56
	Individual Columns of Result File.....	3-57
	Evaluated Tab .....	3-57
	Grid Errors Tab .....	3-58
	Grid Infos Tab.....	3-58
	Sequence Line Tab.....	3-60
	Time Slicing.....	3-60
	Raw Sample Tab .....	3-60
	Raw Reference Tab .....	3-60
	Raw Complete Tab .....	3-61
	Raw Ratios Sample Tab .....	3-62
	Raw Ratios Reference Tab.....	3-62
	Raw Ratios Complete Tab .....	3-63
	Ratios Tab.....	3-63

Interfering Masses .....	3-64
Code Example.....	3-65
<b>Chapter 4 Basic Operations.....</b>	<b>4-1</b>
Leak Check .....	4-2
Water .....	4-3
Air.....	4-3
CO <sub>2</sub> .....	4-3
Dual Inlet Ar Signal .....	4-4
Advanced Leak Checking Procedure .....	4-4
Bakeout of Kiel IV Carbonate Device .....	4-6
Operating the Autosampler .....	4-7
Connecting Vials.....	4-7
Disconnecting Vials .....	4-8
Matching Sample Capillary to Standard Capillary.....	4-9
Cleaning Acid Valve.....	4-14
Disassembling Acid Valve .....	4-14
Brief and Superficial Cleaning Procedure .....	4-18
Thorough Cleaning Procedure .....	4-19
Reassembling Acid Valve.....	4-19
Adjusting Liquid Nitrogen Refill Sensor .....	4-20
Operating Pinch Valve.....	4-21
Troubleshooting.....	4-21
Elementary Handling of Kiel IV Carbonate Device.....	4-22
Vial Test .....	4-23
Principle.....	4-23
Performing a Vial Test .....	4-23
Possible Error Messages during Vial Test .....	4-25
Phosphoric Acid Preparation.....	4-26
Removing Water from Phosphoric Acid.....	4-27
Adding Phosphorous Pentoxide .....	4-27
Handling Sample Vials.....	4-28
Manual Cleaning of Sample Vials .....	4-28
Automatic Cleaning of Sample Vials .....	4-29
<b>Chapter 5 Measurement Procedures for Real Samples.....</b>	<b>5-1</b>
Introduction.....	5-2
Measurement Principle .....	5-3
Placing Sample into Vial .....	5-4
First Alternative.....	5-4
Second Alternative .....	5-5
Preparing Carbonate and IRMS .....	5-6
Procedure.....	5-7
Checking Quality of Result Data .....	5-11

Referencing vs VPDB..... 5-14  
 Remark on the Strange Mathematics of Delta Values... 5-16  
 Ion Correction ..... 5-17  
 Neogloboquadrina Pachyderma (Ehrenberg, 1894)..... 5-19  
 Events during Sample Measurement ..... 5-19  
 Reference Refill ..... 5-22  
 Selecting Correct Amount of Reference Gas..... 5-24

**Chapter 6 Technical Information ..... 6-1**

Spare Parts and Consumables..... 6-3  
 Installation Kit ..... 6-7  
 Acid Valve..... 6-8  
 Trapping Volume ..... 6-9  
 Cooling Unit..... 6-11  
 Microvolume ..... 6-12  
 Valve Unit ..... 6-13  
 Valve Replacement..... 6-14  
 Definitions..... 6-17  
 IAEA Primary Standards ..... 6-17  
 Checking for Internal Leaks ..... 6-18  
 High Vacuum Side..... 6-19  
 Carbonate Side Leaks ..... 6-22  
 Maintenance ..... 6-26  
 Turbo Pumps and Fore Pumps ..... 6-26  
 O Ring Seals ..... 6-26  
 Autocool Unit..... 6-26  
 Programming Information ..... 6-26  
 Board Base Addresses ..... 6-26  
 Device Addresses ..... 6-27  
 Registry Use ..... 6-32  
 Compressed Air Supply..... 6-33  
 Vacuum Schematic ..... 6-35  
 Circuit Diagrams ..... 6-36

**Glossary ..... G-1**

**Index ..... I-1**

# Figures

Kiel IV Carbonate Device - Site Requirements .....	1-2
Finnigan Kiel IV Carbonate Device - Front View .....	2-2
Heating Cabinet - Front View .....	2-3
Front Panel - Upper Part .....	2-4
Close-Up View of Trap 1 and Sensor for Fill Level .....	2-4
Front Panel - Lower Part .....	2-5
Control Panel in Close-Up View .....	2-5
Parts of Front Panel and Rear Panel .....	2-6
Fore Pumps at Rear Panel .....	2-7
Electronics Cabinet .....	2-7
Safety Switch-Off for Acid Valve .....	2-8
Fiberline to Serial Bus Interface .....	2-9
Interface Proximity Switches - Inlet Control Board .....	2-9
Parts within Right Panel .....	2-10
Position of VM2 Fore Vacuum Gauge .....	2-11
Valve 7 and Valve 8 .....	2-11
Connecting Viton Tubes of Acid Valves .....	2-14
Fixing Alligator Crimps to Capillaries .....	2-15
Capillary Heating Transformer .....	2-15
Turbo Pump with Turbo Pump Controller .....	2-17
Trap Arrangement .....	2-17
Parts of a Pneumatic Valve .....	2-18
Pliers to Insert a Pneumatic Valve .....	2-18
Dual Inlet System Valve .....	2-19
Parts of Dual Inlet System Valve .....	2-19
Double Valve Block .....	2-20
Manifold Block with Four Solenoid Valves .....	2-20
Compressed Air Distributor .....	2-21
Schematic of Trap 1 .....	2-22
Microvolume and its Parts .....	2-23
Microvolume Parts to be Inserted into Autocool Unit .....	2-24
Autocool Unit and Interior of Funnel .....	2-24
Cascade with Three Holes .....	2-25
Mounting Instructions for Autocool Unit onto Trap 1 .....	2-25
Dewar and “Lab Boy” .....	2-26
Oven Control .....	2-27
Turret and Correct Vertical Placement of Vials .....	2-29
Turret and Lid for Line 1 and Line 2 Position .....	2-29
Magazine with Cover Plate and Vials .....	2-30
Piston Speed Adjustment .....	2-31
Connections of Pneumatic Levers for Pistons .....	2-31
Adjusting Piston Height .....	2-32
Proximity Switch .....	2-33

Turret Motor and Position Sensing Array .....	2-33
Tool to Readjust Array of Infrared Light Barriers .....	2-34
Schematic of a Standard Liquid Nitrogen Refill Device ...	2-34
Magnetic Valve and Liquid Nitrogen Safety Unit .....	2-35
Example of a Liquid Nitrogen Refill Device .....	2-35
Valves and Distributor for Liquid Nitrogen (Top View) .	2-37
Valves and Distributor for Liquid Nitrogen (Side View) .	2-37
Valves and Distributor for Liquid Nitrogen .....	2-38
Position of Refill Sensor .....	2-38
Liquid Nitrogen Refill Sensor .....	2-39
Reference Gas Refill Units .....	2-40
Manual Valves .....	2-41
Filling Reference Refill Device from External Source - I ..	2-41
Filling Reference Refill Device from External Source - II .	2-42
Schematic of Acid Flow .....	2-43
Acid Reservoir and Manual Valve .....	2-43
Pinch Valve and Drop Counters .....	2-44
Pinch Valve - Mounted .....	2-44
Pinch Valve - Disassembled .....	2-45
Feeding Acid Tubing into Pinch Valve .....	2-45
Mounting Pinch Valve upon Acid Valve .....	2-45
Opening Configurator .....	3-2
Adding a New Configuration .....	3-2
Renaming the New Configuration .....	3-2
Expanding Tree of IRMS .....	3-2
Appending Kiel IV Carbonate Set to IRMS .....	3-3
Kiel IV Carbonate Device Appended to Capillary Port .....	3-3
Visibility of Individual Toolbars .....	3-5
Status Bar .....	3-6
Components of Accessories Bar .....	3-7
Marking or Unmarking Accessories .....	3-8
Dual Inlet Window .....	3-9
Switching Positions of Changeover Valve .....	3-9
Dual Inlet Measurement Loop .....	3-10
Setting Idle Time .....	3-11
Setting Integration Time .....	3-11
Shot Noise Limits of Precision .....	3-12
File Browser .....	3-13
Modifying Result Path .....	3-15
Location of ISL Scripts for Kiel IV Carbonate Device .....	3-15
Info Window .....	3-17
Commands for Using Info Window .....	3-17
Info Properties Window .....	3-18
Enabling Log File .....	3-18
Offline View .....	3-19
Creating a New Method .....	3-21
Instrument Tab - Experiment .....	3-22
Instrument Tab - Isotope MS .....	3-22
Instrument Tab - Peak Center .....	3-23

Instrument Tab - Reference Refill .....	3-23
Peripherals Tab - Dual Inlet System .....	3-24
Peripherals Tab - Background .....	3-24
Peripherals Tab - Pressure Adjust .....	3-25
Peripherals Tab - Time Slicing .....	3-26
Peripherals Tab - Carbonate Device .....	3-26
Evaluation Tab - Cycle .....	3-27
Evaluation Tab - Extended Parameters .....	3-27
Evaluation Tab - Evaluation Type .....	3-28
Evaluation Tab - Standard Parameter (I) .....	3-28
Evaluation Tab - Standard Parameter (II) .....	3-28
Printout Tab .....	3-29
Saving a Method .....	3-30
Creating a New Sequence .....	3-31
Selecting Number of Samples .....	3-32
Sequence Grid .....	3-32
Saving a Sequence .....	3-34
Defining Parameters for Handling Results .....	3-35
Defining Full Path for Results Storage .....	3-36
Defining Parameters for Results Export .....	3-36
Defining Parameters for Results Printout .....	3-37
Properties Box - Comment .....	3-37
Selecting ISL Scripts to be Executed .....	3-37
Standby .....	3-38
Selecting PrimaryStandards.std .....	3-39
Opening Primary Standards Database .....	3-39
Opening Standard Editor .....	3-40
Creating New Standard .....	3-41
New Standard Appearing in Standard Editor .....	3-41
Excel Export Template .....	3-43
Calling Dyn Externals vs. Service Scripts .....	3-44
Basic Functionality Tab .....	3-44
Pressure Fast Adjust Tab .....	3-45
Temp Settings Tab .....	3-46
Process Timing Tab .....	3-47
About Acid Tab .....	3-48
Calling Dyn Externals vs. Service Scripts .....	3-49
Available Service Scripts .....	3-49
Service Script Take Magazine .....	3-49
Service Script Load Magazine .....	3-51
Service Script Acid Drop Test .....	3-52
Service Script Standby and Drop .....	3-53
Service Script Standby and Pump .....	3-55
Kiel IV Carbonate Device in Standby Mode .....	3-55
Actions during Terminate .....	3-56
Raw Tab - Part I .....	3-56
Raw Tab - Part II .....	3-56
Evaluated Tab .....	3-57
Grid Errors Tab .....	3-58

Grid Infos Tab .....	3-59
Sequence Line Tab .....	3-60
Raw Sample Tab .....	3-60
Raw Reference Tab .....	3-61
Raw Complete Tab .....	3-61
Raw Ratios Sample Tab .....	3-62
Raw Ratios Reference Tab .....	3-62
Raw Ratios Complete Tab .....	3-63
Ratios Tab .....	3-63
Interfering Masses Tab .....	3-64
Mass Spectrum of Background Gas Composition .....	4-2
Mass Spectrum of Background Gas Composition .....	4-2
Scan Parameters during Leak Check .....	4-4
Instrument Control as Aid when Leak Checking .....	4-5
Connecting Vials .....	4-7
Disconnecting Vials .....	4-8
Matching Sample Capillary to Standard Capillary .....	4-9
Connecting Transfer Tube to IRMS .....	4-10
Connection for Installation of Shortcut Link Tube .....	4-10
Connecting Transfer Tube to Kiel IV Carbonate Device .....	4-11
Crimping Device at the End of a Capillary .....	4-11
Vacuum Scheme of Kiel IV Carbonate Device .....	4-13
Parts of Acid Valve .....	4-14
Acid Valve in Position .....	4-15
Acid Valve - I .....	4-15
Acid Valve - II .....	4-16
Drop Counter and Electrical Feedthrough .....	4-16
Mounted Spring Plate - Bottom View .....	4-17
Mounted Valve - Bottom View .....	4-17
Valve Flange with Frit .....	4-18
Cleaning Acid Valve .....	4-18
Kiel IV Carbonate Window .....	4-22
ISL Script is Loaded into Method .....	4-24
Selecting ISL Script for Vial Test .....	4-24
Saving Method as a New Vial Test Method .....	4-25
Removing Water from Phosphoric Acid .....	4-27
Checking Specific Gravity of Phosphoric Acid .....	4-28
Oven Rack Containing Some Vials .....	4-29
Finnigan Kiel IV Carbonate Device - Front View .....	5-2
Common Glass Tube .....	5-3
Sample Vial - Part No. 075 4960 .....	5-5
Correct Sample Placement .....	5-5
Flowchart of Carbonate Process Including Timing .....	5-7
Leak Test .....	5-8
Trap Temperatures during Sample Preparation .....	5-9
Expansion Flowchart .....	5-10
Pressure vs. Mass .....	5-11
Signal vs. Pressure .....	5-12
$\delta$ vs. Signal .....	5-13

Calculation Example ..... 5-15

Neogloboquadrina Pachyderma ..... 5-19

Evolution of Measurement Process in Info Window ..... 5-20

Reference Refill Process ..... 5-22

Reference Refill Process - Step 1 ..... 5-23

Reference Refill Process - Step 2 ..... 5-23

Reference Refill Process - Step 3 ..... 5-24

Vacuum Scheme of Device and Dual Inlet ..... 6-2

Capillary for Acid Valve (Part No. 059 8651) ..... 6-3

Oil for Fore Vacuum Pump (Part No. 109 4301) ..... 6-4

Flexible Acid Resistant Tubing (Part No. 056 7830) ..... 6-4

Teflon Gasket for Acid Valve (Part No. 059 8671) ..... 6-4

Gold Gasket (Part No. 054 5270) ..... 6-4

Heatable Capillary (Part No. 067 1182) ..... 6-5

Pinch Valve (Part No. 119 1170) ..... 6-5

Complete Membrane (Part No. 065 3010) ..... 6-5

Gold Stamp (Part No. 065 3041) ..... 6-6

Jacket Ring (Part No. 055 3140) ..... 6-6

Sample Vial (Part No. 075 4960) ..... 6-6

Proximity Switch (Part No. 106 9490) ..... 6-6

Gold Gasket for Microvolume (Part No. 055 1010) ..... 6-7

Copper Shim (Part No. 100 7730) ..... 6-7

Tool to Adjust Turret Readout (Part No. 115 7390) ..... 6-7

Parts of Acid Valve ..... 6-8

Acid Valve (Part No. 106 9450) ..... 6-8

Trapping Volume (Part No. 100 7740) ..... 6-10

Important Parts of Trapping Volume ..... 6-10

Cooling Unit (Part No. 079 2400) ..... 6-11

Microvolume (Part No. 116 3150) ..... 6-13

Microvolume and its Parts ..... 6-13

Valve Unit (Part No. 065 3001) ..... 6-14

Tools for Valve Replacement ..... 6-15

IAEA Primary Standards ..... 6-17

Checking VM1 Pressure Increase of Line 1 ..... 6-18

Leak Check of Pneumatic Valves - Step 1 ..... 6-20

Leak Check of Pneumatic Valves - Step 2 ..... 6-20

Leak Check of Pneumatic Valves - Step 5 ..... 6-21

Leak Check of Pneumatic Valves - Step 6 ..... 6-22

Checking V2 - Step 1 ..... 6-23

Checking V2 - Step 2 ..... 6-23

Checking V1 - Step 1 ..... 6-24

Checking V1 - Step 2 ..... 6-24

Checking V1 - Step 3 ..... 6-25

Adjusting Base Address of Power Distribution Board ..... 6-27

Connections on Inlet Control Board ..... 6-31

Adjusting Sensitivity of Drop Counters ..... 6-32

Compressed Air Supply ..... 6-33

Compressed Air Supply - Schematic ..... 6-34

Kiel IV Carbonate Device - Vacuum Schematic ..... 6-35

## Figures

Circuit Diagram of Device - Control Logic .....	6-36
Circuit Diagram of Device - Mains Connections .....	6-37

# Tables

Kiel IV Carbonate Device - Gas Requirements .....	1-3
Important Parts within Right Panel .....	2-10
Parts of Dual Inlet System Valve .....	2-19
Parameters for Step 1 of Programming .....	2-27
Parameters for Step 2 of Programming .....	2-28
Indications of Dual Inlet Window .....	3-9
Switching Positions of Changeover Valve .....	3-10
Instrument Tab - Experiment .....	3-22
Instrument Tab - Isotope MS .....	3-22
Instrument Tab - Peak Center .....	3-23
Instrument Tab - Reference Refill .....	3-23
Peripherals Tab - Dual Inlet System .....	3-24
Peripherals Tab - Background .....	3-24
Peripherals Tab - Pressure Adjust .....	3-25
Peripherals Tab - Time Slicing .....	3-26
Peripherals Tab - Carbonate Device .....	3-26
Evaluation Tab - Cycle .....	3-27
Evaluation Tab - Extended Parameters .....	3-28
Evaluation Tab - Evaluation Type .....	3-28
Evaluation Tab - Standard Parameter .....	3-29
Printout Tab .....	3-29
Sequence Grid .....	3-32
Defining Full Path for Results Storage .....	3-36
Defining Parameters for Results Export .....	3-36
Defining Parameters for Results Printout .....	3-37
Properties Box - Comment .....	3-37
Selecting ISL Scripts to be Executed .....	3-38
Standby .....	3-38
Parameters of Basic Functionality Tab .....	3-45
Parameters of Pressure Fast Adjust Tab .....	3-45
Parameters of Temp Settings Tab .....	3-46
Parameters of Process Timing Tab .....	3-47
Parameters of About Acid Tab .....	3-48
Explanations of Columns of Result File .....	3-57
Parameters of Grid Infos Tab .....	3-59
Basic Operations via Kiel IV Carbonate Window .....	4-23
Some Spare Parts and Consumables (Part No. 119 1160) .....	6-3
Parts of Installation Kit (Part No. 115 7800) .....	6-7
Parts of Acid Valve (Part No. 106 9450) .....	6-9
Parts of Trapping Volume (Part No. 100 7740) .....	6-10
Important Parts of Cooling Unit (Part No. 079 2400) .....	6-11
Important Parts of Microvolume (Part No. 116 3150) .....	6-12
Important Parts of Valve Unit (Part No. 065 3001) .....	6-14
Parts Needed for Valve Replacement .....	6-14
Definitions .....	6-17

**Tables**

Board Base Addresses .....	6-27
Device Addresses - Complete .....	6-28
DIO Parameters .....	6-29
ADC Parameters .....	6-30
DAC Parameters .....	6-30

# Chapter 1 Preinstallation Requirements


This chapter contains the following topics:

- “Site Requirements” on page 1-2
- “Power Requirements” on page 1-2
- “Gas Requirements” on page 1-2
- “Further Requirements” on page 1-3
- “Reference Sample” on page 1-3

## Site Requirements

The Finnigan Kiel IV Carbonate Device is attached to Thermo Electron isotope ratio mass spectrometers, e.g. a Finnigan DELTA V Plus or a Finnigan DELTA V Advantage, equipped with a Dual Inlet system.

It is placed stand-alone but must be arranged next to the IRMS within a maximum distance of 30 cm. The distance to walls may not be less than 60 cm. The space required is 900 mm width x 900 mm depth. Its height is 1900 mm, and it weighs approximately 100 kg (220 lb). See Figure 1-1.



**Figure 1-1.** Kiel IV Carbonate Device - Site Requirements

## Power Requirements

The Finnigan Kiel IV Carbonate Device will be supplied by the IRMS line distributor. The power consumption of the IRMS will increase by 1.2 kW.

**Note** It is absolutely necessary to run your instrument without disruptions of power supply! Thus, if your local area is susceptible to corrupted power or power disruptions, an uninterruptible power supply (UPS) should be installed in your laboratory. ▲

## Gas Requirements



**Warning** All gas lines should be oil-free and preferably flame-dried. The gas lines or gas tanks should be at a distance of 1 m to 1.5 m to the instrument. ▲



**Warning** All regulators should be oil- and grease-free and be specified for gases of high purity. The supply lines should terminate with 1/8" male Swagelok® type connectors. Thermo Electron (Bremen) recommends to use regulators with an outlet pressure range between 0 and 5 bar (that is, 0 and 73 psi). ▲



**Warning** All compressed air tubing and the air compressor should contain a water and oil trap. Water and oil may fill the compressed air supply and destroy the pneumatic valves! ▲

For operation, the Kiel IV Carbonate Device needs the gases summarized in [Table 1-1](#).

**Table 1-1.** Kiel IV Carbonate Device - Gas Requirements

Gas	Comment
He, N <sub>2</sub>	1 bar as sample vial vent gas
liquid N <sub>2</sub>	For the liquid nitrogen-cooled trap provide approximately 0.5 l of liquid nitrogen per sample.
CO <sub>2</sub>	Used as reference gas. See <a href="#">“Reference Refill”</a> on <a href="#">page 5-22</a> .
Ar	Sometimes, it may be necessary to check the unit for leaks. Therefore, use an argon tank.
compressed air	Supplied by the compressed air distributor of the IRMS. Should be in the range between 2.8 bar and 6 bar (that is, between 40 and 87 psi).

## Further Requirements

- a laboratory to wash and prepare sample vials
- a fume hood, hot plate and stirrer to prepare phosphoric acid<sup>1</sup>
- contaminant-free weighing instruments<sup>1</sup>
- phosphorous pentoxide (500 g)<sup>1</sup>
- acetone p.a. for cleaning vials and weighing instruments<sup>1</sup>
- a suitable micro balance to verify the installation weight specifications<sup>1</sup>
- a liquid nitrogen tank (approximately 30 l per day)<sup>1</sup>

## Reference Sample

To demonstrate the specified precision for small samples, it is required that the customer buys and supplies a sample of NBS-19 from the International Atomic Energy Agency or the National Bureau of Standards. Refer to [“IAEA Primary Standards”](#) on [page 6-17](#).

<sup>1</sup>Keep available within installation distance.




## Chapter 2 Hardware Components

This chapter provides information about the various hardware components of the Finnigan Kiel IV Carbonate Device. It contains the following topics:

- “Layout” on page 2-2
- “Connections” on page 2-12
- “Vacuum System” on page 2-16
- “Valves and Traps” on page 2-17
- “Oven and Oven Control” on page 2-27
- “Autosampler” on page 2-29
- “Liquid Nitrogen Refill Device” on page 2-34
- “Reference Gas Refill” on page 2-39
- “Acid Flow” on page 2-43
- “Pinch Valve” on page 2-44

**Note** For information about spare parts and part numbers, refer to [Chapter 6: “Technical Information”](#). ▲


**Layout** This section describes the parts of the Finnigan Kiel IV Carbonate Device. Figure 2-1 shows the device in front view.



**Figure 2-1.** Finnigan Kiel IV Carbonate Device - Front View

**Upper Section** Inside the heating cabinet, a round plate with holes can be found that holds sample vials. This is called the turret. Two pistons move the vials up or down in order to connect them to the vacuum system. For historical reasons, this connection port is called the acid valve. Phosphoric acid is stored in a borosilicate glass container and fed to the

acid valve via Viton tubing. A pinch valve controls the acid flow. Unused vials can be stored in two vial racks located in the upper part of the thermostated oven. See Figure 2-2.




**Figure 2-2.** Heating Cabinet - Front View

### Lower Section


After the grey side panels have been removed the lower cabinet becomes visible. See Figure 2-3 and Figure 2-5.

In the upper part of Figure 2-3, valve systems can be seen on the right side. They belong to trap 1 shown as **a**. In the middle of this picture the pneumatic levers for the pistons **c** are located, each with a position sensor **b**. Towards the back the turret motor **e** is barely visible. The sensor for the fill level of liquid nitrogen **d** in the dewar is located rightmost in this part of the cabinet.



**Figure 2-3.** Front Panel - Upper Part

Figure 2-4 depicts trap 1 (a) and the sensor for the fill level of liquid nitrogen (d) in close-up view.



**Figure 2-4.** Close-Up View of Trap 1 and Sensor for Fill Level

The lower part of the cabinet contains the fore vacuum pumps and pneumatic control valves that switch on and off compressed air that in turn operates valves and the pneumatic levers for the pistons. A dewar located upon a lab boy can be found here as well. See Figure 2-5. The dewar is used to store a limited amount of liquid nitrogen that is required to operate the two traps.





Figure 2-5. Front Panel - Lower Part



- 1 power indicator
- 2 host connection
- 3 status of turbo pump
- 4 error status of this pump
- 5 "Analyzer Pumps" switch (inoperable)

Figure 2-6. Control Panel in Close-Up View

The front panel allows to:

- control the oven temperature using the Jumo itron 16 temperature controller (2 in Figure 2-7, left). For detailed instructions on how to use the oven controller, see “Oven and Oven Control” on page 2-27.
- switch on the vacuum system via the control panel (Figure 2-6 and 1 in Figure 2-7, left)

**Note** The **Analyzer Pumps** switch, 5 in Figure 2-6, is inoperable with Kiel IV Carbonate Device. ▲

On the rear panel (Figure 2-7, right) the main power switch 5 is located. The connection for vent gas 4 and compressed air 3 can be found here as well.




Figure 2-7. Parts of Front Panel and Rear Panel

Figure 2-8 shows the fore pumps located at the rear panel.



- 1 fore pump for turbo pump
- 2 fore pump for fore vacuum side of Kiel IV Carbonate Device
- 3 compressed air reservoir used during disconnection of vials
- 4 inspection glasses for checking oil levels of fore pumps

**Figure 2-8.** Fore Pumps at Rear Panel



**Figure 2-9.** Electronics Cabinet

Behind the grey panel on the left side the electronics cabinet is located. See Figure 2-9. Access to the electronics parts is restricted to qualified service personnel only. The section is locked and additionally covered by a transparent pane.

**1** shows Inlet Control boards (Figure 6-38, two of them are arranged in stacked order). **2** is the Power Distribution board. **3** depicts the safety switch-off for the acid valve (Figure 2-10). **4** is the fiberline to the serial bus interface (Figure 2-11). **5** shows the interface to connect the proximity switches to the Inlet Control board (Figure 2-12).



**Figure 2-10.** Safety Switch-Off for Acid Valve

The safety switch-off (monoflop) for the acid valve is shown in Figure 2-10. It has been factory-preset to 3 min. The purpose of these monostable relays is to limit the acid flow in case of a computer failure.


**Caution** A defective pinch valve causes acid to overflow the oven section of the Kiel IV Carbonate Device! Even in case of smallest problems with the pinch valve, stop working immediately and exchange it! ▲



**Figure 2-11.** Fiberline to Serial Bus Interface


Figure 2-11 shows the fiberline to the serial bus interface.

Figure 2-12 displays the interface to connect the proximity switches to the Inlet Control board.



**Figure 2-12.** Interface Proximity Switches - Inlet Control Board

**Caution** Proximity switches can be operated up to a maximum temperature of 85 °C. ▲




**Figure 2-13.** Parts within Right Panel

Figure 2-13 and [Table 2-1](#) summarize some important parts to be seen behind the right panel.

**Table 2-1.** Important Parts within Right Panel


No.	Description
1	valve system of trap 1
2	vent valve (magnetic valve for fore vacuum connection, He/N <sub>2</sub> sample vial vent gas) It is closed when the pump is running. The pump controller keeps the valve closed. When the pump controller is switched off, the turbo pump acts as a generator and keeps the valve closed until pump speed decreases to less than 50 %.
3	turbo pump and its power supply (pump controller)
4	fore pump for backup of turbo pump
5	Autocool Unit for trap 1
6	Autocool Unit for trap 2
7	valve system of trap 2
8	Active Pirani Gauge (APG-M)

In the background of Figure 2-13 the fore vacuum gauge VM2 is located beneath the oven at the T piece. Figure 2-14 shows it in close-up view.



**Figure 2-14.** Position of VM2 Fore Vacuum Gauge

Additionally, Figure 2-15 shows V7 and V8 mounted on a distributor that connects to the fore pump, the vial vent gas and the oven section.



**Figure 2-15.** Valve 7 and Valve 8

## Connections

This section provides information about how to connect the Kiel IV Carbonate Device to the IRMS and to the gas supply.

### Connecting Kiel IV Carbonate Device to IRMS

To connect the Kiel IV Carbonate Device to the IRMS, proceed as follows:

1. Remove the side panels of the Kiel IV Carbonate Device.
2. Connect the compressed air connector which is located at the rear side of the Kiel IV Carbonate Device to the IRMS distributor. Use the quick release connection to connect the blue compressed air cable to the compressed air connectors of the IRMS. See Figure 2-5, right.

As the IRMS has four connectors, four screws (wing unions for compressed air, quick release connections) are provided either with the Kiel IV Carbonate Device or with the IRMS itself.

3. Connect the auxiliary gas (e.g. He or N<sub>2</sub>) to the connector at the rear side of the instrument. Set auxiliary gas pressure to less than 0.5 bar.
4. Connect the fiberline cable from this serial bus interface (Figure 2-11) to the rear panel connector of the IRMS.

**Note** If after switch-on the **Host Connection** LED is not on, interchange the grey and blue plugs of the optical fiber. ▲

5. Connect the mains cable to the appropriate connector of the IRMS.
6. Connect the fore vacuum exhaust to an appropriate vent.
7. Connect the reference refill tank capillary to the standard side of the IRMS (V22).
8. If not already active, start Isodat 2.5.
9. Open Instrument Control and close the Changeover Valve. See Chapter 6 of DELTA V Operating Manual.


10. Connect the stainless steel capillary from V3 and V4 (Microvolume of Kiel IV Carbonate Device) to the Changeover Valve of the IRMS by pulling it through the hole in the side panel of the IRMS.
  
11. Switch on **Inlet Pumps** at the front panel of the Finnigan Kiel IV Carbonate Device. After approximately 10 min, the green LED of the pump controller gets on, indicating that the vacuum system is operational.
  
12. Select the configuration for Kiel IV Carbonate Device and ensure that Isodat 2.5 can operate the valves. Refer to “[Accessories Bar and its Components](#)” on [page 3-5](#) for instructions.
  
13. Perform a leak test.
  
14. Set trap 1 and trap 2 to 150 °C. Keep V1, V2, V3, V4 and V5 open. Set oven temperature to 70 °C (Standby mode).
  
15. Prepare phosphoric acid<sup>1</sup> and transfer it to the acid reservoir of the Kiel IV Carbonate Device. See Figure 2-58. Place it inside the oven:
  - a. Take a 50 ml Luer-type syringe.
  - b. Add a silicone tube.
  - c. Put the syringe to the 1/16” tubing.
  - d. Remove any air by sucking it away from the phosphoric acid tubing line using the syringe.

**Note** The Viton tubings and the 1/16”-1/4” straight connectors must be leak-tight. Neither sucking of air nor air release should occur during dropping. ▲

16. Connect the two black Viton tubes of the acid valves to each of the two acid glass ports. See Figure 2-16. Check the acid flow by opening V12 and V22, and open acid valve 1 or 2, respectively.

---


<sup>1</sup>For instructions about phosphoric acid preparation, see “[Phosphoric Acid Preparation](#)” on [page 4-26](#).



**Figure 2-16.** Connecting Viton Tubes of Acid Valves


17. Take a magazine and place vials at positions 1/1 and 2/1.
18. Insert the magazine into the oven and perform **Load Magazine**. Refer to “[Load Magazine](#)” on [page 3-51](#).
19. Take the capillary heating transformer and connect the heating wires as shown in Figure 2-17 and Figure 2-18:
  - a. Connect one 7.7 V input of capillary heating transformer (positive pole, e.g. **2** in Figure 2-18) to the brass contact in center of one capillary (e.g. **2** in Figure 2-17).
  - b. Connect the other 7.7 V input of capillary heating transformer (positive pole, e.g. **3** in Figure 2-18) to the brass contact in center of the other capillary (e.g. **3** in Figure 2-17).

**Note** After the capillary has been wired make sure, the capillary (with or without insulation) has no contact to any plastic surface of tubes, housings, cables etc. to avoid melting or smoldering caused by a hot capillary. ▲



**Figure 2-17.** Fixing Alligator Crimps to Capillaries

20. Switch the capillary heating transformer on. Keep pumping for 12 h, that is over night. In the meantime, leave V1, V2, V3, V4, V5 and V9 open.



- 1 main switch (on/off)
- 2 input of 7.7 V (12 A)
- 3 input of 7.7 V (12 A)
- 4 input of 0 V
- 5 input of 0 V
- 6 automatic fuse  
Press to reset

**Figure 2-18.** Capillary Heating Transformer

21. Perform a vial test. Refer to “Vial Test” on page 4-23.
22. Wait until the acid is warmed up to 70 °C. Open the manual valve of the acid reservoir shown in Figure 2-2.
23. Open V10 and V20 while the two vials are connected and V13, V23 and V7 are open. See Figure 4-22.
24. Look into the vials: in the beginning, acid drops may spray out due to air bubbles. After this initial period however, the acid should drop slowly into the vials. It may take some time until the first drop of acid appears, because acid must fill the Teflon tube first.

**Note** The Viton tubings and the 1/16”-1/4” straight connectors must be leak-tight. Neither sucking of air nor air release should occur during dropping. ▲

For instructions about matching sample capillary to standard capillary, see “Matching Sample Capillary to Standard Capillary” on page 4-9.

## Vacuum System


The vacuum system of the Kiel IV Carbonate Device (see Figure 6-42) is located in the lower section. It consists of two different pump types:

- A turbo pump (TMH 071 P, manufacturer: Pfeiffer, **3** in Figure 2-13 and Figure 2-19) is connected to the trapping section of the Kiel IV Carbonate Device, that is trap 1 and trap 2.
- A rotary vane pump (DUO 2.5, manufacturer: Pfeiffer, **1** in Figure 2-8) provides the fore vacuum for the turbo pump.
- A second rotary vane pump of the same type (DUO 2.5, manufacturer: Pfeiffer, **2** in Figure 2-8) is used during connection and disconnection of sample vials.

The fore vacuum generated by this second fore pump of the Kiel IV Carbonate Device is monitored by an Active Pirani Gauge (**8** in Figure 2-13, APG-M, manufacturer: Edwards<sup>1</sup>) located in the lower cabinet.

---

<sup>1</sup>See Instruction Manual of BOC Edwards and [www.bocedwards.com](http://www.bocedwards.com).




- 1 electrical connection to Power Distribution board
- 2 LEDs showing pump controller status
- 3 vent valve

Figure 2-19. Turbo Pump with Turbo Pump Controller

## Valves and Traps

This section provides information about the valves and traps of the Kiel IV Carbonate Device. See Figure 2-20.



- 1 trap 1 assembly
- 2 trap 2 (Microvolume) assembly

Figure 2-20. Trap Arrangement

## Pneumatic Valves

**Layout** A cylinder on top is actuated by compressed air. Its gold-made plunger then presses a membrane underneath and thus tightens. Gas transfer is then impossible. When no compressed air is present, the cylinder is not actuated. Its plunger will not press, and a spring assembly resets the membrane and thus does not tighten. Gas transfer is possible.

**Parts of a Pneumatic Valve** Figure 2-21 depicts the parts of a pneumatic valve.



Figure 2-21. Parts of a Pneumatic Valve

**Inserting a Pneumatic Valve** Figure 2-22 shows the pliers **1** to properly insert a pneumatic valve **2**. Hold the pneumatic valve tight by jamming it within the pliers (left in Figure 2-22). Then shove the outer sleeve above the pneumatic valve (right in Figure 2-22).





Figure 2-22. Pliers to Insert a Pneumatic Valve

## Valves of Dual Inlet System

The Dual Inlet system is operated by pneumatic valves with a nominal closing pressure of 4 bar. Even though they are all made of stainless steel, after long-term operation they might be worn nevertheless. All the valves base upon the same construction principle. Figure 2-23 shows such a valve with its high-vacuum side opened. Open the pneumatic valve and grease the plunger without disassembling it.

## Parts of a Valve



The stainless steel membrane 4 is turned down by 180° and then laid onto the gold-made gasket 2, that is with the plunger 5 oriented downwards.

- 1 valve block, made of stainless steel
- 2 gasket (gold), seals valve block against stainless steel membrane
- 3 Swagelok connector as gas inlet; laterally welded on the valve block
- 4 stainless steel membrane with valve plug that closes the valve. See also Figure 2-22.
- 5 plunger, made of gold; fits exactly to the knife edge
- 6 knife edge, located in center of valve block; also gas exhaust
- 7 hole, acting as gas inlet to the valve connected to 3 via an internal bore

**Figure 2-23.** Dual Inlet System Valve



**Figure 2-24.** Parts of Dual Inlet System Valve

The parts of a Dual Inlet system valve are also made of stainless steel. They are depicted in Figure 2-24 and summarized in [Table 2-2](#).

**Table 2-2.** Parts of Dual Inlet System Valve\*

No.	Description
1	valve block (made of stainless steel)
2	stainless steel membrane with valve plug that closes the valve (a sleeve not to be seen in Figure 2-24 is attached to its rear side)
3	actuator for compressed air (usually lying within rear side) Refer to <a href="#">"Pneumatic Valves"</a> on <a href="#">page 2-18</a> .
4	covering cap
5	Four screws to fasten covering cap

\*See Figure 2-24.

### Arrangement in Valve Blocks

Compressed air is either supplied by an optional compressor attached to the IRMS or by a user supplied pressure air line. The metal valves are equipped with gold gaskets and gold seats acting on knife edges. Up to six valves are machined into one monoblock, thus considerably reducing the volume in plumbing as well as possible leak of the installation.

This type of valve block is used throughout all inlet modules. For plumbing the valve blocks are fitted with special 1/4" Swagelok connectors. Compressed air is fed to the pneumatically operated valves by solenoid valves. These are controlled by dedicated electronics linked to the computer via a data bus.





Figure 2-25. Double Valve Block

### Manifold Block with Solenoid Valves

Four of the solenoid three-way valves are located on a manifold block. The solenoid valves are operated with a voltage of 24 V. The voltage is supplied by the Dual Inlet board.




- 1 four LEDs reveal individual switching status (on/off)
- 2 connector to Dual Inlet board

Figure 2-26. Manifold Block with Four Solenoid Valves

The solenoid valves are normally open (with the exception of TubeCracker, Part No. 108 2840). The working condition is signaled by a red LED located on the board. The actuators for compressed air transform a signal A into another signal B: they switch an electrical signal generated at the Dual Inlet board into a compressed air signal. Thereby, compressed air is provided which forcefully switches the actual valves.

**Note** In case of a power failure, all solenoid valves open automatically (with the exception of TubeCracker, where they close automatically). Thus, the pneumatic valves in the entire Dual Inlet system close avoiding its contamination. ▲

## Compressed Air Distributor




**Figure 2-27.** Compressed Air Distributor

The compressed air distributor, Figure 2-27, is part of the Dual Inlet system. Therefore, it is missing if no Dual Inlet system is available. The compressed air connections of the distributor are all equivalent ones. Compressed air enters at **1** in Figure 2-27 and is then distributed to all compressed air valves.

The number of compressed air valves of the system depends on which particular options for the Dual Inlet system are available. If only few compressed air valves must be connected to the compressed air distributor, one or more blind plugs **2** allow to close the unused connections tightly.

## Trap 1

This trap consists of a trapping volume valve block and an Autocool Unit. See Figure 2-28. This trapping volume has two tubes, which are stuck into each other. The inner tube must be connected to trap 2. The other end is open and hangs inside the outer tube. See Figure 6-19. The outer tube is connected via valves 1 and 12 (22) to the carbonate preparation vial.



**Figure 2-28.** Schematic of Trap 1\*

\*The arrow points in counterflow direction.

By means of the Autocool Unit, temperature can be set individually, if required, within a range of about -196 °C to +150 °C. Before any analysis of carbonate the computer sets the trap to +150 °C in order to remove any impurities. After the cleaning procedure, the Autocool Unit cools the trap to -196 °C. As soon as acid drops into the vial which contains carbonate, the reaction takes place and CO<sub>2</sub>, H<sub>2</sub>O, N<sub>2</sub> and O<sub>2</sub> are present in different concentrations. Since the trapping volume is at -196 °C, the CO<sub>2</sub> and H<sub>2</sub>O gases get frozen inside the outer tube. This means that they are isolated and separated from other produced gases.

After pumping O<sub>2</sub>, N<sub>2</sub> as non-condensable gases, Autocool sets the trap temperature to -120 °C to -115 °C. At this temperature, CO<sub>2</sub> as gas in a vacuum is released into the inner tube for later use. The release temperature of -120 °C to -115 °C is extremely important. If temperature is too high, water can get trapped by trap 2. If temperature is too low, CO<sub>2</sub> might be retained. The punched arrow of the trap does not show the gas direction, but counterflow direction. It indicates the installation direction. See “Lower Section” on page 2-3.


**Note** If the trap finger has to be exchanged, proper gold gaskets and flow direction of the inner and outer trapping volume must be guaranteed! Any scratch or the wrong direction influences trapping efficiency and incorrect CO<sub>2</sub> trapping, that is isotope fractionation. ▲

## Microvolume (Trap 2)

In Dual Inlet mode, measurement is performed out of a volume for smallest amounts, the so-called “Microvolume”. Figure 2-29 shows the Microvolume and its parts<sup>1</sup>. It consists of a trapping volume, a valve block (valve 3 and 5), an Autocool Unit and a separate capillary. This capillary leads directly to the Changeover Valve of the IRMS.

**1** reminds of a specially designed valve block. The trapping volume **2**, a drilled tube with a clamped gold seal on top, comprises a pre-set, included volume of 50 µl. Its lower part is cooled to the temperature of liquid nitrogen, e.g. -196 °C so that the gas to be measured will be frozen therein.

**3** is a copper-made heat transfer tube that surrounds the trapping volume. The surface of the tube has been gold-plated. K indicates the connection for the stainless steel capillary. V on its rear side indicates the Swagelok socket for the vacuum connection.



**Figure 2-29.** Microvolume and its Parts\*


\*The numbers refer to Table 6-6.

The total volume in front of the capillary crimp including trapping volume and the capillary is around 145 µl. Due to the viscous flow conditions which require a pressure of at least 15 mbar for CO<sub>2</sub> gas in front of the capillary, a sample of 5 to 50 mbar/ml has to be concentrated to a small volume. The concentration in a Microvolume is performed by freezing CO<sub>2</sub> gas from trap 1 using liquid nitrogen and expanding it by subsequent heating. By means of the Autocool Unit, the temperature can be set individually within a range of about -190 °C and +150 °C, if required.

<sup>1</sup>“Microvolume” is a collective term for all the parts displayed in Figure 2-29.

## Hardware Components

### Valves and Traps




**Figure 2-30.** Microvolume Parts to be Inserted into Autocool Unit

## Autocool Unit

Figure 2-31 shows the dismantled Autocool Unit. Under operating conditions, its funnel will be completely immersed into liquid nitrogen.

The funnel contains a heater to be turned on and off. If a current flows through it, bubbles will ascend into the tube and thus transport liquid nitrogen upwards to the cascade. The transport direction of liquid nitrogen is indicated by the arrow in Figure 2-31.




**Figure 2-31.** Autocool Unit and Interior of Funnel

The cascade with its three stages is depicted in Figure 2-32. Liquid nitrogen is transported towards the uppermost cascade and then drips through holes to the cascade below. Thereby, the unit is cooled to the temperature of liquid nitrogen, e.g.  $-196\text{ }^{\circ}\text{C}$ .


As the cascade contains a heater and a temperature sensor, the temperature can be raised above this temperature and adjusted there very precisely within a wider range.

The biggest hole, **1** in Figure 2-32, takes up the heat transfer tube (**3** in Figure 2-29). Gilding its cupreus surface facilitates pulling it out of **1** after longer periods of operation for servicing or repair.

The heater cartridge fits into hole **2**. The smallest hole **3** contains a temperature sensor resistor (Pt 100).



**Figure 2-32.** Cascade with Three Holes



**Figure 2-33.** Mounting Instructions for Autocool Unit onto Trap 1

In order to remove the cooling/heating cascade with funnel:

1. Loosen screws 1 and 2 in Figure 2-33 and the electronic connection to the electronic board. Trap numbers are indicated on the cable.
2. Push the clamp to the left (or right). Then, gently push down the Microvolume bar.

**Warning** Do not loose the Teflon 30 mm washer! It must be reattached at the same position. ▲





**Warning** Careless handling of liquid nitrogen might cause personal injury including frostbite. Wear protective clothing when operating this equipment including protective gloves and face shield. ▲

The temperature of the Autocool Unit which cools trap 1 or trap 2 can be set by a routine in Isodat 2.5's Instrument Control. If, e.g. the temperature level is set to  $-80\text{ }^{\circ}\text{C}$ , the heater works against the temperature of the liquid nitrogen in order to keep the set temperature.

A temperature between  $-196\text{ }^{\circ}\text{C}$  and  $+150\text{ }^{\circ}\text{C}$  can be set. Trap 1 or trap 2 fit into a thermal contact attached to the lid of the dewar. The dewar contains liquid nitrogen up to a certain level. Fitted to the contact pipe are an electrical heater element, a temperature sensor and cascade arrangement of three small bowls. All parts of the assembly are made of a high thermal conductivity material and are placed in close thermal contact to each other. This achieves a quick changeover from one temperature to another.

To heat trap 1 or trap 2 to a defined temperature, the heater element is activated, and the temperature sensor controls the heating phase. To cool trap 1 or trap 2, another electrical heater element immersed in liquid nitrogen is activated and causes evaporation as well as agitation.

Above the latter heater element, a funnel-shaped hood of standpipe is positioned which leads to the uppermost bowl of the cascade. This arrangement enables about one drop of liquid nitrogen per second to be carried by the stream of evaporated nitrogen. Small holes in the bottom of the bowls enable a constant trickle of liquid nitrogen back into the dewar, and the continuous flow of liquid nitrogen rapidly cools down trap 1 or trap 2. By suitable balancing of the liquid nitrogen flow and heating trap 1 or trap 2, any temperature within the temperature range of  $\pm 2\text{ }^{\circ}\text{C}$  can be attained.



A dewar serves as reservoir for liquid nitrogen. It is adjusted beneath the trap assembly using a "lab boy".

**Figure 2-34.** Dewar and "Lab Boy"

## Oven and Oven Control

The Jumo itron 16 temperature controller allows to control the temperature of the oven. Notice its three keys indicated by the arrows in Figure 2-35.



Key	Function
P key	for programming. Values will be accepted automatically after 2 s
Arrow Up key	to increase a particular value
Arrow Down key	to decrease a particular value

**Figure 2-35.** Oven Control

### Programming - Step 1

1. Press the P key and hold it for 2 s.
2. Pass through the menu until Y.0 is displayed.
3. Again, press the P key and hold it for 2 s.
4. Set the parameters according to [Table 2-3](#).

**Table 2-3.** Parameters for Step 1 of Programming

Parameter	Value	Explanation
C111	003	temperature sensor, that is transducer type (e.g. NiCr-Ni, K)
C112	1	number of decimal places and temperature unit (e.g. 1 and °C)
C113	33	controller type (e.g. double setpoint)
C115	0	ramp function, that is ramp function is off
C116	0	outputs on fault, that is 0 % minimum output limiting Y.2 is effective
SPL	0	lower setpoint limiting
SPH	100	upper setpoint limiting
OFFS	0	process value correction

For details, refer to Jumo itron 16 temperature controller manual.

## Programming - Step 2

1. Again, press the P key and hold it for 2 s.
2. Press the Arrow Up key and Arrow Down key to change values.
3. Set the parameters according to [Table 2-4](#).

**Table 2-4.** Parameters for Step 2 of Programming

Parameter	Value	Explanation
Pb.1	0.4	proportional band 1
Pb.2	0.2	proportional band 2
d.t	2	derivative time [s]
r.t.	11	reset time [s]
CY.1	2.3	cycle time 1 [s]
CY.2	2.3	cycle time 2 [s]
db	0.0	contact spacing
HyS.1	1.0	differential 1
HyS.2	1.0	differential 2
Y.0	0	working point [%]
Y.1	85	maximum output [%]
Y.2	0	minimum output [%]
d.F	1.1	filter time constant [s]

For details, refer to Jumo itron 16 temperature controller manual.


## Alternative - Automatic Programming

Let the temperature controller program itself automatically. Thereby, you don't need to specify all the parameters mentioned above on your own. For details, refer to Jumo itron 16 temperature controller manual.

## Tune Scan after Programming

**Note** After both programming steps have been finished, perform an Autotune. This requires to set the ramp function C115 to 0 before you press the Arrow Up key and the Arrow Down key simultaneously! During Autotune, the controller displays the setpoint value and “tune“ alternating. ▲

## Autosampler



**Figure 2-36.** Turret and Correct Vertical Placement of Vials

## Installing and Removing Magazine

**Caution** Be careful not to destroy the acid dropper capillary and the drop counter spring during installation of the removable tray! ▲

The whole magazine turret consists of the magazine part containing the vials and a cover plate with two holes and a locator to install it inside the oven the specific way shown below. See Figure 2-36, Figure 2-37 and Figure 2-38.



**Figure 2-37.** Turret and Lid for Line 1 and Line 2 Position

The turret can be removed from the oven to load/unload it with vials. For this purpose, a special stand comes with Kiel IV Carbonate Device.



**Warning** Be sure to put the stand in a safe place before you remove the hot turret from the oven! ▲


In order to remove the turret, stop any sequence and press **Take Magazine** in Instrument Control. After loading the turret, you can put it back to the oven. Then, press **Load Magazine** in Instrument Control.

## Adjusting Magazine Position

To adjust the position of the vials in the magazine to the piston position, perform as follows:

1. Take away the cover of the magazine and place the magazine inside the oven.
2. Loosen the screws 1-4. Move the magazine. Tighten the screws.
3. Rotate the magazine until the vials are aligned to the piston.
4. Take away the magazine and tighten the screws 1-4 again.
5. Rotate the magazine several times.
6. Perform a **Take Magazine** procedure and a **Load Magazine** procedure to check positioning again.

Moving the support rod can align the cover plate to the magazine. The support rod is fixed by two screws, which are located at one end.




**Figure 2-38.** Magazine with Cover Plate and Vials


## Adjusting Piston Speed

The speed of the piston to move the vials up- or downwards can be adjusted as follows (see Figure 2-39 and Figure 2-40):

- Turn the screw **clockwise** to **reduce** the speed of the piston movement and **counterclockwise** to **increase** it.
- The screw at the **upper** end of the piston is used to move the piston **downwards!** The screw at the **bottom** of the piston is used to move the piston **upwards!**



**Figure 2-39.** Piston Speed Adjustment



- 1 compressed air lever for piston line 1
- 2 compressed air lever for piston line 2
- 3 magnetic position sensor for piston with LED to indicate position
- 4 flow limiter
- 5 compressed air connection

**Figure 2-40.** Connections of Pneumatic Levers for Pistons


Figure 2-40 displays the connections of the pneumatic levers for the two pistons.

## Adjusting Piston Height

Inside the piston, a spring assembly is located to ensure proper pressure when connecting vials. It may be necessary to adjust the overall height to ensure proper operation of the connect process.

For this purpose, the piston height can be adjusted. The piston itself is screwed into the piston compressed air lever and secured by a second nut. See Figure 2-41.

This position is difficult to access. The screw can be found on top of the pneumatic lever, at the bottom of the piston itself.



- 1 compressed air lever
- 2 safety nut
- 3 bottom of piston


**Figure 2-41.** Adjusting Piston Height

A red PTFE vial bottom lever is located inside the piston. If phosphoric acid has spoiled into the piston, this lever can latch vials. This leads to a hardware failure.

**Note** In this case, remove the entire piston and clean it thoroughly, also on its inside. ▲

## Proximity Switch

The proximity switch used to detect the presence of vials is located inside the acid valve. This switch is a rod which at one end contains a coil as shown. The electronic switch contains no mechanical parts. Once the piston moves the vial to the acid valve, the vial pushes a metal U-shape spring upwards close to proximity switch. See Figure 2-42. The induction of the coil changes, which means the vial is connected.



- 1 window to check sensor status
- 2 adjustment screw
- 3 electrical connection of sensor

Figure 2-42. Proximity Switch

### Turret Motor





Figure 2-43. Turret Motor and Position Sensing Array

The motor to drive the turret is located beneath the oven, together with the position sensor for the autosampler section. The setup consists of an array of infrared light barriers and a plate with holes that is binary coded. Depending on the position of the rotating plate, more or less of the light barriers are closed resulting in a code that represents the position of the turret. See Figure 2-43.




In case the position sensor shown in Figure 2-43 has eventually shifted, use the tool (Part No. 115 7390) to readjust its position.

**Figure 2-44.** Tool to Readjust Array of Infrared Light Barriers

## Liquid Nitrogen Refill Device

This section contains information to be read before operating the refill device. Read the manufacturer's handling instructions carefully as well.




**Figure 2-45.** Schematic of a Standard Liquid Nitrogen Refill Device



**Warning** The refill device contains an extremely cold liquid. Careless handling might cause personal injury including frostbite. Wear protective clothing including protective gloves and a face shield. ▲



**Warning** Do not overfill or tilt the refill device and prevent spills. Use the refill device only in well-ventilated areas. Poor ventilation causes suffocation. Keep in mind: Safety first! ▲



- 1 liquid fill and decant valve
- 2 vent valve
- 3 pressure gauge
- 4 fill level display
- 5 connection to transfer tube
- 6 magnetic valve for refill (open/close)
- 7 5 bar LN2 overpressure gauge for overpressurized arm next to 6
- 8 central O ring seal to close dewar

**Figure 2-46.** Magnetic Valve and Liquid Nitrogen Safety Unit



**Figure 2-47.** Example of a Liquid Nitrogen Refill Device

All liquid nitrogen reservoir tanks contain some kind of level indicator, a pressure gauge to read the internal pressure, safety pressure relieve valves and a mechanism to build up pressure. This mechanism can usually be switched off (pressure raising valve).

Furthermore, there are connections to decant the liquid phase or the gas phase, and usually there is a connection to vent the reservoir. Finally, there is a pressure regulator to adjust the reservoir pressure.

For proper operation of the refill device, a pressure of about 0.7-1.4 bar (10-20 psi) is required inside the tank to ensure liquid nitrogen flow, if the decant and refill valves are opened. For this purpose, all liquid nitrogen storage devices have a pressure raising valve and a pressure limiter.

**Note** If the central O ring seal, **8** in Figure 2-46, is not closed properly, the Kiel IV Carbonate Device will not operate, because no pressure raise is possible! ▲

## **Warning Note for Liquid Nitrogen Supply**

Concerning the automatic refill device, you may have either a liquid nitrogen tank of your own or a tank delivered by Thermo Electron optionally (30 l or 90 l).

As a part of the tank, a **liquid fill and decant valve** (**1** in Figure 2-46) is also included. It is a **manual cutoff valve** to close and open the tank. During tank transport it must be closed. After the tank has been connected to the gas line, the valve must be opened to provide nitrogen for the entire system. Therefore, never close it during operation!

The **solenoid valve** (that is **magnetic valve**, **6** in Figure 2-46) is controlled by Isodat 2.5. Its status depends on the processes inside the device: when nitrogen is required, it will be opened allowing new nitrogen to flow along the tube. When no more nitrogen is needed, it will be automatically closed. The manual cutoff valve may still be open.



**Warning** The tube between solenoid valve and manual cutoff valve may never be closed! When liquid nitrogen inside the closed tube is warmed, pressure will increase considerably! ▲



**Warning** If the manual cutoff valve of the liquid nitrogen tank and the solenoid valve are closed simultaneously, the liquid nitrogen distributor as well as the gas line to the tank could burst and seriously injure operating staff! ▲


Therefore, together with the distributor a pressure control valve is delivered by Thermo Electron. This safety valve serves to reduce excess pressure and thereby prevents the tube from bursting. It must be installed by a technician. In case of an **upgrade**, perform as follows:

1. Unscrew the distributor's dummy plug.
2. Screw in the pressure control valve.


In case of a **new system**, the pressure control valve has already been screwed onto the distributor.




**Warning** Always mount the 5 bar safety valve between solenoid valve and manual cutoff valve! Never operate the instrument without safety valve, that is never unscrew it! In case the 5 bar safety valve is damaged, renew it immediately! ▲



**Figure 2-48.** Valves and Distributor for Liquid Nitrogen (Top View)



**Figure 2-49.** Valves and Distributor for Liquid Nitrogen (Side View)




**Figure 2-50.** Valves and Distributor for Liquid Nitrogen

## **LN2 Transfer from Refill Device into Dewar**


The refill device provides a constant level of liquid nitrogen in the dewar for the Autocool Unit. It is equipped with the necessary safety devices, valves and pressure gauges required to safely handle liquid nitrogen.

A solenoid-operated refill valve controls the transfer line to the dewar of the Autocool Unit. This refill valve must be directly connected to liquid fill and decant valve of the refill device. A level sensor installed in the dewar of the Autocool Unit activates the refill device.

The refill sensor used for liquid nitrogen refill can be found besides the trap assemblies. See Figure 2-51 and Figure 2-52.



**Figure 2-51.** Position of Refill Sensor



**Figure 2-52.** Liquid Nitrogen Refill Sensor

### Functional Principle of Refill Sensor

The refill sensor is made up from a volume filled with silica gel, nickel wool and air that is connected leak-tight to a pressure gauge. If the volume is at room temperature, the inside of the tube is at atmospheric pressure. If the tube is emerged in liquid nitrogen, the nitrogen inside the tube will condense on the surface of the silica gel, and the pressure will drop below 10 mbar. This pressure change is recorded by the gauge and measured by Isodat 2.5.

Whenever the temperature of the traps is set to less than 20° C, the state of the sensor will be checked and, if necessary, liquid nitrogen refill will be initiated.

### Possible Problems

- Liquid nitrogen is refilled only on request (via Isodat 2.5).
- If the sensor is defect or improperly adjusted, overflow may occur.
- If the sensor reacts too sensitive, cooling to temperature of liquid nitrogen is no longer possible. Change the sensor position by pulling it out of the dewar.
- Setting the -120 °C level seems to be problematic, if the sensor position is too high.

## Reference Gas Refill

When working with the Finnigan Kiel IV Carbonate Device, a reference gas refill (synonymously called standard gas refill) is necessary to avoid running out of reference gas during measurements.



**Warning** When installing CO<sub>2</sub> reference gas tanks, keep in mind that standard high pressure tanks for CO<sub>2</sub> contain a liquid phase that is subject to fractionation when temperature changes. These tanks must be stored at constant temperature to obtain stable isotope values for your reference gas. ▲




**Warning** CO<sub>2</sub> from a high pressure gas tank is not suitable! ▲

Reference refill provides the reference gas supply to the inlet system. See Figure 2-53 and Figure 2-54. It is a hardware option and consists of a metal tank of approximately 5 l with a manual valve connected via a capillary (see Figure 6-7) to one of the inlet ports on the standard side. With the reference refill selected, the standard side of the Dual Inlet system is completely pumped out before it is filled again for the next measurement.



**Figure 2-53.** Reference Gas Refill Units

Figure 2-53 shows two reference gas refill units, the left one, **1** with a capillary for normal operation and the right one, **2** with a bigger tube (e.g. 1/4”) suitable for baking and refilling operation. Both capillary and tube are included in the Reference refill option.



**Figure 2-54.** Manual Valves

Figure 2-54 shows manual valves. The left one, **1** denotes the manual valve to close the refill tank in open position. The right one, **2** denotes the manual valve in closed position. The reference refill parameters can be set at “Reference Refill” on page 3-23.


### Filling Reference Refill Device from External Source



**Figure 2-55.** Filling Reference Refill Device from External Source - I

## Hardware Components

### Reference Gas Refill




**Figure 2-56.** Filling Reference Refill Device from External Source - II


Whereas in Figure 2-55 the connections are outlined, Figure 2-56 shows the switching of the valves. To fill the Reference Refill Device from an external source, proceed as follows according to Figure 2-56:

1. Evacuate the Reference Refill device over night using the high vacuum pump.
2. If you fill from a high pressure tank with pressure regulator, adjust the manometer of the CO<sub>2</sub> tank to 0.6 bar. In order to avoid contamination, flush the manometer for about 5 min before connecting the line to the extern left.
3. Evacuate the entire vacuum line first using the fore vacuum pump and then using the high vacuum pump. See arrows in Figure 2-56.
4. Open the tank connection and allow gas from the tank to fill the Reference Refill device. Equilibrate for 15 min before breaking the connection. In order to obtain an ideal pressure of about 0.6 bar (8.5 psi) in the Reference Refill device, adjust about 0.6 bar (8.5 psi) at the pressure regulator of the CO<sub>2</sub> tank.

## Acid Flow



**Figure 2-57.** Schematic of Acid Flow




**Figure 2-58.** Acid Reservoir and Manual Valve

## Hardware Components

### Pinch Valve

Figure 2-57 schematically shows the acid flow of the instrument. For its upper part, see also Figure 2-58. For its lower part, see also Figure 2-59.



- 1 Viton tubing from reservoir
- 2 pinch valve
- 3 "acid valve"
- 4 drop counter connection

**Figure 2-59.** Pinch Valve and Drop Counters

## Pinch Valve

The pinch valve is shown in Figure 2-60 (mounted), Figure 2-61 (disassembled), Figure 2-62 (with acid tubing) and Figure 2-63 (connected to acid valve).



- 1 drop counter assembly
- 2 drop counter electrical connection

**Figure 2-60.** Pinch Valve - Mounted



**Figure 2-61.** Pinch Valve - Disassembled



**Figure 2-62.** Feeding Acid Tubing into Pinch Valve



**Figure 2-63.** Mounting Pinch Valve upon Acid Valve



## Chapter 3 Isodat 2.5

This chapter outlines the following topics:

- “Creating a Kiel IV Carbonate Device Configuration” on page 3-2
- “Acquisition Mode” on page 3-4
- “Accessories Bar and its Components” on page 3-5
- “Creating a New Method” on page 3-20
- “Structure of Methods for Kiel IV Carbonate Device” on page 3-21
- “Creating a New Sequence” on page 3-31
- “Standards” on page 3-38
- “Excel Export” on page 3-42
- “Dyn External” on page 3-44
- “Service Scripts” on page 3-49
- “Interpreting Results” on page 3-56
- “Time Slicing” on page 3-60
- “Interfering Masses” on page 3-64

# Creating a Kiel IV Carbonate Device Configuration


**Note** Do not create a Kiel III Carbonate configuration! ▲

To create a new Kiel IV Carbonate Device configuration:




Start Isodat 2.5's Configurator.

**Figure 3-1.** Opening Configurator



Click the **Add Configuration** button.

**Figure 3-2.** Adding a New Configuration




Replace **New Configuration** by a significant name, e.g. **Kiel IV Carbonate**.

The new configuration contains the IRMS you previously specified, e.g. Delta V Plus.


Click on its + symbol.

**Figure 3-3.** Renaming the New Configuration




The tree of the IRMS will be expanded. The ports **Source** and **Capillary** become visible.

**Figure 3-4.** Expanding Tree of IRMS



**Figure 3-5.** Appending Kiel IV Carbonate Set to IRMS

- a. Click on the **Dual Inlet Sets** tab.
- b. Among the Dual Inlet Sets in the right pane, mark the one that contains the Kiel IV Carbonate Device, that is **Dual Inlet+Kiel IV Carbonate+Reference Refill**.
- c. Drag&drop it to **Capillary** port in the left pane.



The entire Kiel IV Carbonate Device has been appended to the capillary port.

Clicking on the + signs unfolds the tree showing the individual components of the device.

Close the Configurator. All settings will be saved automatically.

**Figure 3-6.** Kiel IV Carbonate Device Appended to Capillary Port

## Acquisition Mode

This section outlines Acquisition Mode. For detailed information, refer to:

- ISODAT *NT Operating Manual*, Part No. 109 2481
- *ISODAT NT Operating Manual - Upgrade to Version 2.0*, Part No. 115 4991

## Starting Acquisition Mode



Start Isodat 2.5 by a double-click.




Start Acquisition Mode.

You are now able to run any application which gives you full control over the automated measurement.


## Activating Toolbars

**Note** It is recommended to check first, whether the following toolbars (that is dialog bars) are activated. Proceed as follows. ▲




1. Move your cursor to the title bar **Isodat Acquisition**.

2. Right-click on it.



3. Select **Properties**.




4. Select the **Bars** tab.
5. Select the toolbars to be displayed. We recommend to primarily select **Status** bar, **Basic** bar, **Accessories** bar and **Information** bar.
6. Confirm by **OK**. The bars will appear in the **Acquisition** window

**Figure 3-7.** Visibility of Individual Toolbars \*

\*The individual bars mentioned in Figure 3-7 are described in detail in the ISODAT NT Operating Manual, Part No. 109 2481.

## Accessories Bar and its Components




1. In order to display the **Accessories** bar at all, mark the corresponding checkbox at **Bars** tab.

**Note** It is important that you have already created a configuration that contains the Kiel IV Carbonate Device, e.g. **Kiel IV Carbonate**. Refer to “[Creating a Kiel IV Carbonate Device Configuration](#)” on page 3-2. ▲

## Isodat 2.5

### Accessories Bar and its Components



The **Status** bar displays


- the actual configuration (e.g. **Kiel IV Carbonate**)
- the actual Gas Configuration (e.g. **CO2**).

2. Select just this configuration.

**Figure 3-8.** Status Bar

The **Accessories** bar displays the selected configuration together with its configured features (as previously defined in the Configurator). See Figure 3-9.

For detailed information about the components of the **Accessories** bar, refer to the *ISODAT NT Operating Manual*, Part No. 109 2481.



The **Accessories bar** contains information about:

- 1 Quick Access buttons
  - 2 Readout of high voltage, magnet current and ion source pressure (MS)
  - 3 MS State of heaters
  - 4 Voltage settings in the ion source (Focus Delta)
  - 5 Dual Inlet window
- See "Dual Inlet Window" on page 3-8.
- 6 Kiel IV Carbonate window
- See "Kiel IV Carbonate Window" on page 3-8.
- 7 Object Properties
  - 8 A lot of running scripts are displayed (ISL Scripts)
  - 9 File Browser


**Figure 3-9.** Components of Accessories Bar

## Changing Visibility of its Components

Change visibility of individual components of the Accessories bar as follows:



1. **Right-click** on an arbitrary title bar (e.g. **Kiel IV Carbonate**).
2. Click the **Administrate Panels** button.



3. Mark the information to be displayed on the Accessories bar, e.g. Kiel IV Carbonate.
4. Unmark the information not to be displayed on the Accessories bar, e.g. ISL Scripts.
5. Confirm by **OK**.


**Figure 3-10.** Marking or Unmarking Accessories

## Kiel IV Carbonate Window

The **Kiel IV Carbonate** window is shown as **6** in Figure 3-9. It will be discussed in detail at [“Elementary Handling of Kiel IV Carbonate Device”](#) on page 4-22.

## Dual Inlet Window

The **Dual Inlet** window is shown in Figure 3-11 and as **5** in Figure 3-9.




**Figure 3-11.** Dual Inlet Window

**Table 3-1.** Indications of Dual Inlet Window\*

No.	Indication	Refer to
1	actual pressure of left side of Dual Inlet system (in mbar)	pressure transducer, 1 in Figure 4-2
2	actual pressure of right side of Dual Inlet system (in mbar)	pressure transducer, 1 in Figure 4-2
3	optical reminder of the crimp position of a capillary	
4	actual fore vacuum pressure in mbar measured by a Pirani gauge (as fore vacuum gauge of Dual Inlet system)	
5	volume proportion of left bellow (in %)	
6	volume proportion of right bellow (in %)	
7	Dual Inlet system fore pump	"Dual Inlet System Fore Pump" on page 2-30
8	Dual Inlet system turbo pump	"Dual Inlet System Turbo Pump" on page 2-27
9	Additional valve	"Additional Valve" on page 2-30

\*Refer to Figure 3-11.

### Operation of Changeover Valve




**Figure 3-12.** Switching Positions of Changeover Valve

Changeover Valve can be switched to three different positions described in Figure 3-12 and Table 3-2:

**Table 3-2.** Switching Positions of Changeover Valve

No.	Designation	Comment
1	Changeover Left	The left side capillary is opened to the ion source while the right side capillary is opened to Dual Inlet system turbo pump.
2	Changeover Right	The right side capillary is opened to the ion source while the left side capillary is opened to Dual Inlet system turbo pump.
3	Changeover Closed	Both capillaries are pumped by Dual Inlet system turbo pump.


Usually, the Changeover Valve is controlled by Isodat 2.5 in order to accomplish a Dual Inlet measurement. A flowchart of this basic measurement is shown in Figure 3-13.




**Figure 3-13.** Dual Inlet Measurement Loop

Two important time constants can be adjusted in the Dual Inlet method:

- Idle time (Figure 3-14, is **Pre Delay** in Figure 3-13) and
- Integration time (Figure 3-15)



**Figure 3-14.** Setting Idle Time




**Figure 3-15.** Setting Integration Time

### Shot Noise Limits of Precision in Dual Inlet Measurements

Figure 3-16 displays how precision (one  $\sigma$  value) varies when integration time and amplitude are changed. Higher amplitudes and longer integration times result in enhanced measurement precision (that is low standard deviation of the  $n$  repetitions with  $n$  selected at **Number of Cycles** in Figure 3-14, e.g.  $n=8$ ).

Thus, it can be used to select a reasonable integration time for a given measurement and to calculate the precision that can be expected.



**Figure 3-16.** Shot Noise Limits of Precision

The diagram contains the results of a calculation for the shot noise (statistical noise) on a cup, taking into account the integration time, the resistor of the cup and the signal height in that cup. The sample calculation is taken out for the middle cup and thus is true for  $\delta^{13}\text{C}$ . For  $\delta^{18}\text{O}$  roughly multiply the results by 1.4.

The same mathematics is used with the calculations that are contained on the **All Products** CD supplied by our marketing department.

When we talk about precision here, we call it "internal precision". This number is reported as "Standard Deviation" in the output grid. Standard deviation  $\sigma$  ("Std. deviation" column in Figure 3-76) is given by:

$$\sigma = \frac{1}{\sqrt{n-1}} \cdot \sqrt{\sum_n (\delta - \delta_{\text{mean}})^2}$$

where n denotes the number of repetitions.

**Note** This is not the number reported as "standard error" although in discussion both are often mixed up. ▲

The "Standard Error" is a number that is generally smaller because it takes into account the repetitions of the measurement (number of cycles - usually set to 8). The standard error is also reported in the output grid and represents the error in determining the average of a distribution (in

our case the average  $\delta$  value of the  $n$  repetitions of individual measurements). Standard error  $\sigma_e$  (“Std. error“ column in Figure 3-76) is related to standard deviation  $\sigma$  as follows:

$$\sigma_e = \frac{1}{\sqrt{n}} \cdot \sigma$$

## File Browser

The File Browser, also called File Browser bar comprises several tabs. See 9 in Figure 3-9 and Figure 3-17.




Figure 3-17. File Browser

### Methods Tab

- Methods provide the complete description of a single measurement.
- Methods can be programmed or changed by the user.

Refer to “Creating a New Method” on page 3-20.

**Caution** Take the displayed methods only as a guideline. Do not use them for measurements! For measurements, always create your own methods! ▲

### Sequences Tab

- Sequences contain the description of a sequence of single measurements (methods).
- Sequences can be programmed or changed by the user.
- Different sequences have been predefined covering all basic measurements (in the **Examples** folder of the **Sequences** tab).

Refer to “Creating a New Sequence” on page 3-31.

**Caution** Take the displayed sequences only as a guideline, but do not use them for measurements! For measurements, always create your own sequences! ▲

**Caution** You must create and save a new method and a new sequence on your own!

The predefined methods and sequences delivered by Thermo Electron in the **Examples** folders are only example files. They only show guidance through helpful default values, but must never be used for measurements!

Never overwrite an example file with a method or sequence created on your own! Depending on your software version, these examples may not work properly. ▲

### **Export Tab**


- Edit voluminous amounts of Kiel IV Carbonate Device acquisition data for your own data systems using export templates (cf. LIMS).
- Refer to **Excel Export** in *ISODAT NT Operating Manual*; Part No. 109 2481.
- Use Isodat 2.5's Result Workshop to select and display particular aspects of your acquisition data.

Refer to **Result Workshop** in *ISOAT NT Operating Manual - Upgrade to Version 2.0*; Part No. 115 4991.

### **Results Tab**

- Provides access to all previously acquired measurement results.
- Gives an overview of all results.
- Is empty prior to the first measurement.

**Note** To easily transfer and store data at your place of choice (e.g. on a drive where data security is guaranteed), change the result path by a right-click and then select **Set Path**. The basic path is automatically installed. For reasons of data security, we recommend using this feature frequently. ▲




**Figure 3-18.** Modifying Result Path

From now on, all method, sequence and result files will be stored at a different location. See Figure 3-18.

**ISL Tab** Refer to ISODAT NT Operating Manual - Upgrade to Version 2.0; Part No. 115 4991.

**Note** All acquisition scripts are usually named **acquisition.isl**, no matter to which application they belong (e.g. GasBench II, ConFlo III or Kiel IV Carbonate Device).

However, they are stored in separate folders (e.g. in a particular Kiel IV Carbonate Device folder at C:\Thermo\Isodat NT\Global\ISL\Kiel IV Carbonate). ▲



**Figure 3-19.** Location of ISL Scripts for Kiel IV Carbonate Device

**Calibrations Tab**

- shows the mass calibrations for the IRMS

**Scans Tab**

- Instrument scans can be saved here.

**Search Tab**

- allows to find any result files of data acquisitions by pressing the **File Search** button.

- Like a file manager, it displays the results of a file search and allows to move files.

### Browser Tab

- If a Result Workshop document is open, this tab shows the objects that can be imported (e.g. methods, sequences, results).
- A file manager that allows browsing to an arbitrary directory of your choice, even to a root of a harddisk drive.
- As with other file managers, files and folders can be created, moved or deleted.

### Info Window


In the **Info** window, information during the process of data acquisition appears online. It displays the content of the info log files as well.

### Activating Info Window


Ensure that the **Info** window will indeed be displayed as follows:



Open **Acquisition** mode.



Press **Options** button.




Choose **Bars** tab.

Mark **Information** checkbox.


Confirm by **OK**.

After the **Info** window has been activated, it will appear beneath the **Accessories** bar. See Figure 3-20.



Right-click somewhere into the **Info** window.

**Figure 3-20.** Info Window



**Properties:** Log file activation and filter handling for **Info** window

**Print:** Print online information.

**Clear Buffer:** Clear buffer of **Info** window.

**Figure 3-21.** Commands for Using Info Window

### Properties Command



If you choose **Properties**, the window shown in Figure 3-22 will appear.




Figure 3-23

Figure 3-22. Info Properties Window

In general, all changes are effective globally, that is for all other Isodat 2.5 configurations.

**Note** During an acquisition, import and sharing of the log file are deactivated. A notation indicating a sharing violation appears. Importing a log file is only allowed during idle acquisition. Opening an ASCII log file in an appropriate program (e.g. MS Editor™) is only allowed, if Isodat 2.5 is entirely closed. ▲

### Enabling Log File




Figure 3-23. Enabling Log File\*


\* upper part of Figure 3-22

- A log file will be permanently saved at the assigned file location. Default is C:\Thermo\Global\bin\ Info Window Log Files\Workspace\Info.txt. See **1** in Figure 3-23.
- Choose the size of the log file at **3**.

**Note** The log file is stored in ASCII format. This allows to save huge amounts of information. Therefore, it is not necessary to activate the filter for the log file. See **4**. ▲

- Use the filter explained below for the online printout.

### Print Command



Use this command to print the log file.

Then, select your printer and the printer settings you prefer.

**Note** During an acquisition, move the mouse pointer onto the **Acquisition** bar in order to reactivate the acquisition. ▲

Now, you can click the **Print** button to print the total current "online information" directly while the acquisition is running.

### Offline View

Info	Detail
[02/10/06 13:37:48] 0 files saved.	
[02/10/06 13:37:48] Result Path	C:\Finnigan\Isodat NT\Global\User\Dual Inlet System\Results\ACQ-Results\
[02/10/06 13:37:51] Sample ISL Script 0	0
[02/10/06 13:37:52] ResetRatioPlotControls	
[02/10/06 13:37:53] Prepare Run - Start	
[02/10/06 13:37:53] Pump COV Sample Side 3 [Sec]	
[02/10/06 13:38:02] Set Trap 1 to 150 [°C]	
[02/10/06 13:38:02] Set Trap 2 to 150 [°C]	
[02/10/06 13:38:03] Set Trap 1 to 150 [°C]	
[02/10/06 13:38:54]	514.508
[02/10/06 13:38:54] <Action Thread Done>	
[02/10/06 13:38:54] The running Script 'Acquisition.isl' returns an Error - 'Script terminated by User'	The running Script 'Acquisition.isl' returns an Error - 'Script terminated by User'
[02/10/06 13:38:54] Method	
[02/10/06 13:38:54] Acquisition Warning I	
[02/10/06 13:38:56] Sequencer	Terminate Sequence I
[02/10/06 13:38:56] Sequencer	Sequence finished I

**Figure 3-24.** Offline View

## Creating a New Method

Isodat 2.5's **Acquisition** mode allows fully automated isotope ratio determination of carbon (CO<sub>2</sub>) and oxygen of carbonate samples. All parameters relevant for data acquisition of a sample are stored in a method.


**Note** For an extensive description of the options of the method definition, refer to ISODAT *NT Operating Manual* (Part No. 109 2481). This section describes only the entries that are specific for operating the Kiel IV Carbonate Device. ▲

**Caution** You must create and save a new method on your own! The predefined methods delivered by Thermo Electron in the **Examples** folder are only example files. They only show guidance through helpful default values, but must never be used for measurements! Never overwrite an example file with a method created on your own! Depending on your software version these examples may not work properly. ▲

The following steps are needed to create a new method.



1. Open **Acquisition** mode.
2. Select a configuration for Kiel IV Carbonate Device applications, e.g. **Kiel IV Carbonate**.
3. Select the appropriate Gas Configuration for the intended measurement type, e.g. CO<sub>2</sub>.
4. Press the **New** button.



5. Click on the **Method** icon.

6. Confirm by **OK**.

7. Proceed with “[Structure of Methods for Kiel IV Carbonate Device](#)” on page 3-21.

**Figure 3-25.** Creating a New Method

## Structure of Methods for Kiel IV Carbonate Device

All Kiel IV Carbonate Device methods are organized by the following tabs:

- **Instrument** tab
- **Peripherals** tab
- **Evaluation** tab
- **Printout** tab

**Note** In **Evaluation** tab and **Printout** tab, the currently active Gas Configuration is indicated. For example, Evaluation@CO<sub>2</sub> alludes to CO<sub>2</sub>, whereas Evaluation@N<sub>2</sub> alludes to N<sub>2</sub>. ▲

The following values are a guideline explaining the parameters for a CO<sub>2</sub> method using a Finnigan Kiel IV Carbonate Device and a Reference Refill.

As soon as a Kiel IV Carbonate Device configuration has been created, a predefined and stored method named **Kiel\_Carbo.met** exists.

### Instrument Tab

#### Experiment

Figure 3-26 shows the **Experiment** part of **Instrument** tab.

## Isodat 2.5

### Structure of Methods for Kiel IV Carbonate Device

**Figure 3-26.** Instrument Tab - Experiment

**Table 3-3.** Instrument Tab - Experiment

No.	Parameter	Description
a	Comment	Per default, this field is empty. You can type in comments about Method, Acquisition Script, Time Events, etc.
b	Gasconfiguration	Select the appropriate entry, e.g. CO2. Usually, the default entry can be accepted. Refer to the Status bar in <a href="#">"Accessories Bar and its Components"</a> on <a href="#">page 3-5</a> .
c	Acquisition Script	Select an appropriate acquisition script by a click on the  button. Acquisition.isl is the default entry and can usually be accepted. It controls the acquisition cycle.  To edit the acquisition script, press the  button.

**Caution** An acquisition script should only be edited by users trained on script editing, debugging and error tracking. Otherwise, potential errors within scripts, which are due to editing, may be hardly discovered afterwards. ▲

## Isotope MS

**Figure 3-27.** Instrument Tab - Isotope MS

**Table 3-4.** Instrument Tab - Isotope MS

No.	Parameter	Description
a	Integration Time [s]	time during which the current in the cups is integrated to form a data point triplet

## Peak Center

**Figure 3-28.** Instrument Tab - Peak Center

**Table 3-5.** Instrument Tab - Peak Center

No.	Parameter	Description
a	Predelay [s]	waiting time between activation of reference gas and start of peak center cycle
b	Postdelay [s]	waiting time between end of peak center cycle and start of data acquisition
c	Cup	Select the Peak Center cup, e.g. cup 3 as narrow center cup in a triple collector.

## Reference Refill

Some important parameters concerning the vacuum in the Reference Refill can be preset according to Figure 3-29 and [Table 3-6](#).

**Figure 3-29.** Instrument Tab - Reference Refill

**Table 3-6.** Instrument Tab - Reference Refill\*

No.	Parameter	Description
a	Pump Overlay Time	capillary pump out time of Reference Refill tank
b	Refill Time	gas flow time from Reference Refill tank into inlet port of standard bellow
c	FV Threshold	minimum pressure of standard bellow including valves and tubes evacuated with fore pump before pumping with turbo pump
d	HV Pump Time	pump time of bellows - including valves and tubes - with turbo pump

\*For further details, refer to "Reference Refill" on [page 5-22](#).

## Peripherals Tab

By means of the elements of **Peripherals** tab, you determine the properties of some peripheral options in the active method.

### Dual Inlet System

Here, you determine some properties for the Dual Inlet system.

## Isodat 2.5

### Structure of Methods for Kiel IV Carbonate Device

**Figure 3-30.** Peripherals Tab - Dual Inlet System

**Table 3-7.** Peripherals Tab - Dual Inlet System

No.	Parameter	Description
a	Reference	where the reference gas is available, that is in the left or right bellow
b	Sample	where the sample is available, that is in the left or right bellow
c	Number of Cycles	Measure e.g. 8 times sample and 8 times standard.
d	Idle Time [s]	waiting time after changing from sample to standard side and vice versa before integrating the ion intensities of m/z 44, m/z 45 and m/z 46
e	FV Threshold [mbar]	This value is used as a reference to the background vacuum reading of the Dual Inlet. The value will be checked prior to each measurement. Exceeding this value causes a fatal script error: a failure with pumps or leak exists. In this case, check the hardware!
f	HV Pump Time [s]	no function in this configuration
g	FV Pump Time [s]	no function in this configuration

## Background

**Figure 3-31.** Peripherals Tab - Background

**Table 3-8.** Peripherals Tab - Background

No.	Parameter	Description
a	Pre Delay [s]	waiting time (after the Changeover Valves 31 and 33 have been closed) until the real mass spectrometer background will be measured

**Table 3-8.** Peripherals Tab - Background, continued

No.	Parameter	Description
b	Integration Cycles	number of repetitions of the background determination The integration time selected above is used each time.

### Pressure Adjust

**Figure 3-32.** Peripherals Tab - Pressure Adjust**Table 3-9.** Peripherals Tab - Pressure Adjust

No.	Parameter	Description
a	On Cup	Always select the cup where m/z 44 is measured.
b	Delay Time [s]	waiting time after changing from sample to standard gas (or vice versa) before matching standard bellow to sample ion intensity Prior to each press adjust determination, this time will elapse as a predelay. It appears as "Equilibration Time" in the Info Window.
c	Tolerance [mV]	maximum acceptable ion intensity difference between sample and standard after matching (e.g. 50 mV meaning $\pm 25$ mV) If tolerance is exceeded, the press adjust will be repeated.
d	Master	The standard bellow at the right side must be adjusted to the level of sample ion intensity minus signal up $\pm$ tolerance. For carbonate applications, always take sample (e.g. if left = 4200 mV, match the right bellow to 4200 mV - 100 mV $\pm$ 25 mV).
e	Signal Up [%]	Match the ion intensity of standard gas less than e.g. 100 mV before closing V 25 and starting the Dual Inlet acquisition. After V 25 is closed, the signal will increase. The exponential relationship to the pressure in the tubing before the crimp and V 25 is approximately considered to be a linear one.

### Time Slicing




**Figure 3-33.** Peripherals Tab - Time Slicing

**Table 3-10.** Peripherals Tab - Time Slicing

No.	Parameter	Description
a	Slices	The single integration time specified in "Instrument" tab (see Figure 3-27) will be replaced by this number of individual integrations.
b	Slope Threshold	
c	Outlier Test	

### Carbonate Device



**Figure 3-34.** Peripherals Tab - Carbonate Device

**Table 3-11.** Peripherals Tab - Carbonate Device

No.	Parameter	Description
a	VM2 Leak Threshold [μbar]	After connecting the vial and pumping it (e.g. line 1, V 7 and V 13 open, V 12 closed), the pressure at the vacuum gauge VM2 must be below this threshold value. Otherwise, an error message will appear and preparation will start with the next sample.
b	VM1 Leak Threshold [μbar]	Maximum pressure rise accepted in any sample vial during leak test. If one of these vials shows a leak, the system skips this sample and starts the next preparation. Line 1 and line 2 are treated separately: here, the fourfold value is accepted. If one of these lines shows a leak, the system stops the whole sequence.
c	Acid Temperature [°C]	The oven temperature (controlled by the Jumo itron 16 temperature controller) must be set to this value. The temperature controller inside the heating cabinet (not the Pt 100 resistor of the Jumo itron 16) acts as the temperature set point!

**Table 3-11.** Peripherals Tab - Carbonate Device, continued

No.	Parameter	Description
d	Number of Acid Drops	Enter the desired number of acid drops per vial.
e	Reaction Time 1 [s]	Reaction time of carbonate. Time starts as soon as the last acid drop is injected.
f	Reaction Time 2 [s]	Pumpout time of non-condensable gases
g	Transfer Time [s]	CO <sub>2</sub> transfer time from trap 1 to trap 2
h	VM1 Expansion [μbar]	If the CO <sub>2</sub> gas pressure released from trap 1 exceeds this value, the gas will be systematically expanded and pumped until a pressure beneath this threshold is achieved.

**Note** Transfer time and CO<sub>2</sub> gas pumpout time are set at “[Process Timing Tab](#)” on [page 3-47](#). ▲

## Evaluation Tab


### Cycle

**Figure 3-35.** Evaluation Tab - Cycle**Table 3-12.** Evaluation Tab - Cycle

No.	Parameter	Description
a	Outlier Test	As outlier test, select “None” or “Sigma”. The “Sigma” outlier is varying criteria of rejection according the gaussian normal distribution. It simply rejects the values outside a barrier of 1 to n times the standard deviation.

**Note** The Kiel III Carbonate Device can be used only with Isodat 2.0. This feature is used to import Dual Inlet measurement files that have been recorded with Isodat 2.0. ▲

### Extended Parameters

**Figure 3-36.** Evaluation Tab - Extended Parameters

**Table 3-13.** Evaluation Tab - Extended Parameters

No.	Parameter	Description
a	Ion Correction Location	New ion corrections can be created by modifying predefined ones and saving them as new ones. They can be imported by Isodat 2.5. Use the folder C:\Thermo\Isodat NT\Global\NSL\Ion Correction for selecting the predefined and storing the new ion corrections.

### Evaluation Type




**Figure 3-37.** Evaluation Tab - Evaluation Type

**Table 3-14.** Evaluation Tab - Evaluation Type


No.	Parameter	Description
a	Evaluation Type	Select an appropriate ion correction for CO2 data evaluation from the list: "None", "CO2_SSH" (default) or "CO2_Craig".

Press the  button to add own scripts for ion corrections.

### Standard Parameter



**Figure 3-38.** Evaluation Tab - Standard Parameter (I)



**Figure 3-39.** Evaluation Tab - Standard Parameter (II)


**Table 3-15.** Evaluation Tab - Standard Parameter

No.	Parameter	Description
a	Ref. Name	<ol style="list-style-type: none"> <li>1. Select a Ref. Name from your standard database set in the Standard Editor, e.g. "Haus2" (Figure 3-38) or</li> <li>2. Edit the related <math>\delta</math> values (Figure 3-39). In this case, "User Defined" will be shown at "Ref. Name".</li> <li>3. New standards can be created in the Standard Editor.</li> </ol>

**Note** For correct reporting of sample  $\delta$  values, the data entered in this field must resemble the true isotopic values of your reference gas. ▲

### Printout Tab

In **Printout** tab, the use of printout templates is controlled. See Figure 3-40.

**Figure 3-40.** Printout Tab**Table 3-16.** Printout Tab

No.	Parameter	Description
a	Single	Selects a print template from the Result Workshop for an individual printout per sample.
b	Sequence	Selects a print template from the Result Workshop for a reduced printout per sample within a sequence summary.

### Saving a Method

After you met all your decisions throughout the tabs of the method, you must save it.


**Caution** You must create and save a new method and a new sequence on your own! The predefined methods and sequences delivered by Thermo Electron in the "Examples" folders are only example files. They only show guidance through helpful default values, but must never be used for measurements! ▲

**Caution** Never overwrite an example file with a method or sequence created on your own! Depending on your software version these examples may not work properly. ▲

To save a method, proceed as follows:

1. Do one of the following:

### Save Command




Click on the **Save** button to save a method previously created on your own.

### Save As Command



Click on the  arrow and choose **Save as...** to optionally choose a new name and folder for the currently active method (single document). See Figure 3-41.



**Figure 3-41.** Saving a Method

**Note** Notice that the particular folder is shown that contains the currently active method. See Figure 3-41. ▲

**Caution** Choose the folder above the **Example** folder, not the **Example** folder itself! This ensures not to mix or even overwrite the predefined example method with your own method. ▲

2. Give the method a significant name, e.g. similar to the sequence it corresponds to. Keep the extension .met.



3. Confirm by **Save**.

## Save All Command



Click on the  arrow and choose **Save All** to save all currently active Isodat 2.5 documents (e.g. methods, sequences, result files, Result Workshop files).

They will be stored without changing names and folders.

## Creating a New Sequence


After creating and saving a method (see “[Creating a New Method](#)” on page 3-20), a sequence must now be created as follows.

**Caution** As with methods, you must create and save a new sequence on your own! The predefined sequences delivered by Thermo Electron in the **Examples** folder are only example files. They only show guidance through helpful default values, but must never be used for measurements! ▲

**Caution** Never overwrite an example file with a sequence created on your own! Depending on your software version these examples may not work properly. ▲




1. Press the **New** button.



2. Click on the **Sequence** icon.

3. Confirm by **OK**.

**Figure 3-42.** Creating a New Sequence



4. Specify the number of samples, e.g. 10.
5. Confirm by **OK**.

**Figure 3-43.** Selecting Number of Samples

The sequence grid, Figure 3-44, contains all information about the individual samples bundled together in the sequence. See [Table 3-17](#).

**Note** Individual columns (e.g. **Line** or **Sample**) can be copied from example sequences. ▲

**Figure 3-44.** Sequence Grid

**Table 3-17.** Sequence Grid

Parameter	Icon	Description
Row		Each row refers to an individual sample.
Peak Center		<p>Marking it  allows performing a peak center procedure prior to measuring the particular sample. This ensures the peak to be in the middle of the cup.</p> <p>As this standard procedure is time-consuming, save a lot of time by omitting some peak centers. The device is sufficiently stable to operate during a certain time period without a peak center.</p>

**Table 3-17.** Sequence Grid, continued

Parameter	Icon	Description
Background		Marking it  allows performing a background scan.  COV will be closed. The instrument is idle for "Background Pre Delay". See a in Figure 3-31. Then, the current in the three cups will be recorded.
Press Adjust		Marking it  allows performing a press adjust prior to measurement.  As press adjust is essential, always mark "Press Adjust"!
Reference Refill		Marking it  allows performing a Reference Refill. Refer to "Reference Refill" on page 5-22. The frequency of the Reference Refill during a sequence depends on the sample amount to be analyzed, on the reference gas amount and on the fill level obtained during a given time period.
Line		Select line 1 or line 2 of the autosampler.
Sample		Select the position of the autosampler's turret.
Weight [mg]		optional, used for sample weight
Identifier 1, 2		optional, mostly used to identify the particular sample
Comment		optional, add an arbitrary comment concerning the particular sample.
Preparation		optional, add an arbitrary comment concerning sample preparation.
Method		important; the IRMS method edited at "Creating a New Method" on page 3-20 can be selected here from the pulldown list. By selecting it here, you determine the particular IRMS method to be used indeed during measurement. Without a selection from the pulldown list, no measurement will take place. Instead, the error message "No valid method found in sequence grid" will appear.



Fill Grid with Data

**Note** After you typed data in only one cell of the sequence grid, easily fill each of its columns: right-click the column and choose the **Fill Grid with Data** command. ▲


## Saving a Sequence

As done with a method (see "Saving a Method" on page 3-29), after defining the new sequence, you must save it before it will start. Proceed as follows:

**Caution** The predefined sequences in the **Examples** folder are only example files. They only show guidance through helpful default values, but must never be used for measurements! ▲


**Caution** Never overwrite a sequence example file with a sequence you created! Depending on your software version, these examples may not work properly. ▲

### Save Command




Click the **Save** button to save a sequence previously created on your own.

### Save As Command



Click on the  arrow and choose **Save as...** to optionally choose a new name and folder for the currently active sequence (single document).




**Figure 3-45.** Saving a Sequence

**Note** Notice that the particular folder is shown that contains the currently active sequence. See Figure 3-45. ▲


**Caution** Choose the folder above the **Example** folder, not the **Example** folder itself! This ensures not to mix or even overwrite the predefined example sequence with your own sequence. ▲

6. Give the sequence a significant name, e.g. similar to the method it corresponds to. Keep the extension .seq.



7. Confirm by **Save**.

### Save All Command



Click on the  arrow and choose **Save All** to save all currently active Isodat 2.5 documents (e.g. methods, sequences, result files, Result Workshop files). They are stored without changing names and folders.


### Starting a Sequence

To start the sequence, proceed as follows:




1. Press the **Start** button.

2. Define Parameters for Results Export, Printout and Sequence Scripts as follows. See Figure 3-46.



**Figure 3-46.** Defining Parameters for Handling Results

## Results



**Figure 3-47.** Defining Full Path for Results Storage

**Table 3-18.** Defining Full Path for Results Storage

No.	Parameter	Description
a	Folder Name	Define a folder for results storage. To choose another folder than the proposed one, browse via  .
b	File Name	Within the defined folder, define a file for results storage.

Use the **Pre** and **Post** buttons to automatically create meaningful folder names and file names (e.g. **ACQ-Results** as file folder and **Date** as a postfix will result in a folder **ACQ-Results\_060725**, if the acquisition was started the 25<sup>th</sup> of July 2006).

Folder and path for storage of single result data will be set at the **Results** tab of the File Browser. If no entry is made at **Folder Name** (a in Figure 3-47), the result data will be stored directly at the **Results** tab without creating a particular folder.

## Export



**Figure 3-48.** Defining Parameters for Results Export

**Table 3-19.** Defining Parameters for Results Export

No.	Parameter	Description
a	Format	Define the format of result data to be exported via the .wke template. Choose between "None", "Excel", "Lotus" and "ASCII". 
b	Export File Name	Name the export file.

**Note** Text strings exceeding 256 characters in one row will be truncated when exporting result data in Excel format (.xls). However, in case of export in Lotus format (.wk1) or ASCII format (.csv), no truncation of information will happen. ▲

### Printout




**Figure 3-49.** Defining Parameters for Results Printout

**Table 3-20.** Defining Parameters for Results Printout

No.	Parameter	Description
a	Yes/No	Decide, whether you want a printout (Yes ) or not (No).
b, c		If you want a printout, choose between:  one printout per sequence (b) or one printout per sample (c)

### Properties




**Figure 3-50.** Properties Box - Comment

**Table 3-21.** Properties Box - Comment

No.	Parameter	Description
a	Comment	Type an arbitrary comment applied to all result files in this sequence.
b	Measure only Selection	Mark, if only specific vials are to be measured, e.g. 1-7. The checkbox is active, if a selection in the sequence list has been made.  For Kiel IV Carbonate Device, only samples in connected lines (that is without an interception between them) can be measured by "Measure only Selection". Scattered runs however, are not possible and will lead to an error!

### Sequence Scripts



**Figure 3-51.** Selecting ISL Scripts to be Executed

**Table 3-22.** Selecting ISL Scripts to be Executed

No.	Parameter	Description
a	Pre Script	Select an ISL script (*.isl) to be executed before the sequence.
b	Post Script	Select an ISL script (*.isl) to be executed after the sequence. Useful for adding the "Standby & Drop" service script after a sequence

### Standby



**Figure 3-52.** Standby

**Table 3-23.** Standby

No.	Parameter	Description
a	Interface Standby after Sequence	Advanced users may add additional interfaces which can be put on idle by ISL script language.

3. Finally, confirm by **OK**. The measurement will be started.  
Refer to ["Measurement Procedures for Real Samples"](#) on page 5-1.

## Standards

Standards can be set in


- a. Isodat **Generic Editor** (recommended only for advanced users) or in
- b. Isodat **Standard Editor**.

## Using Generic Editor

### Determination of Primary Standards


1. Browse to the folder **C:\Thermo\Isodat NT\Global\bin**.
2. Open the Isodat **Generic Editor** by a double-click on the file **IsodatEditor.exe**.
3. Click on its **Open** button.





**Figure 3-53.** Selecting PrimaryStandards.std

4. Browse to the folder **C:\Thermo\Isodat NT\Global\Databases**.  
Select the file **PrimaryStandards.std**.  
Confirm by **Open**.  
The database containing the Primary Standards will be opened.



**Figure 3-54.** Opening Primary Standards Database

5. In the left pane, mark **PDB, Mineral**.  
The primary ratios will be shown.  
Check your primary ratios in the right pane.

## Using Standard Editor




In **Acquisition Mode**, click on the arrow right to the **Editors** icon.



Select **Standard Editor**.

The Standard Editor will be opened as displayed in Figure 3-55. It is a database that contains the defined **Working Standards**.



**Figure 3-55.** Opening Standard Editor

The Standard Editor contains a description of the referred isotopic standards. It allows to:



create new Working Standards that will be added to the database



add  $\delta$  values that will be added to an existing Working Standard



delete individual  $\delta$  values or entire Working Standards from the database

To create a new Working Standard:

**a** Type in a significant name for the standard.

**b** Type in the name of the gas or select it from the pulldown menu.

**Figure 3-56.** Creating New Standard

Name	Gas	Delta	Value	Primary Standards	Prim. Std. Ratio
		d 32O2/29N2	0.000	None	1.0
		d 40Ar/28N2	0.000	None	1.0
		d 32O2/40Ar	0.000	None	1.0
		d 44CO2/28N	0.000	None	1.0
		d 44CO2/40A	0.000	None	1.0
SO-SO2	SO,SO2				
		d 33S/32S	0.000	CDT	0.007878
		d 34S/32S	0.000	CDT	0.0450045
		d 17O/16O	0.000	VSMOW	0.0003799
		d 18O/16O	0.000	VSMOW	0.0020052
oliver	CO2				
		d 18O/16O	0.000	PDB, Mineral	0.002079
		d 13C/12C	0.000	VPDB	0.0111802
Haus 3	CO2		0.000	VSMOW	0.0020052
		d 13C/12C	0.000	VPDB	0.0111802


**Figure 3-57.** New Standard Appearing in Standard Editor

The new Working Standard has been included into the database and therefore appears in the Standard Editor. See Figure 3-57.

**Note** As based on Primary Standard Ratios, in the **Value** column of the new Working Standard, 0.000 is displayed. ▲

Primary Standards are based either upon calibrated international Standards (if available) or upon their relation to own Standards (if unavailable).

**Note** Primary Standard Ratios cannot be changed in Standard Editor! This is only possible in the database of the Primary Standards. ▲

A rectangular button with a thin blue border and the text "Save & Close" in a sans-serif font.

To save your changes, finally click **Save & Close**.

## Excel Export

Figure 3-58 shows a simple Excel Export template created for Kiel IV Carbonate Device using the Excel Export Editor. It can be used as an example for creating a customized export template.

**Note** In all cases, the following columns shown in Figure 3-58 should be exported. ▲

For more information about Excel Export, refer to *ISODAT NT Operating Manual - Upgrade to Version 2.0*; Part No. 115 4991.





Figure 3-58. Excel Export Template

## Dyn Externals

These hardware parameters are factory-set and are used in different parts of the process to perform carbonate analysis with the Finnigan Kiel IV Carbonate Device. Generally, it is possible to perform all kinds of carbonate analysis without changing these parameters. Since the unit can be connected to any Dual Inlet IRMS, we recommend to adjust the following parameters in order to optimize the individual application to be performed.

Two different groups of hardware parameters, **Dyn Externals** and **Service Scripts**, can be set in Kiel IV Carbonate window as follows.




- When **right**-clicking on a free spot within the device window, **Service Scripts** will be called. Refer to “[Service Scripts](#)” on [page 3-49](#)
- When **left**-clicking **directly** on the symbol, **Dyn Externals** will be called and can be edited. They are arranged in five tabs and will be discussed now.

**Figure 3-59.** Calling Dyn Externals vs. Service Scripts

### Basic Functionality Tab

**Basic Functionality** tab deals with variables that control the connect/disconnect procedure. See [Figure 3-60](#) and [Table 3-24](#).




**Figure 3-60.** Basic Functionality Tab

**Table 3-24.** Parameters of Basic Functionality Tab

Parameter	Comment
a	delay time for piston movement
b	pressure threshold for fore vacuum section
c	pressure threshold for auxiliary gas (used, when vials are disconnected)

**Pressure Fast Adjust Tab**     **Pressure Fast Adjust** tab deals with signals for bellow adjustment. See Figure 3-61 and Table 3-25. Refer to “Reference Refill” on page 5-22 as well.




**Figure 3-61.** Pressure Fast Adjust Tab

**Table 3-25.** Parameters of Pressure Fast Adjust Tab

Parameter	Comment
Smallest Useful Signal to Set [mV]	minimum ion intensity acceptable to perform measurement (default: 500) corresponds to the signal after Reference Refill with bellow at 100 % and depends on refill time
Largest Signal Achievable [mV]	maximum ion intensity achievable with the current amount of gas in the bellow Default: 6500 for IRMS (Delta <sup>plus</sup> , Delta <sup>plus</sup> XL, MAT 252); preamplifier dynamic range: 10 V. If the ion intensity of m/z 44 exceeds this value, Isodat 2.5 waits for the signal to drop before an attempt to adjust pressures is made. Default 36000 for IRMS (Delta <sup>plus</sup> XP, MAT 253); preamplifier dynamic range: 50 V. If the ion intensity of m/z 44 exceeds this value, Isodat 2.5 waits for the signal to drop before an attempt to adjust pressures is made.

**Temp Settings Tab**    **Temp Settings** tab is used to determine the temperatures for the traps. See Figure 3-62 and Table 3-26.




**Figure 3-62.** Temp Settings Tab

**Table 3-26.** Parameters of Temp Settings Tab

Parameter	Comment
Heatout Temperature [°C]	Temperature used during pumpout and cleaning cycles of the two traps (before and after each carbonate measurement in order to remove any impurities) Default: 150
CO <sub>2</sub> Freeze Temperature [°C]	must be set to -170 °C
Close V4 Temperature [°C]	initial pump out and cool temperature of trap 1. Default: -20 During leak test period, valves 12 (22), 1, 2, and 4 are open. If leak test is successful, Isodat 2.5 starts cooling trap 1 in order to trap the produced gases of the prepared carbonate. As soon as temperature of trap 1 reaches -20 °C, valve 4 is closed.
CO <sub>2</sub> Release Temperature [°C]	Default: -120 <b>Be extremely careful in changing this value!</b> If setting it to lower values, ensure that all gas is released during the equilibration time.
Start Transfer Temperature [°C]	initial temperature to transfer CO <sub>2</sub> gas from trap 1 to trap 2. As soon as trap temperature reaches -15 °C, valve 3 opens and transfer of CO <sub>2</sub> from trap 1 into trap 2 starts. Default: -15 It has no effect on isotope fractionation that the CO <sub>2</sub> expands into the crimped part of the capillary at -15 °C. At -120 °C, the CO <sub>2</sub> will be frozen into trap 2. The CO <sub>2</sub> loss through the crimping hole is negligible and conditioning of the capillary head is started.
Measurement Temperature [°C]	temperature of trap 2 just before starting an acquisition. Default: 30
Standby Temperature [°C]	at this temperature, traps 1 and 2 are set after finishing an acquisition. Default: 30

**Process Timing Tab**    **Process Timing** tab is used to set some time parameters for the traps. See Figure 3-63 and Table 3-27.




**Figure 3-63.** Process Timing Tab

**Table 3-27.** Parameters of Process Timing Tab

Parameter	Comment
Time allowed for pumping trap(s) [s]	pump time after setting the traps to pump temperature. Default: 120
Expansion Pump Time [s]	pump time during expansion cycles where part of the CO <sub>2</sub> is pumped away. Default: 180
Expansion Equilibration Delay [s]	waiting time after expansion of CO <sub>2</sub> gas before next action. Default: 180
Leak Test Rise Time [s]	will always be used to calculate the leak rate (in min <sup>-1</sup> ). Default: 30
Sample Side Pump Time [s]	pump time of the volume between valves 16, 15 and 32 of Dual Inlet system. Is used before each run. Default: 3

**About Acid Tab**    **About Acid** tab deals with parameters used in acid handling.




**Figure 3-64.** About Acid Tab

**Table 3-28.** Parameters of About Acid Tab

Parameter	Comment
Acid temperature tolerance [°C]	If oven temperature varies more than this value, acquisition stops. An error message will occur. Default: 1. Better is 3, because opening the door of the heating cabinet before the run is finished will lead to a fatal error. Using $\pm 3$ °C, this can be avoided.
Acid dropping timeout [s]	If no acid drop can be counted within this time (default: 120), the system stops indicating "fatal error". If at least one acid drop can be counted, the current sample will be measured. Afterwards, the system stops indicating a fatal error.
Pre process drop interval [s]	In order to avoid acid contamination due to rigid $H_3PO_4$ /carbonate reaction (phosphoric acid might well up into the steel body of the acid valve), a delay time after each dropping can be applied. Default: 10
Number of acid drops in drop test	Number of drops during an Acid Drop Test. Default: 5. See Figure 3-69.
Standby drop interval [s]	If the system is in Standby and Drop mode (e.g. after finishing or stopping an acquisition), every interval time (default: 3600), one acid drop will be injected into the waste vials (2/1 and 2/2). During the remaining time, the system stays in pump mode.
Standby drop: maximum number of drops	Default: 500. E.g. after 100 drops, the 2.5 ml vial is half-full.

## Service Scripts



- When **left-clicking directly** on the symbol, **Dyn Externals** will be called and can be edited. Refer to “[Dyn Externals](#)” on [page 3-44](#).
- When **right-clicking** on a free spot within the device window, **Service Scripts** will be called. They will now be discussed.

**Figure 3-65.** Calling Dyn Externals vs. Service Scripts




The following five service scripts will be outlined:

- Take Magazine
- Load Magazine
- Acid Drop Test
- Standby and Drop
- Standby and Pump

**Figure 3-66.** Available Service Scripts


## Take Magazine

Service Script **Take Magazine**



**Figure 3-67.** Service Script Take Magazine


To remove the magazine, proceed as follows:



1. Select **Take Magazine**.



2. If vials are connected, the adjoining information appears.
3. Confirm by **Yes**.




If, for some reason, the original position of the vials actually connected to the acid valve is not stored in the registry, a request appears: **Check the turret, and name an empty position for the vials to put.**

4. Confirm by **OK**.
5. The magazine will be moved to the respective position and vials will be taken away from the acid housing
6. After moving to position 1 automatically, the message **Ready to take magazine. OK** will appear.
7. Open the oven door and take the magazine out of the oven very carefully.
8. Close the door to maintain the oven temperature constant.

## Load Magazine

Service Script **Load Magazine**



**Figure 3-68.** Service Script Load Magazine


After inserting the magazine into the oven, this test should be performed to keep the lines free from water and other impurities. The following actions take place:

1. The left piston moves up vial 1/1 (first vial/line one).
2. Valve 7 opens and the vial is pumped. The computer waits until the pressure in the vial, measured by the pressure meter VM2, is below the acceptable value. See Figure 3-60. Then, the piston moves down.
3. The right piston moves up vial 2/1 (second vial/line one).
4. Valve 7 opens and the vial is pumped. The computer waits until the pressure in the vial, measured by the pressure meter VM2, is below the acceptable value. See Figure 3-60. Then, the piston moves down.
5. Leave the valves 7, 13 and 23 open.

**Note** If no vial is available or the associated proximity switch gives no response, an error message appears. ▲

## Acid Drop Test

Service Script **Acid Drop Test**



**Figure 3-69.** Service Script Acid Drop Test\*

\* See Table 3-28.

The Acid Drop Test checks the function of acid valves and counters. After connecting and measuring the VM2 value, valve 10 opens and injects one drop of acid to line 1. Then, valve 20 is opened, and one drop is injected to line 2.




**Warning** Do not perform the Acid Drop Test, if samples are loaded into the vial positions L1/P2 and L2/P2! ▲

This action is repeated until the requested number of drops is reached. See Table 3-28. To stop the action, select the properties again and click on **Stop Acid Drop Test**.

**Note** The vials 1/1 (vial one/line one) and 1/2 (vial one/line two) are so-called pump-position vials. This means that they are not sample prepared vials. They are used to keep the system clean and at vacuum. ▲

## Standby & Drop


### Service Script "Standby and Drop"




**Figure 3-70.** Service Script Standby and Drop

If the Kiel IV Carbonate Device is not used for more than 5 h, we recommend selecting the **Standby and Drop** function to keep the acid lines in flow. Otherwise, the acid may clog the lines or the magnetic valves 10 or 20 may squeeze the acid tubing. This function is placed as a postscript automatically after finishing an acquisition (see Figure 3-51). Standby and Drop can be performed manually as well.


**Note** Make sure that vials are available in position 1/1 (vial one/line one); 1/2 and 2/1 and 2/2 in the magazine prior to starting Standby and Drop. ▲



1. If vials are connected, the adjoining information appears.
2. Confirm by **Yes**.



3. The pistons take the vials. The magazine moves to position one and loads the vials of position 1/1 and 2/1, and the valves 7, 13 and 23 are opened.
4. The pumping period is started. After a time defined by the user (Standby Drop Interval, see Figure 3-6 and Table 3-2), vial 1/1 is taken away, and the magazine moves to position 2. Vial 2/1 (vial two/line one) is connected and pumped. After pumping, acid is dropped and vial 2/1 is removed into the magazine.
5. The same procedure takes place for vial 2/2 (vial 2/line 2).



To stop the function, at **ISL Scripts** select the script **StandbyAndDrop.isl**.

Then, click **Terminate Script**.  
This will terminate the executed script.

## Standby and Pump

Service Script **Standby and Pump** 1. Insert at least two clean vials in position 1/1 and 2/1.




2. The script sets Kiel IV Carbonate Device in Standby mode, that is, it closes valves 12 and 22 and opens valves 1, 2, 3, 4 and 5.

3. Leave valves 7, 13, 23 open.

4. Bring the vials of position 1 to respective acid valve housing and keep pumping.

**Figure 3-71.** Service Script Standby and Pump



**Figure 3-72.** Kiel IV Carbonate Device in Standby Mode

Figure 3-72 shows the Kiel IV Carbonate Device in **Standby** mode.

### Terminate

When **Terminate** appears within the above-mentioned service scripts, the process displayed in Figure 3-73 is executed.

Terminate




Figure 3-73. Actions during Terminate

## Interpreting Results

### Raw Tab

Raw <CO2>	Evaluated <CO2>		Grid - Errors	Grid - Infos	Sequence Line							
	44 Sample [mV]	45 Sample [mV]	46 Sample [mV]	44 Reference [mV]	45 Reference [mV]	46 Reference [mV]	rR 45CO2/44	rR 46CO2/44	R 45CO2/44	R 46CO2/44	rd 45CO2/44 [per mil] vs. CO2 I ah Tank	rd 46CO2/44 [per mil] vs. CO2 I ah Tank
Pre				1859.576	2229.242	2625.168						
1	1833.946	2211.448	2620.630	1806.475	2165.569	2550.127	1.205841	1.428957	0.012060	0.004007	5.885	12.238
2	1786.200	2153.937	2552.424	1755.019	2103.892	2477.417	1.205876	1.428969	0.012060	0.004007	5.917	12.277
3	1739.410	2097.538	2485.593	1704.928	2043.912	2406.787	1.205890	1.428988	0.012060	0.004008	5.910	12.287
4	1693.851	2042.652	2420.460	1656.660	1985.975	2338.518	1.205921	1.428968	0.012060	0.004008	5.934	12.283
5	1649.445	1989.025	2356.944	1609.554	1929.494	2271.962	1.205875	1.428931	0.012060	0.004008	5.917	12.299
6	1606.289	1936.954	2295.258	1564.098	1874.989	2207.808	1.205856	1.428919	0.012060	0.004008	5.910	12.305
7	1564.147	1886.160	2234.942	1520.109	1822.236	2145.665	1.205871	1.428858	0.012060	0.004007	5.932	12.270
8	1523.382	1836.941	2176.869	1477.230	1770.858	2085.122	1.205830	1.428970	0.012060	0.004008	5.897	12.367

Figure 3-74. Raw Tab - Part I

d 45CO2/44 [per mil] vs. VPDB VSMOW	d 46CO2/44 [per mil] vs. VPDB VSMOW	d 13C/12C [per mil] vs. VPDB	d 18O/16O [per mil] vs. VSMOW	d 17O/16O vs.	AT% 13C/12C [%]	AT% 18O/16O [%]
6.299	-2.927	6.862	-2.945	6.296	1.113161	0.199531
6.330	-2.889	6.894	-2.907	6.316	1.113196	0.199538
6.323	-2.879	6.886	-2.898	6.321	1.113188	0.199540
6.348	-2.883	6.913	-2.901	6.319	1.113216	0.199539
6.330	-2.867	6.893	-2.886	6.327	1.113195	0.199542
6.324	-2.861	6.886	-2.880	6.330	1.113187	0.199544
6.345	-2.895	6.910	-2.914	6.313	1.113214	0.199537
6.311	-2.800	6.870	-2.818	6.363	1.113169	0.199556

Figure 3-75. Raw Tab - Part II

## Individual Columns of Result File

Table 3-29 explains the individual columns of the result file shown in the **Raw** tab (Figure 3-74 and Figure 3-75).

**Table 3-29.** Explanations of Columns of Result File

Number	Column	Explanation
1, 2, 3	Sample [mV]	The signal readings during the given integration time (see Figure 3-27) are reported for all masses selected in the Gas Configuration and for sample and reference gas separately (columns 1-3 vs. columns 4-6). The values shown are corrected for background.
4, 5, 6	Reference [mV]	
7, 8	rR	The signal values are used to create amplifier-corrected raw molecular ratios for all selected ratio traces (defined in Ratio Editor for the Ratio Group used in the Gas Configuration).
11	rd45CO2/44CO2	Amplifier-corrected elemental ratios. The molecular $\delta$ value for the reference gas selected in the method (see Figure 3-38) is calculated from the elemental isotopic composition given in the same location. This is done by utilizing a reverse ion correction based on the algorithm selected in the "Evaluation Type" entry of the selected method (see Figure 3-37).
12	rd46CO2/44CO2	
13	d45CO2/44CO2	With the molecular ratios between sample and reference and the known molecular ratio of the reference gas, the molecular $\delta$ value of the sample can now be calculated. Refer to "Excel Export" on page 3-42.
14	d46CO2/44CO2	
15	d13C/12C	From the derived molecular $\delta$ values, elemental isotopic compositions are calculated using the ion correction selected (see Figure 3-37) in the method.
16	d18O/16O	
17	d17O/16O	

## Evaluated Tab

The **Evaluated** tab contains statistical information about the eight repetitions of a single measurement, e.g. their mean, standard deviation and standard error (that is, error of mean determination). Standard error depends on the number of repetitions whereas standard deviation does not. See columns in Figure 3-76.

The **Outlier** column reveals the number of outliers found during calculation of mean. The outliers found after running the outlier test will be highlighted red in the **Raw** tab. All outliers are then eliminated before a new, "corrected" mean, standard deviation and standard error will be calculated without them. These new, "corrected" values will finally be displayed in **Evaluated** tab.

Raw <CO2>	Evaluated <CO2>	Grid - Errors	Grid - Infos	Sequence Line
	Mean	Std. Deviation	Std. Error	Outlier
d 45CO2/44CO2	6.326	0.016	0.006	0
d 46CO2/44CO2	-2.875	0.036	0.013	0
d 13C/12C	6.889	0.018	0.006	0
d 18O/16O	-2.894	0.036	0.013	0
d 17O/16O	6.323	0.019	0.007	0

**Figure 3-76.** Evaluated Tab

$\delta_{\text{mean}}$  ("Mean" column in Figure 3-76) is given by:

$$\delta_{\text{mean}} = \frac{1}{n} \cdot \sum_n \delta$$

where n denotes the number of repetitions.

Standard deviation  $\sigma$  (“Std. deviation“ column in Figure 3-76) is given by:

$$\sigma = \frac{1}{\sqrt{n-1}} \cdot \sqrt{\sum_n (\delta - \delta_{\text{mean}})^2}$$


Standard error  $\sigma_e$  (“Std. error“ column in Figure 3-76) is related to standard deviation  $\sigma$  as follows:

$$\sigma_e = \frac{1}{\sqrt{n}} \cdot \sigma$$

## Grid Errors Tab

If errors occur during sample measurement, the related error messages will be shown in **Grid Errors** tab, Figure 3-77.

**Examples:** Acid cannot be dropped, acid has run out, sample vials cannot be found, sample vials are not tight, liquid nitrogen has run out.



**Figure 3-77.** Grid Errors Tab

## Grid Infos Tab

During the entire preparation process of the carbonate sample prior to measurement start, parameters are saved for informational purposes. They are summarized in **Grid Infos** tab. See Figure 3-78 and [Table 3-30](#).

After their determination, measurement will start. The results of the measurement will finally be shown in Figure 3-74 and Figure 3-75.

Raw <CO2>	Evaluated <CO2>	Grid - Errors	Grid - Infos	Sequence Line
Acid: 71.2 [°C]				
LeakRate [µBar/Min]: 142				
P: VM 1 Pre-Reaction : 45				
Total CO2 : 704				
Nr of Exp.: 0				
CO2 after Exp.: 706				
Peak Center found at [59483]				
Background: BGD 2006/Feb/24 - ,6.99 mV,4.58 mV,7.65 mV				
PressAdjust: L: 2294.8 R: 2218.5 ( Master Capillary )				

**Figure 3-78.** Grid Infos Tab

**Table 3-30.** Parameters of Grid Infos Tab

Parameter	Description
Acid [°C]	temperature of carbonate oven when acid is being dropped. Default: 72 °C
Leak Rate [µbar/min]	leak rate measured at sample vial after connecting it. Must be less than 400 µbar/min! Otherwise, sample preparation stops (no acid will be dropped. An error message will be displayed in “Grid Errors” tab). Maintenance: high leak rate indicates a) a leak in the line to V4 and V2 (both closed) b) density of H <sub>3</sub> PO <sub>4</sub> is beneath the range 1.92-1.95 g*cm <sup>-3</sup> due to water therein.
P VM1 Pre-Reaction [µbar]	pressure read off above trap 1 (VM1) after leak test, before acid will be dropped should be less than 100 µbar!
Total CO <sub>2</sub> [µbar]	overall amount of CO <sub>2</sub> , measured above trap 1 prior to possible expansion and transfer to trap 2 (that is, Microvolume). Correlates to sample size: small samples yield a smaller value than larger ones (Figure 5-9).
Nr. of Exp.	number of expansion cycles performed. Correlates to sample size (larger gas amounts yield more expansions than smaller ones). Should be small. Above 70 µg of sample, an expansion will be performed. The parameter “VM1 Expansion [µbar]” in the “Peripherals” tab of the method defines the VM1 pressure that must be reached to start an expansion (the system automatically dilutes the sample). See h in Figure 3-34 and in Table 3-11.
CO <sub>2</sub> after Exp. [µbar]	CO <sub>2</sub> pressure after expansion, measured at VM1. Should correlate with sample signal
Peak Center found at [DAC steps]	If you previously marked “Peak Center” in the sequence grid, a peak center will be performed before the gas is measured. See Figure 3-44 and Table 3-17. After the Peak Center procedure, this high voltage-dependent value (in DAC steps) reveals, where the peak center has been found.* It can be converted into a high voltage value and should be the same for all measurements to be performed. DELTA V: 1 DAC step equals 21.9 mV MAT 253: 1 DAC step equals 6.55 mV
Background [mV]	If you previously marked “Background” in the sequence grid, the result of the background measurement <sup>†</sup> will be displayed here. See Figure 3-44 and Table 3-17. If no background is measured, the time and date of the last background determined will be reported. The value should be small. Default: each of the three values (for m/z 44, m/z 45 and m/z 46) must be less than 10 mV. Commonly, background is stable and will be used only in the beginning of a sequence.
Press Adjust	As press adjust is essential, always mark “Press Adjust” in the sequence grid! The result of the press adjust will be displayed here. See Figure 3-44 and Table 3-17.

\* 65535 DAC steps at 10 kV (in case of a MAT 253) or at 3.3 kV (in case of a DELTA IRMS)

† That is, no gas will be let into the ion source, changeover valve is closed.

## Sequence Line Tab

In **Sequence Line** tab, the original content of the sequence line is displayed. See Figure 3-79.

Raw <CO2>		Evaluated <CO2>	Grid - Errors		Grid - Infos		Sequence Line						
Row	Peak Center	Background	Pressadjust	Reference Refill	Line	Sample	Weight [mg]	Identifier 1	Identifier 2	Analysis	Comment	Preparation	Method
32	1	0	1	0	2	17		17/2		620			Carbo Kiel IV.met

Figure 3-79. Sequence Line Tab

## Time Slicing

This feature has been introduced with Isodat 2.5. In Time Slicing mode, instead of a single, long integration, a series of shorter integrations will be performed. This is similar to run the IRMS in Continuous Flow mode. The number and duration of these integrations is determined by the number of time slices and the integration time selected in the method.

## Raw Sample Tab

**Raw Sample** tab, Figure 3-80, presents the results of the individual integrations for the sample side. This corresponds to columns 1-3 in Figure 3-74.




Figure 3-80. Raw Sample Tab

## Raw Reference Tab

**Raw Reference** tab, Figure 3-81, presents the results of the individual integrations for the reference side. This corresponds to columns 4-6 in Figure 3-74.




Figure 3-81. Raw Reference Tab

### Raw Complete Tab

**Raw Complete** tab presents the results of the individual integrations for sample and reference combined. This corresponds to columns 1-6 in Figure 3-74.

Figure 3-82 nicely shows the signal decrease of sample and reference gas due to viscous CO<sub>2</sub> sample gas flow in the two different sample and reference gas capillaries.





Figure 3-82. Raw Complete Tab

## Raw Ratios Sample Tab


**Raw Ratios Sample** tab, Figure 3-83, displays the ratio of m/z 45 to m/z 44 and the ratio of m/z 46 to m/z 44 for the sample only.

**Note** If a Dual Inlet valve failure appears, the signal will not decrease constantly. ▲



**Figure 3-83.** Raw Ratios Sample Tab

## Raw Ratios Reference Tab




**Figure 3-84.** Raw Ratios Reference Tab

**Raw Ratios Reference** tab, Figure 3-84, displays the ratio of m/z 45 to m/z 44 and the ratio of m/z 46 to m/z 44 for the reference only.


### Raw Ratios Complete Tab

**Raw Ratios Complete** tab, Figure 3-85, displays the ratio of m/z 45 to m/z 44 and the ratio of m/z 46 to m/z 44 for sample and reference. As failures become obvious here, it provides a tool for indication of a stable sample measurement.



**Figure 3-85.** Raw Ratios Complete Tab

### Ratios Tab



**Figure 3-86.** Ratios Tab

**Ratios** tab, Figure 3-86, corresponds to columns 7 and 8 in Figure 3-74. It contains the result after subsumption of all slices (the outliers have been mathematically eliminated).

## Interfering Masses

Interfering masses can be implemented in acquisition scripts as follows:

1. In an arbitrary acquisition script, include the following line in the section after the comment “place your includes here”:

Include "lib\InterferingMass\_lib.isl".

If you save a method that contains an acquisition script with this line and load it again, you will notice that a new tab has been added to the method tabs. If you click on this tab, the method page shown in Figure 3-87 will appear.

The screenshot shows a software window titled "Interfering Masses" with a tabbed interface. The active tab is "Interfering Masses". Below the title bar, there is a section titled "Interfering Masses" with the instruction "Type the requested mass or zero to skip". There is a checkbox for "Master Mass" with the note "(scales relative to first Mass)". Below this, there is a table with five rows, each representing a mass measurement. Each row has four columns: "Mass" (with a text input field containing "0"), "Channel" (with a dropdown menu showing "Cup 2"), and "Delay [s]" (with a text input field containing "0").

Mass	Channel	Delay [s]
0	Cup 2	0

**Figure 3-87.** Interfering Masses Tab

Here, you can select up to five masses to be determined during execution of the interfering masses function. The result of each measurement (mass, cup and intensity) will be reported within the **Extended** tab of the result file.

2. But first, a call to this function needs to be implemented in the acquisition file. To do this, fill in the following line at an arbitrary position within the “main” section of the acquisition file:

call DiMeasureInterferingMasses(2,2000)

The position of this 'function call' of course determines when the interfering masses will be measured during workflow. In the code example below, which is valid for the Kiel IV Carbonate Device, this point in time is immediately after the accomplishment of the sample measurement.

## Code Example

```
main()
{ call InitScript();
_LoadDynExternals("Carbonate Device\Dynexternals.iso");
if (_GetSequenceInfo(IS_FIRST_SAMPLE_RUNNING,TRUE))
{

    _RegSetProfileNumber("TASK_LOCKING", "PREP_READY",0);
    _RegSetProfileNumber("TASK_LOCKING", "T1_FREE",1);
    _RegSetProfileNumber("TASK_LOCKING",
    "HVPUMP_FREE",1);
    _RegSetProfileNumber("TASK_LOCKING",
    "DROP_WHEN_STANDBY",0);
    _RegSetProfileNumber("TASK_LOCKING",
    "ACID_IN_SAMPLE",0);
    _RegSetProfileNumber("TASK_LOCKING",
    "SEQUENCE_ABORT_FATAL",1);
    _RegSetProfileNumber("TASK_LOCKING",
    "FIRST_PREPARATION_DONE",0);
    call PrepareRun();

}

_StartFirstPreparation("\Carbonate Device\Instrument
Control\PrepareSample.isl");

string sSeqText=""; sSeqText=_GetSequenceText("Line","Empty");
number nLine=_strtod(sSeqText);
sSeqText=_GetSequenceText("Sample","Empty"); number
nSample=_strtod(sSeqText);

call MeasureOneSample(nLine,nSample);
call DiMeasureInterferingMasses(2,2000);

}
```

3. In some cases, it will be necessary to control the Changeover Valve to measure the interfering masses while the COV is closed or in a specific position. To do this, include one of the following lines in your code immediately before calling the interfering masses function:

```
call ChangeOverLeft();  
call ChangeOverRight();  
call ChangeOverClose();
```

4. Finally, it is possible to call the interfering masses function several times during an acquisition script. The set of masses and cups selected will remain the same, but all results are reported as expected in additional lines of the **Extended** section.

## Chapter 4 Basic Operations

This chapter contains the following topics:

- “Leak Check” on page 4-2
- “Bakeout of Kiel IV Carbonate Device” on page 4-6
- “Operating the Autosampler” on page 4-7
- “Matching Sample Capillary to Standard Capillary” on page 4-9
- “Cleaning Acid Valve” on page 4-14
- “Adjusting Liquid Nitrogen Refill Sensor” on page 4-20
- “Operating Pinch Valve” on page 4-21
- “Troubleshooting” on page 4-21
- “Elementary Handling of Kiel IV Carbonate Device” on page 4-22
- “Vial Test” on page 4-23
- “Phosphoric Acid Preparation” on page 4-26
- “Handling Sample Vials” on page 4-28

# Leak Check




Figure 4-1. Mass Spectrum of Background Gas Composition \*

\*for Finnigan DELTA V




Figure 4-2. Mass Spectrum of Background Gas Composition \*

\*for Finnigan MAT 253

To check whether the IRMS is ready to operate, close the inlet valve and run a mass scan from 3000 magnet steps to 12000 magnet steps. It should look more or less like Figure 4-1 or Figure 4-2, respectively.

The mass scan shows the composition of the background gas in the ion source region and informs about the amount of gases present. Try to identify the following patterns and compare them with the maximum values below (magnet step values for DELTA series):

**Water** Water in the background gas is indicated by the following properties:

- Contains ions of  $m/z$  16,  $m/z$  17 and  $m/z$  18.
- Appears at magnet current values approximately between 5300 steps and 6000 steps for DELTA series at 3 kV.
- Peak intensity after a couple of days of continuous operation should be max. 1 V for DELTA V and less than 500 mV for MAT 253.
- Intensity ratio of the three peaks is 1:2:4.
- CO<sub>2</sub> reference gas can be checked at normal sensitivity (e.g. 4 V gives a background level of  $m/z$  18 as the normal background).
- If the  $m/z$  18 signal of the CO<sub>2</sub> reference gas is too high, the CO<sub>2</sub> reference gas tank must be baked out and refilled again.

**Air** Air in the background gas is indicated by the following properties:

- Contains ions of  $m/z$  28,  $m/z$  32 and  $m/z$  40.
- Appears around magnet current values of approximately 7800 steps, 8500 steps and 9700 steps, respectively.
- Maximum intensity for  $m/z$  40 is 30 mV.
- Intensity ratio of the three peaks is 4:1:0.7.

**CO<sub>2</sub>** CO<sub>2</sub> in the background gas is indicated by the following properties:

- Contains ions of  $m/z$  28 (CO) and  $m/z$  44 (CO<sub>2</sub>).
- Appears around magnet current values of approximately 7800 steps and 10300 steps, respectively.
- Intensity of  $m/z$  44 must be less than 50 mV.

- The CO portion can easily be confused with nitrogen from air.

If air appears in the spectrum, check the IRMS for leaks, e.g. by using argon from a tank. In case of a too high water level, heat out the IRMS using the source heaters for at least 12 h. If a high water level is present in the ion source, usually some air is leaking into the IRMS as well.

Once this check has been performed within the given limits, open the inlet valve and repeat the mass scan. If air appears in the spectrum again, check all gas connections at the Kiel IV Carbonate Device for air leaks. Do not forget to check all connections under excess pressure as they may leak, too. The best way to find leaks in the excess pressure section is to use a standard soap solution (e.g. SNOOP®) which is applied to the connectors. Small bubbles appear when gas is leaking.

## Dual Inlet Ar Signal

The Dual Inlet argon signal should be very small, that is less than 20 mV on m/z 46 as the center cup ( $3 \cdot 10^9 \Omega$ ). When measuring small sample amounts (that is, about 10 µg), an argon background of less than 10 mV on the highest amplified cup is recommended. This is mandatory for all high vacuum lines.

## Advanced Leak Checking Procedure

An advanced leak checking procedure that can be used to search for the smallest leaks is described in this section.

1. Make sure that the IRMS uses the **CO2** configuration.
2. Select m/z 40 in middle cup where the amplification is set to  $3 \cdot 10^{10}$ .




Figure 4-3. Scan Parameters during Leak Check

3. Switch on the **Tune Scan** option in Instrument Control and start an acquisition in order not to miss events. Use a reasonable big acquisition time, e.g. 1200 s and zoom the time axis to display a window of 100 s. See Figure 4-3.

4. Make sure that you can look at the screen while performing the following leak test.
5. Use a gentle flow of argon, so small that it is hardly detectable when directed towards the lips. Spray all fittings and flanges using a plastic syringe tip as gas outlet.

**Note** During this entire procedure, **no** raise in the actual argon level should be visible! To precisely locate a possible leak, it is required to use an argon flow as small as possible. ▲

6. Directly touch all fitting surfaces and welded connections and move the gas outlet slowly over them.




Figure 4-4. Instrument Control as Aid when Leak Checking

**Note** Using high flows during this procedure may produce false signals due to argon moving around the IRMS in an unpredictable way. See Figure 4-4. ▲

Use Instrument Control as an aid when checking for leaks by monitoring m/z 40: while spraying argon to the vicinity of a leak, it will be detected by a sharp rise of argon intensity. See Figure 4-4.

## Bakeout of Kiel IV Carbonate Device

Upon installation and after some time, if precision deteriorates, the tubing of the Kiel IV Carbonate Device can be baked out in order to vacuum-clean the tubing. This procedure should be performed with the vacuum system switched on and the Kiel IV Carbonate Device in “Standby“ valve position. The following steps are required:

1. Heating the **capillary** to the IRMS by using the designated 7 V capillary heating transformer (see Figure 2-18 and especially step 19 on [page 2-14](#)):
  - a. Attach the middle steel point to 7 V.
  - b. Put the grounding to both capillary ends. At 7 V, temperatures above 180 °C are obtained.

**Note** Bake out capillaries only, if all other parts are baked out at the same time. ▲

2. Heating the **valve blocks** associated with traps 1 and 2. This can be performed by using a heat gun capable of heating to a maximum of 120 °C. The grease and O ring seals in the valve blocks are suited for this temperature. From now on, both valve blocks of trap 1 and trap 2 can be heated by a heating cartridge<sup>1</sup>.
3. Heating the **Micovolume** using an open flame. The tip of the Microvolume finger can be heated up to 600 °C using a suitable torch or a cigarette lighter.



**Warning** Be careful to heat only the finger and not the surrounding plastics or the valve block! ▲


---

<sup>1</sup>In former times, either a hair blower or heating bands (usable up to 450 °C) have been used instead of a heating cartridge. With these, it is easy to destroy the pneumatic part of the valve by uncontrolled heating. The valve contains an O ring seal and grease not suitable for temperatures above 150 °C!

4. Heating the **reference tank** before refilling with CO<sub>2</sub>.
  - a. Attach a 1/4" tubing to the Swagelok connector.  
See Figure 2-54.
  - b. Attach the other end of the 1/4" tubing to the internal standard inlet. If internal standard inlet is not available, valve ports of internal sample inlet or external right or left inlet ports can be used.
  - c. Bake out the reference tank using the heating bands over night, that is for 12 h.
  - d. Open the manual valve (for those manual parts which cannot be heated using heating bands).

## Operating the Autosampler

### Connecting Vials



**Figure 4-5.** Connecting Vials

The connection of a vial to the acid valve (see Figure 4-5) is a straightforward process, if no errors occur. Two ISL Dyn Externals variables control this process.

The first one is a delay that should be set according to the movement speed of your piston (see Figure 3-60). The second is a pressure threshold,  $p_x$  in Figure 3-60. It is used to mark an upper level to the vacuum quality that you want to accept in the vial section of the Kiel IV Carbonate Device.

Correct function and positioning of the proximity switch in the acid valve that is used to detect the presence of a vial is important for proper operation of the connect algorithm as well. Refer to “Proximity Switch” on page 2-32.

## Disconnecting Vials

Disconnecting Vials





Figure 4-6. Disconnecting Vials

Disconnecting a vial (see Figure 4-6) requires the fore vacuum section of the Kiel IV Carbonate Device to be leaktight and the vent gas to have sufficient pressure.

The first point requires that the connect algorithm works properly. For the second point, it is necessary for the vent gas pressure to be approximately 0.5 bar. In this case, the gauge VM2 shows a pressure of approximately 14 mbar. After venting again, the pressure is checked with V12 and V22 closed. This level can be adjusted via the ISL hardware parameters (see 5 in Figure 3-60).

## Matching Sample Capillary to Standard Capillary

This section describes how to match the sample flow from the Kiel IV Carbonate Device to the standard flow of the IRMS. As an example, see Figure 4-7. Flows must be matched to ensure equal voltage readings during the course of the measurement.



**Figure 4-7.** Matching Sample Capillary to Standard Capillary

The volume between V25, V26, V34 and V33 plus the capillary at the standard side of the IRMS has almost the same size as the volume between V3 and V5 plus Microvolume and the capillary of the Kiel IV Carbonate Device.

This means that matching the two signals from the two capillaries for an arbitrary pressure setting will result in an equal flow for both capillaries and, subsequently, for the parallel voltage drop for all possible signal heights.

The following procedure shows how to match the two capillaries:

1. Switch off the turbo pump of the Kiel IV Carbonate Device

**Note** In case of Kiel III Carbonate Device, the turbo pump can continue operating.

2. Close V1, V2, V3, V4 and V5 on the Kiel IV Carbonate Device.

## Basic Operations


### Matching Sample Capillary to Standard Capillary

**Note** In case of Kiel III Carbonate Device, V9 protects V4 from being vented. Therefore, only close V3, V5 and V9. Open V4.


3. Remove the flexible tubing which connects the valve combinations VC5 and VC3 (trap 1 and trap 2, respectively).
4. Connect an 1/8" stainless steel tube on the left sample port of the Dual Inlet (at V11). See arrow in Figure 4-8.



**Figure 4-8.** Connecting Transfer Tube to IRMS




**Figure 4-9.** Connection for Installation of Shortcut Link Tube



**Figure 4-10.** Connecting Transfer Tube to Kiel IV Carbonate Device

5. Connect the other side of the above tube to V5.  
See Figure 4-9 and Figure 4-10.
6. Now that the Microvolume (trap 2) of the Finnigan Kiel IV Carbonate Device is connected to the Dual Inlet, open both bellows to 100 %. Carefully pump all tubes including sample and reference of Dual Inlet.
7. Close V15 of sample side and let 50 mbar of CO<sub>2</sub> flow in at both sides of the bellows. V25, V24, V23, V13, V11, V14, V4, V2 and V3 are open.
8. Squeeze the carbonate capillary until an acceptable tolerance (that is, less than 50 mV) is reached. See Figure 6-7 and Figure 4-11.



**Figure 4-11.** Crimping Device at the End of a Capillary

## Basic Operations

### Matching Sample Capillary to Standard Capillary

9. Match the signal intensity of the carbonate capillary to the signal intensity of the reference side (usually the right side, see Figure 4-12). Frequently check the result by switching the Changeover Valve to the respective side. Allow for signal settling times in the range of minutes.
  
10. After reference matching, the lines must be baked out including all tubings and the entire inlet system (IRMS and Kiel IV Carbonate Device).

Figure 4-12 illustrates the complete plumbing with the link tube to perform capillary matching.





Figure 4-12. Vacuum Scheme of Kiel IV Carbonate Device

## Cleaning Acid Valve

In order to clean the acid valve, it must first be **dismantled**. After the **cleaning** procedure, it must, of course, be **reassembled**.

## Disassembling Acid Valve



**Figure 4-13.** Parts of Acid Valve \*

\*For the numbers and a more detailed description, refer to “Acid Valve” on page 6-8 and to Table 6-3.

**Note** At regular intervals, approximately every two or three months, the acid capillary might be clogged by crystallized acid. Disassembling and cleaning the acid valve avoids clogging. ▲

Always disassemble and clean the acid valve, if:

- a drop is indicated all along at trap 1 or trap 2 and therefore a “vial connect” failure appears either randomly or permanently.
- the leak rate is too high and therefore phosphoric acid gets inside the acid valve.

To disassemble the acid valve, proceed as follows:

1. Unmount the pinch valve by opening the two screws on each side. See Figure 4-14.



Figure 4-14. Acid Valve in Position



**Warning** Be sure to remove the valve together with the acid tubing in order not to spoil acid onto the pinch valve plunger! ▲

2. Remove the screw from the proximity sensor and pull out the sensor.
3. Remove the acid valve holder. See Figure 4-15.




Figure 4-15. Acid Valve - I\*

\*View from above with pinch valve mounted

4. Mark the three parts of the acid valve so that they can be tightened in the same position later on.

5. Loosen the four screws and remove the upper part with the acid capillary. See Figure 4-16.



**Figure 4-16.** Acid Valve - II\*

\*View from above without pinch valve mounted

**Note** Always renew the Teflon gasket, **12** in [Table 6-3](#). ▲


6. Remove the drop counter insert and electrical connection by opening the Swagelok connection at the electrical feedthrough. See Figure 4-17.



**Figure 4-17.** Drop Counter and Electrical Feedthrough

**Caution** Be careful not to destroy the acid dropper capillary and the drop counter spring during installation of the removable tray! ▲

7. Remove the spring plate from the bottom of the common mounting plate to get access to the screws that hold the acid valve in place. See Figure 4-18.



**Figure 4-18.** Mounted Spring Plate - Bottom View


8. Remove the bottom screws and the remaining Swagelok connection to the valve block in order to completely remove the acid valve. See Figure 4-19.



**Figure 4-19.** Mounted Valve - Bottom View

**Note** The Swagelok connection contains a removable stainless steel frit, Part No. 115 7670. See Figure 4-20. ▲

**Basic Operations**  
Cleaning Acid Valve




**Figure 4-20.** Valve Flange with Frit

All parts can be cleaned in deionized water and dried in acetone or in an oven at 70 °C.

**Note** Reassemble the acid valve in reverse order and thoroughly check for leaks! ▲

**Note** Use Apiezon H for sealing the Viton rings. See [www.apiezon.com](http://www.apiezon.com). ▲

### Brief and Superficial Cleaning Procedure



**Figure 4-21.** Cleaning Acid Valve

We recommend to clean the acid valve after every magazine that has been measured. Use a synthetic wiper.

**Note** The ends of the drop counter and the acid dropping capillary must be aligned. Otherwise, acid dropping stops (surface tension of the two wires or drops at the side wall of the vial). Dropping at the side wall will avoid acid reaction. ▲

To briefly and superficially clean the acid valve:

1. Use a lint-free wiper. Roll it up approximately as thick as a match.
2. Clean the inner part of the acid valve as shown at Fig. A in Figure 4-21.
3. Clean the Viton ring.

### **Thorough Cleaning Procedure**

**Note** During the cleaning procedure, be careful not to lose any of the Teflon rings of the drop counter! ▲

1. Put all parts into a beaker containing deionized water and heat it above 80 °C for about 30 min.
2. Clean the parts twice or more times.
3. Put the parts into acetone for fast drying.
4. Put the parts into the heating cabinet for drying.

### **Reassembling Acid Valve**

To reassemble the acid valve, first attach the acid capillary before installing the drop counter.

**Note** During reassembly, exchange the Teflon washer positioned beneath the first thick metal washer! ▲

### **Attaching Acid Capillary**

1. Lay the metal body onto a table.
2. Put the acid capillary into the metal body.

## Basic Operations

### Adjusting Liquid Nitrogen Refill Sensor

3. Put the two thick metal washers which fix the dropping capillary.
4. Tighten all screws which fix the two washers and the drop counter.

**Note** The hole of the proximity switch must be in-line! ▲

### Installing Drop Counter

1. The grounding (“shift“) must be drawn inwards through the Teflon body.
2. Consequently, the wedged tip can be tightened with pressure into the small hole of the drop counter.
3. Put the drop counter into the full metal acid valve body from below.

**Note** The counter tip must be opposite to the grounding tip! This will avoid a short circuit to the full metal body of the acid valve. ▲

## Adjusting Liquid Nitrogen Refill Sensor

If you encounter problems with setting temperatures below room temperature, the liquid nitrogen refill sensor must be adjusted. Before you actually touch the sensor itself:

1. Check that the liquid nitrogen tank is pressurized between 0.5 bar and 1.5 bar and that it contains liquid nitrogen. See Figure 2-47.

**Note** If the tank head pressure exceeds 1.5 bar, the consumption of liquid nitrogen during the liquid nitrogen refill process will increase dramatically. ▲

2. Check the electrical connection to the liquid nitrogen refill valve and the status of the manual valve in series with it.
3. Make sure that the resistor beneath the funnel can be heated and that no ice has been formed inside the capillary.

**Note** Frozen water will prevent the liquid nitrogen from rising into the cooling cascade! ▲

4. Only then, think of moving the sensor: move it carefully, only 1 mm at a time. If cooling is insufficient, move the sensor up. If heating is insufficient, move the sensor down.

## Operating Pinch Valve

The pinch valve operates at 24 V. Without an impressed voltage, it will therefore be closed. The connection operates with grounding. The pinch valve controls the acid flow and operates with viscous  $\text{H}_3\text{PO}_4$  up to 80 °C. The acid flow will be stopped at medium vacuum, that is approximately below  $7\text{-}9 \times 10^{-3}$  mbar.



**Warning** The supplied pinch valve is not acid-resistant. If any phosphoric acid is attached to the crimp of the pinch valve, the acid will degrade the pinch valve! The pinch valve must be immediately exchanged to avoid flooding the tubing of the Kiel IV Carbonate Device with acid! ▲

If no drop appears, the outer diameter of the Viton tubing is too wide. In this case, the tubing must be cut by a few centimeters. When the correct outer diameter is given, it will operate again. Open the acid valve using pliers. Now, phosphoric acid may flow through Viton tubing and acid valve.

Prior to a run a carbonate sequence, an acid drop test using carbonate-free position 2 test vials must be performed. Therefore, acid flow from the container through tubing and acid valve as well as vial connection and vial disconnection can be tested. If any of these tests fails, the underlying reason must be found and solved.

Acid flow can be stopped as a consequence of:

- improperly operating pinch valve
- crimped acid tubing because of too long idle time during closed pinch valve.  
Open the tube and avoid any phosphoric acid to drop on the crimp of the pinch valve. Cut or open the tubing using pliers in order to obtain acid flow. If no acid flows, the acid reservoir might be under vacuum or closed.
- a leak in the glass-to-tubing connectors.  
Close the straight connectors slowly and smoothly.

## Troubleshooting

**Note** From time to time, take a look at the checklist shown below. It outlines the performance achievable by the system. Check all mentioned items. ▲

1. A basic test must be performed, that is testing the IRMS alone.
2. Only if all previous items are performed to specifications, carry out the measurement.

## Basic Operations

### Elementary Handling of Kiel IV Carbonate Device

3. Finnigan Kiel IV Carbonate Device and IRMS must be free of leaks.
4. Bake out the transfer capillaries, traps and valves whenever possible to avoid contamination.
5. Check the drop counters for proper built-up of drops. Spraying of drops must be avoided!
6. Remove acid from O ring seals and drop counters using a lint-free cloth such as Kimwipes®.
7. Use new O ring seals on the acid valve.

## Elementary Handling of Kiel IV Carbonate Device




Figure 4-22. Kiel IV Carbonate Window

Figure 4-22 and Table 4-1 summarize basic operations to be executed via the **Kiel IV Carbonate** window.

It allows direct access to all Kiel IV Carbonate Device hardware components: autosampler turret control, trap temperature control, liquid nitrogen refill control, individual valves, vial connection and vial disconnection.

You can set or reset hardware components at any time, even during an acquisition. Click on a graphical object to operate the specific device.

**Table 4-1.** Basic Operations via Kiel IV Carbonate Window

No. in Figure 4-22	Comment
1	Move magazine to another position
2	Move piston up or down
3	Connect or disconnect vial. Use this button to check the connect algorithm and disconnect algorithm.
4	Set temperature of trap 1 or trap 2. Temperatures between -150 °C and 150 °C can be adjusted. Temperatures below -150 °C are not regulated.
5	Open or close valve. Refer to Figure 6-27 and Figure 6-42 for valve numbers.
6	Show available ISL scripts (via right-click somewhere)
7	Liquid nitrogen control. Press once to fill to 50 % fill level
8	Oven lamp (manual switch)
9	Oven control. Provides power to the Jumo itron 16 temperature controller. Refer to <a href="#">"Oven and Oven Control"</a> on <a href="#">page 2-27</a> .

## Vial Test

### Principle

The vial test monitors basic operations (connect/disconnect procedure on each position of the turret) and controls the usability of valves and vials for carbonate analysis.

A vial test method is a separate method without acid dropping and reduced timing. The leak rate will be reported for each vial. Thus, the vial test indicates, whether the individual vials ensure the needed quality.

We recommend performing a vial test when the Kiel IV Carbonate Device is set for first time of operation, after changing any hardware such as new pistons, magazine and after getting a new set of vials.

### Performing a Vial Test


The vial test is performed as follows:

1. Take a magazine and place a vial into each position.
2. Insert the magazine into the oven.


## Basic Operations

### Vial Test

3. Establish a vial test method as follows:



**Figure 4-23.** ISL Script is Loaded into Method




Browse to the folder **Kiel IV Carbonate** (see Figure 4-23).

Mark the ISL script **vialtest.isl**.

The ISL script will be loaded into the method (see Figure 4-23).

**Figure 4-24.** Selecting ISL Script for Vial Test



Save the method as a new method, e.g. as **Kiel IV vialtest.met**.

**Figure 4-25.** Saving Method as a New Vial Test Method

4. Use a new sequence to be run with the new method **Kiel IV vialtest.met** on each position (refer to “[Creating a New Sequence](#)” on [page 3-31](#)).  
During the connect/disconnect procedure, the following steps will be applied vial by vial:
  - a. Take a vial and connect it to the acid valve.
  - b. Measure the proximity switch, when the vial is connected.  
Refer to “[Proximity Switch](#)” on [page 2-32](#).
  - c. Open V7 and wait until the pressure at VM2 has fallen below 1000 mbar.
  - d. Bring the piston down and measure the pressure at VM2.  
Refer to “[Adjusting Piston Height](#)” on [page 2-32](#).
  - e. Open V13 and wait until the pressure at VM2 has fallen below 200 mbar.
  - f. When the vacuum is within an acceptable range (see **VM1 Leak Threshold** as **b** in [Figure 3-34](#) and [Table 3-11](#)), continue with the next vial.

### Possible Error Messages during Vial Test

During a vial test, the following error messages may occur:

- Vial is not connected.
- Vial is defective.

- Vial is missing.
- Connect procedure failed. In this case, test the connect/disconnect procedure and readjust the spring plate. See Figure 4-18.
- Proximity switch fails. See “Proximity Switch” on page 2-32.
- VM2 pressure is too high, that is a leak is present.
- Trap 1 and/or trap 2 are not operating properly. In this case, check the Autocool Unit (refer to “Autocool Unit” on page 2-24) and the liquid nitrogen dewar (refer to “LN2 Transfer from Refill Device into Dewar” on page 2-38).

## Phosphoric Acid Preparation

Phosphoric acid,  $\text{H}_3\text{PO}_4$ , is prepared from "Puranal" grade orthophosphoric acid ( $\geq 85\%$ ) and "Puriss" grade phosphorous pentoxide or trade names of equivalent purity. Inside a fume hood, one "Winchester" (that is, a 2.5 l package) of phosphoric acid is poured into a 5 l beaker positioned on a magnetic stirrer's hotplate. Use a magnetic stir bar (PTFE).



**Warning** Gloves and a face mask must be worn whenever handling phosphorous pentoxide! Goggles are not sufficient! ▲



**Warning** Between the additions and during the final cooling stage, the beaker is kept covered with cling film. ▲


**Note** A useful thermometer or stirring rod can be obtained by enclosing the thermometer in a large piece of heavy-walled Pyrex tubing with its bottom sealed off (that is shaped like a test tube). ▲

**Note** An alternative way of phosphoric acid preparation has been described by J. Burman, O. Gustafsson, M. Segl and B. Schmitz: A simplified method of preparing phosphoric acid for stable isotope analyses of carbonates. Rapid Commun. Mass Spectrom. 2005, **19**, 3086-3088. ▲

## Removing Water from Phosphoric Acid

Figure 4-26 shows the apparatus used to remove water and absorbed gases from phosphoric acid before phosphorous pentoxide is added. It avoids a rigid reaction with phosphorous pentoxide. Furthermore, the apparatus can be used to regularly dewater prepared 105 %  $\text{H}_3\text{PO}_4$ <sup>1</sup>.

**Note** Procedures and considerations about water content of phosphoric acid have been described by E.A. Wachter and J.M. Hayes: Exchange of oxygen isotopes in carbon dioxide–phosphoric acid systems. Chem. Geol. (Isotope Geoscience Section) **52** 365–374 (1985). ▲



- 1 500 ml water trap
- 2 acid
- 3 hotplate and magnetic stirrer
- 4 vacuum pump
- 5 dewar for liquid nitrogen

**Figure 4-26.** Removing Water from Phosphoric Acid

## Adding Phosphorous Pentoxide

It normally takes about 2 kg of phosphorous pentoxide to obtain the required final specific gravity of 1.92  $\text{g}/\text{cm}^3$  (or greater) at 25 °C. This quantity of phosphorous pentoxide is added gradually over a period of 2-3 h while constantly stirring and heating to a temperature of about 80 °C. The powder initially forms gelatinous lumps, but will gradually be dissolved. The complete process can take 4-5 h.

<sup>1</sup>The density of 105 %  $\text{H}_3\text{PO}_4$  is 1.921  $\text{g}/\text{cm}^3$  at 25 °C. Refer to [www.innophos.com](http://www.innophos.com).

If the phosphoric acid supplied by Thermo Electron (Bremen), Part No. 111 2640, has been dewatered before, addition of approximately 30 g of phosphorous pentoxide is sufficient.



**Warning** Take care during the initial stage of adding phosphorous pentoxide: the reaction can be vigorous as the powder contacts the relatively “wet” acid! ▲

The stirrer’s hotplate is switched off allowing the acid to cool down to room temperature before checking its specific gravity. If it is less than  $1.92 \text{ g/cm}^3$ , the acid must be reheated and more phosphorous pentoxide needs to be added. Finally, the acid which should be about 3 l after phosphorous pentoxide addition, is stored in bottles until required. Use Parafilm® to seal the screw cap.



**Figure 4-27.** Checking Specific Gravity of Phosphoric Acid


## Handling Sample Vials

### Manual Cleaning of Sample Vials

The sample vials used for carbonate measurements should be free of organic and inorganic contaminations before they are loaded with carbonate. To clean them, perform the steps described below:

1. Put the used vials into diluted phosphoric acid.

2. Repeatedly, rinse the vials with distilled water of high quality using a washing bottle (the quality of the distilled water is important).
3. Rinse the vials with acetone using a washing bottle, as well. This helps to dry the vials faster. Removal of dissolved inorganic carbon from the distilled water is also ensured.
4. Dry the vials in the Kiel IV Carbonate Device oven rack or in a drying chamber at 72 °C for 2.5 h. Cover them with aluminum foil to protect them against contamination. See Figure 4-28.



**Figure 4-28.** Oven Rack Containing Some Vials

## **Automatic Cleaning of Sample Vials**

Used sample vials can be cleaned automatically in a laboratory dishwasher made of stainless steel. The dishwasher must be connected to deionized water of high quality. Commonly, the procedure outlined below is performed:

1. The vials are cleaned using an alkaline detergent at 85 °C.
2. An acidic neutralization using citric acid or acetic acid is performed.
3. A complete automated rinsing with deionized water of high quality at 85 °C is performed.  
If the dishwasher is only connected to normal deionized water, perform an extra manual cleaning using deionized water of high quality.
4. The drying procedure can be enhanced by rinsing the vials in acetone or methylchlorine. This ensures removal of residual water that may contain acid-soluble minerals as well.

**Note** As a disadvantage, acetone always contains residues. Therefore, the vials should be dried upside down. ▲

5. Dry the vials in the Kiel IV Carbonate Device oven rack or in a drying chamber at 72 °C for 2.5 h. Cover them with aluminum foil to protect them against contamination. See Figure 4-28.


# Chapter 5 Measurement Procedures for Real Samples

This chapter treats the following topics:

- “Introduction” on page 5-2
- “Placing Sample into Vial” on page 5-4
- “Preparing Carbonate and IRMS” on page 5-6
- “Procedure” on page 5-7
- “Checking Quality of Result Data” on page 5-11
- “Referencing vs VPDB” on page 5-14
- “Reference Refill” on page 5-22

## Introduction

Before starting any sample preparation and performing measurements with the Finnigan Kiel IV Carbonate Device (Figure 5-1), it is important to get an overview of its operation procedure.



**Figure 5-1.** Finnigan Kiel IV Carbonate Device - Front View


It is assumed that the user is not only familiar with clean operating procedures and sample preparation, but also has working experience with IRMS, Microsoft Windows and Isodat 2.5.

Working with the Finnigan Kiel IV Carbonate Device means to operate a complex system running many processes at the same time. Therefore, to become an overview of the system, read this Operating Manual carefully prior to starting work. It is recommended to read the text below and procedure seriously as well as to compare the hardware layout to the system.

In this chapter, simultaneous measurement of  $^{13}\text{C}$  and  $^{18}\text{O}$  isotope ratios in calcite, aragonite (that is, mainly  $\text{CaCO}_3$ ) or dolomite (that is,  $\text{CaMg}(\text{CO}_3)_2$ ) will be covered. The latter is subject to a lot of discussion, and results should be treated carefully. The carbonate species reacts with phosphoric acid yielding  $\text{CO}_2$  that carries an image of the isotopic value of the carbonate ion  $\text{CO}_3^{2-}$ .

## Measurement Principle

Several methods and apparatus are available for the measurement of stable carbon and oxygen isotopes in carbonates. Isotope ratio determination of the above elements is usually carried out using  $\text{CO}_2$  produced in the reaction between carbonates and phosphoric acid. The common method employs a glass tube with a side arm. See Figure 5-2.



**Figure 5-2.** Common Glass Tube

The carbonate samples are loaded in the end of the glass tube. The side arm contains phosphoric acid. Reaction takes place after tilting the glass tube in order to bring acid to the sample. This method requires a stable temperature, e.g. 25 °C (or 50 °C to accelerate the reaction). The reaction time is between 2 h for calcite (sample amount above 5 mg) and 48 h for dolomite (sample amount above 10 mg).

After the reaction is finished, the glass tube contains  $\text{CO}_2$ ,  $\text{O}_2$ ,  $\text{N}_2$ ,  $\text{H}_2\text{O}$  and even  $\text{SO}_2$ , if the sample is not pure  $\text{CaCO}_3$ . The  $\text{CO}_2$  gas for the isotope ratio determination must be isolated and separated from other produced gases. This action takes place by trapping at different temperatures and pumping out  $\text{O}_2$  and  $\text{N}_2$  as non-condensable gases.

The Finnigan Kiel IV Carbonate Device is a fully automated preparation device for the precise and accurate determination of oxygen and carbon isotope ratios in carbonates, which are widely used in many fields of geology. One of the major applications is the estimation of ocean paleotemperatures from  $^{18}\text{O}/^{16}\text{O}$  ratios in marine microfossils.

Ultimate precision and accuracy are prerequisites for such work, because a temperature difference of 1 °C leads to a  $\delta^{18}\text{O}$  difference of about 0.2 ‰.

The Finnigan Kiel IV Carbonate Device has been developed in close cooperation with leading academic researchers with the first unit introduced in 1982. It incorporates a number of design changes which enhance ease of use, diminish the trace metal content of the acid and decrease both the costs of acquisition and operation, without changes to the underlying operation principles.

The system is designed for throughput at the highest level of precision and accuracy. Loading and exchanging the carousel are the only manual interactions with the system. Exchange of one autosampler carousel against a newly loaded one takes only minutes. Therefore, the system can run almost continuously in fully unattended operation. Each sample measurement takes about 15 min. Sample throughput of more than 10,000 samples per year is reported by a number of laboratories and can in fact be regarded as a routine.

## Placing Sample into Vial

### First Alternative

Most of the sample material used with the Kiel IV Carbonate Device is rather fragile in handling. Foraminifera and carbonate powders tend to charge electrostatically and stick to glass walls. Therefore, it is necessary to follow certain rules when loading this kind of samples into the sample vials (Part No. 075 4960). See Figure 5-3 and Figure 5-4.

We recommend using a small piece of weighing paper during weighing and in all filling steps. It can be folded to form a kind of boat.

**Note** The boat used for carbonate measurements should be clean and free of organic and inorganic contaminations. ▲

To place the sample into the vial, perform the following steps:


1. Carefully place the sample in the middle of the boat.  
See Figure 5-4 for correct sample placement.
2. Carefully lead the boat to the vial bottom in horizontal position.
3. Knock the vial in vertical position several times to ensure that all sample material is located at the bottom of the vial.

4. Take the boat out of the vial.
5. Use a new boat for each sample.

**Note** If weighing of the sample is required, weigh after step 1 and weigh the empty boat again after step 4. Often, parts of the sample material will not be released from the boat. ▲



**Figure 5-3.** Sample Vial - Part No. 075 4960



**Figure 5-4.** Correct Sample Placement

## Second Alternative

Another possibility to place a sample of known weight into the sample vial is widely used in paleoclimatological research.

You need a grained standard with a closely selected grain size with known average weight (e.g. NBS-19; grain weight here is 7  $\mu\text{g}$ ). Further, you need a small brush where you remove all hairs but six or eight, and a mirror or a glass plate on a black surface and good lighting.

Place a small number of grains on the glass plate and pick up one grain using the brush. Move the brush into a sample vial and gently knock the brush at the glass wall. The grain will fall into the vial and stay there. Repeat this until a suitable number of grains is placed within the vial.

Then, knock the vial on your desk to force the grains to the bottom. Place the vial in the turret and continue with the next sample.


## Preparing Carbonate and IRMS

1. Always use clean vials. Make sure that during sample preparation, no dust or other impurities fall into the vial.
2. Place the sample exactly at the bottom of the vial.
3. Do not leave the filled or clean vials outside the oven for a long time.
4. After placing the magazine within the oven make sure that the door is properly closed. Do not start the measurement immediately. Even if the oven shows the correct temperature, it takes at least 15 min until the vials are at the oven temperature.
5. Make sure that the oven temperature does not vary more than 0.5 °C within 15 min.
6. Make sure that clean vials without a sample are placed at position 1/1 and 2/1 (pump position).
7. In order to check the quality of your measurements or for calculation of your samples, place a few standard carbonates with known values in each 7<sup>th</sup> or 9<sup>th</sup> location.
8. If the Kiel IV Carbonate Device has not been used longer than 12 h, the first two measurements may not have the same precision as usual. Put more samples than necessary at this position.
9. Make sure that liquid nitrogen tank is properly filled and the manual valve (next to the magnet valve) is open. Take care that the liquid nitrogen filling tube is below the coil of the magnet valve.
10. Make sure that acid drop counter and O ring seal are clean. Refer to [“Brief and Superficial Cleaning Procedure”](#) on page 4-18.
11. Stay near the units at least during the first measurement and watch the measurement procedure.
12. Watch the liquid nitrogen refill unit during the first filling procedure of liquid nitrogen.

- Define the sequence parameters. Finally, start the acquisition.  
Refer to “Starting a Sequence” on page 3-35.

## Procedure

### Carbonate Process



**Figure 5-5.** Flowchart of Carbonate Process Including Timing

The Finnigan Kiel IV Carbonate Device consists of a temperature controlled oven, an acid tank, pneumatic valves, an autosampler containing the glass vials (magazine), gas cleaning facilities (trap 1) and a sample trapping arrangement (trap 2).

The reaction region is housed in a precision temperature-controlled oven and consists of an autosampler with 48 sample vials or thimbles. The thimbles can be loaded with as little as 10 µg of sample. The vials are made of special glass, allowing visual inspection of the sample and easy cleaning prior to re-use.

The reaction region also houses the reservoir of concentrated phosphoric acid, which is dispensed through two metal-free, acid-resistant dosing valves of new design.

The sample preparation magazine is a round tray with 48 holes arranged in two concentric rows - each with 24 holes. The inner row is indicated as line 1, and the outer row as line 2. A glass-made sample preparation vial can be placed into each hole.

The vials at position 1 are called pump vials and are indicated as vial 1/1 and vial 1/2 (this means vial 1/line 1 and vial 1/line 2). These vials are not used for sample measurement.

The vials at position 2 (2/1 and 2/2) have two functions. Initially, they are sample preparation vials and, after the samples of these vials are measured, they are secondly used as waste vials.

Once the sample containers are loaded and the sequence grid is filled out, all subsequent steps are fully automated. Prior to reaction and measurement, each container is evacuated. Now, a leak test is performed. The steps of a leak test are denoted in Figure 5-6.




Figure 5-6. Leak Test

Only if no leak is present, a precisely controlled amount of acid is added to the sample. The acid is exposed to metal surfaces at no time, with the exception of a short final steel capillary and a steel wire for droplet generation and counting. Thereby, contaminations due to the acid delivery system as well as cross-contaminations are minimized.

Thus, it is possible to chemically characterize the trace element content of carbonate samples by analyzing the spent acid amount, e.g. by ICP-OES or ICP mass spectrometry.

The CO<sub>2</sub> formed in the reaction of carbonate with acid is transferred into the trapping and gas cleaning system. It consists of a temperature-controlled trap with associated valves, an ultra high vacuum system, a pressure readout and a Microvolume.

With the temperature-controlled trap at temperature of liquid nitrogen, in a first step, CO<sub>2</sub> is quantitatively removed from the reaction region (along with some water) and frozen within trap 1. The CO<sub>2</sub> is then transferred into the Microvolume (that is trap 2) at -90 °C, whereas water is retained within trap 1. See Figure 5-7.

Trap Temperatures during Sample Preparation





Figure 5-7. Trap Temperatures during Sample Preparation

This CO<sub>2</sub> transfer is the basic principle of the Kiel IV Carbonate Device, because it ensures identical performance with varying sample sizes. Since most precise results require similar pressures in both sample and reference line of the Dual Inlet system, the total amount of CO<sub>2</sub> available is first detected by reading the CO<sub>2</sub> pressure at the trap when at -90 °C. A precision pressure transducer is therefore mounted at the position of the trap. The data system then decides which portion of the available CO<sub>2</sub> must be transferred into the Microvolume in order to

achieve the inlet pressure desired for measurement. If there is too much CO<sub>2</sub>, consecutive gas expansions between two volumes take place until the amount of CO<sub>2</sub> is just right to be transferred.

### Expansion Flowchart



**Figure 5-8.** Expansion Flowchart

This procedure assures maximum precision and is free of isotope fractionation, if the associated times are chosen above 2 min. See Figure 3-63.

It also eliminates any waiting times which would otherwise be necessary to allow the inlet pressure to drop to the desired level. Following this, the water is removed from the gas cleaning and trapping system by baking the trap and pumping all valves and gas lines.


Simultaneously, the Microvolume is warmed up, and the CO<sub>2</sub> flows from the Microvolume into the Changeover Valve of the Dual Inlet system. The Finnigan Kiel IV Carbonate Device is connected to the Changeover Valve via a dedicated capillary. This design minimizes the interaction of CO<sub>2</sub> with metal surfaces, a cause of memory effects and cross-contamination. It does not compromise any other inlet system that might be connected to the IRMS. Since all the valves are of all-metal gold-sealed design and all gas lines are thoroughly pumped, memory effects and contaminations are excluded.

**Note** A gas tank that contains a liquid phase like CO<sub>2</sub> requires absolute temperature stability! ▲

At this point in time, the Dual Inlet press adjust takes place. It intends to match the sample signal to the signal from the reference bellow. Refer to Figure 3-32 and Table 3-9 for details on press adjust. Immediately after the press adjust the Dual Inlet measurement starts. See Figure 3-13.

## Checking Quality of Result Data

The data present in the result document, and here especially the data collected in the **Grid Infos** tab, allow extended diagnostics on the collected data. Extract the information from the result document using an export template. Refer to **Excel Export** in *ISODAT NT Operating Manual*; Part No. 109 2481.



**Figure 5-9.** Pressure vs. Mass

Figure 5-9 shows a typical response curve for VM1. On the x-axis, the sample weight is denoted. In this special case, single grains of NBS-19 were counted and not actually weighed. This explains the pattern.

On the y-axis, the pressure displayed by VM1 for the respective sample is shown. The point close to 0/0 is faked to yield a nice regression curve for this example.

Use this plot type to check the completeness of reaction between carbonate sample and acid. This plot allows to judge the efficiency of trapping CO<sub>2</sub> in trap 1 as well. If any of the above parts of the preparation process fails, the gas amount (denoted by the VM1 pressure) will be smaller than expected for the sample weight used.

If, on the other hand, a leak is present, the displayed pressure at VM1 will be higher than expected. This allows to precisely check the dynamic behavior of the valves (V1, V2, V4, V12 or V22).




Figure 5-10. Signal vs. Pressure

Figure 5-10 explains the relation between the amount of gas measured in VM1 after expansion vs. the signal for m/z 44 achieved in the IRMS.


Use a plot like this to control complete gas transfer from trap 1 to the IRMS source. All points must be located along a quadratic correlation curve. Otherwise, severe problems in gas transfer are present.

Several reasons may contribute to gas loss in the section between valve V2 and the IRMS. Valve V5 is located closely to the liquid nitrogen dewar. During refill cycles, it may cool down and get stuck. In this case, you notice complete gas loss for a given sample.

Another possible reason for a virtual gas loss is the Press Adjust process. If the amount of reference gas available is too small to adjust for a given beam intensity, the algorithm may fail and decide to wait for the sample signal to drop to the reference level. Be sure to set a suitable value for the “largest signal achievable“ in the external parameters for the Kiel IV Carbonate Device. See Figure 3-61.

Figure 5-11 explains the relation between  $\delta$  value and signal. Ideally, for both isotopes inspected, the correlation curve should be a line with no slope. The scatter of the data points around this curve should always be within the Standard Gaussian distribution.

Unfortunately, in reality the scatter increases below a certain sample size. For details, see the advertised specifications. The real correlation curve exhibits some general behavior:



**Figure 5-11.**  $\delta$  vs. Signal\*

\*Both the range of erroneous  $\delta$  values and the linearity region are general properties, that is they occur in every carbonate measurement ( $^{13}\text{C}$  and  $^{18}\text{O}$ ).

Below a certain voltage, usually around 1.5 to 2 V, the  $\delta$  value drops considerably below the real value. This lower limit of the lower range depends largely on the Press Adjust and Reference Refill parameters.

When using Reference Refill, the bellow is completely open (that is 100 %) and gas will be filled in during a certain time period. After the gas has been filled in, that is Reference Refill is finished, a certain pressure exists in the bellow corresponding to a certain signal height.

In case much gas was filled in, no smaller signal heights (that is smaller reference sample amounts) can be adjusted via the Press Adjust algorithm as the bellow is already completely open. Therefore, any smaller sample amounts will not have a correct Press Adjust.

The measurement will be carried out in spite of the incorrect Press Adjust and must consequently yield erroneous  $\delta$  values. These are displayed in the left part of Figure 5-11.

This can be counterbalanced by filling a smaller gas amount into the bellow (that is less than 100 %) via a changed Reference Refill. Press Adjust is thus corrected via Reference Refill.

The linearity region (see right part of Figure 5-11) will always exhibit a marginal positive or negative slope about 0.01-0.02 ‰/V or even smaller. In case the standard deviation within the linearity range is extremely small, this linearity error can be counterbalanced as well. This will add further precision and accuracy.

## Referencing vs VPDB

All carbonate  $\delta$  values must be referenced to the international standard VPDB (Vienna Pee Dee Belemnite), the successor of PDB as PDB is exhausted. However, VPDB with  $\delta^{13}\text{C} = 0$  and  $\delta^{18}\text{O} = 0$  as one would expect, does not exist. Instead, standards exist which are related to this virtual, that is unreal definition. See [Table 6-26 on page 6-17](#).

**Note** See “Reference and intercomparison materials for stable isotopes of light elements”. In: IAEA-TECDOC-825, IAEA, ed., Vienna, 1995. ▲

At present, there are a couple of primary standards available from IAEA and NIST, respectively, with given  $\delta$  values for  $^{18}\text{O}$  and  $^{13}\text{C}$ . To determine the actual  $\delta$  value of a sample relative to VPDB, measure standard and sample under the same conditions and perform the following procedure:

- Determine the  $\delta$  value of your working standard.
- Calibrate versus known standards supplied by IAEA or NBS.
- Use a primary standard to determine the  $\delta$  value of the reference gas.
- With x meaning working standard and z denoting VPDB, the following equation is valid. Refer to “[Remark on the Strange Mathematics of Delta Values](#)” on [page 5-16](#):


$$\delta_z^x = \frac{\delta_y^x \cdot \delta_z^y}{1000} + \delta_y^x + \delta_z^y$$

and

$$\delta_z^y \neq \delta_y^z$$

with:

- x working standard
- y gas
- z absolute standard (that is VPDB)



**Figure 5-12.** Calculation Example\*

\*SHK designates Solnhofen limestone.

Figure 5-12 depicts an example for obtaining  $\delta$  values specified against VPDB starting from measured and corrected  $\delta$  values.

1. Determine absolute  $\delta$  value of the primary standard.  
In this example:

$$\delta_{\text{PDB}}^{\text{SHK}} = 1.750$$

2. Invert the measured value for primary standard versus gas used:

$$\delta_{\text{Gas}}^{\text{SHK}} = 1.908$$

Thus:

$$\delta_{\text{SHK}}^{\text{Gas}} = -1.903$$

3. Determine the absolute  $\delta$  value of the gas used today with the aid of the equation:

$$\delta_z^x = \frac{\delta_y^x \cdot \delta_z^y}{1000} + \delta_y^x + \delta_z^y$$

Thus:

$$\delta_{\text{PDB}}^{\text{Gas}} = -0.157$$

4. Use this value and any measured sample  $\delta$  vs. reference gas to calculate the  $\delta$  value of sample vs. PDB with the aid of:

$$\delta_z^x = \frac{\delta_y^x \cdot \delta_z^y}{1000} + \delta_y^x + \delta_z^y$$

Thus:

$$\delta_{\text{PDB}}^{\text{NBS 19}} = 1.767$$

In this case, the result is incorrect.

### Remark on the Strange Mathematics of Delta Values

The  $\delta$  definition:

$$\delta_y^x = \left( \frac{R_x}{R_y} - 1 \right) \cdot 1000$$

with:

$\delta_y^x$   $\delta$  value of x against y  
 $R_x$  raw ratio of x (that is  $A_{13}/A_{12}$ )

can be rearranged:

$$\frac{R_x}{R_y} = \frac{\delta_y^x}{1000} + 1$$

As x and y are only arbitrary notations and thus can be interchanged, an analogous equation for  $\delta_x^y$  can be written:

$$\frac{R_y}{R_x} = \frac{\delta_x^y}{1000} + 1$$

Considering reciprocity:

$$\frac{R_y}{R_x} = \frac{1}{(R_x/R_y)}$$

Combination of both equations yields the relationship between  $\delta_y^x$  and  $\delta_x^y$  we were aiming at:

$$\frac{\delta_x^y}{1000} = \frac{1}{\frac{\delta_y^x}{1000} + 1} - 1$$

This shows indeed:

$$\delta_x^y \neq \delta_y^x$$

The  $\delta$  definition results in the following rule when calculating a  $\delta$  value with an intermediate result, which is always the case when referencing to a gas or a working standard:

$$\delta_z^x = \left(\frac{R_x}{R_z} - 1\right) \cdot 1000 = \left(\frac{R_y \cdot R_x}{R_y \cdot R_z} - 1 - \frac{R_y}{R_z} + \frac{R_y}{R_z} - \frac{R_x}{R_y} + \frac{R_x}{R_y} - 1 + 1\right) \cdot 1000$$

$$\delta_z^x = \left(\frac{R_y \cdot R_x}{R_y \cdot R_z} - \frac{R_y}{R_z} - \frac{R_x}{R_y} + 1 + \frac{R_y}{R_z} + \frac{R_x}{R_y} - 1 - 1\right) \cdot 1000$$

$$\delta_z^x = \left(\left(\frac{R_x}{R_y} - 1\right) \cdot \left(\frac{R_y}{R_z} - 1\right) + \left(\frac{R_x}{R_z} - 1\right) + \left(\frac{R_y}{R_z} - 1\right)\right) \cdot 1000$$

$$\delta_z^x = \frac{\delta_y^x \cdot \delta_z^y}{1000} + \delta_y^x + \delta_z^y$$

This equation has been used above (special case: working standard x, absolute standard z, that is VPDB). See [“Referencing vs VPDB” on page 5-14](#).

## Ion Correction

The isotopic composition of a compound A is expressed by its  $\delta$  value,  $\delta_A$ ,

$$\delta_A = 10^3 \cdot \left(\frac{R_A}{R_{ST}} - 1\right)$$

Here,  $\delta_A$  is given in ‰.  $R_A$  denotes the isotope ratio of compound A and  $R_{ST}$  the defined isotope ratio of a standard sample.

## Examples 1. Carbon<sup>1</sup>

<sup>1</sup>As usual, the index SA refers to sample and the index ST to standard.

$$\delta^{13}\text{C} = 10^3 \cdot \frac{(^{13}\text{C}/^{12}\text{C})_{\text{SA}} - (^{13}\text{C}/^{12}\text{C})_{\text{ST}}}{(^{13}\text{C}/^{12}\text{C})_{\text{ST}}}$$

## 2. Oxygen (measured as CO<sub>2</sub>)<sup>1</sup>

$$\delta^{18}\text{O} = 10^3 \cdot \frac{(^{18}\text{O}/^{16}\text{O})_{\text{SA}} - (^{18}\text{O}/^{16}\text{O})_{\text{ST}}}{(^{18}\text{O}/^{16}\text{O})_{\text{ST}}}$$

For CO<sub>2</sub>, the IRMS measures the 45/44 and 46/44 ratios. M/z 45 comprises <sup>13</sup>C<sup>16</sup>O<sub>2</sub> as well as <sup>12</sup>C<sup>17</sup>O<sup>16</sup>O so that the 45/44 ratio is different from the ratio <sup>13</sup>C/<sup>12</sup>C by a correction regarding the ratio <sup>17</sup>O/<sup>16</sup>O in the sample or standard. Therefore,

$$R_{45} = \frac{^{13}\text{C}^{16}\text{O}_2 + ^{12}\text{C}^{17}\text{O}^{16}\text{O}}{^{12}\text{C}^{16}\text{O}_2} = \frac{^{13}\text{C}}{^{12}\text{C}} + 2 \cdot \frac{^{17}\text{O}}{^{16}\text{O}} = R_{13} + 2 \cdot R_{17}$$

In this equation, the definitions

$$R_{45} = \frac{\text{signal (m/z 45)}}{\text{signal (m/z 44)}}$$

and

$$R_{13} = \frac{\text{isotopic abundance (m/z 13)}}{\text{isotopic abundance (m/z 12)}}$$

(similarly for R<sub>17</sub> etc.) are used.

Ion correction routines must be applied to the measured ratios in order to account for the additional ion species contributing to the measured ratios. Likewise, other ion species must be subtracted from the 46/44 ratio for oxygen in CO<sub>2</sub>, the 65/64 and 66/64 ratios for sulfur, and the 34/32 ratio for elemental oxygen<sup>1</sup>.

The <sup>17</sup>O correction applied in Isodat 2.5 follows the suggestions given by Santrock and coworkers<sup>2</sup> with:

<sup>1</sup>Refer to H. Craig: Isotopic standards for mass spectrometric analysis of carbon dioxide. *Geochimica Cosmochimica Acta*, **12** (1957) 113-140.

$$R_{17} = K \cdot R_{18}^a$$

using  $a = 0.516$  and  $K = 0.0099235$ .

The  $\delta$  values of the working standard against an international standard must be known.

## Neogloboquadrina Pachyderma (Ehrenberg, 1894)

Neogloboquadrina Pachyderma is the most abundant planktonic foraminifer of high latitudes. As with any planktonic foraminifer, it avoids low-salinity and shallow waters. The left-coiled morphotype prevails at lowest temperatures and occurs throughout the Arctic Ocean.




Figure 5-13. Neogloboquadrina Pachyderma

## Events during Sample Measurement

To illustrate the primary steps, Figure 5-5 shows a flowchart of the events that belong to the measurement of one sample. The running software relies on multithreading to save precious time during a full run. Thus, preparation and measurement are separate functions that are interlocked by registry variables.

When following the sequence on screen, it is sometimes difficult to judge which step exactly takes place at this moment. Therefore, it is helpful that you familiarize with the output presented in the **Info** window. See Figure 5-14. Here, various information is presented in order to judge on the evolution of the measurement process. Also, if error messages appear, the previous informative “user infos” can be used to judge on the origin of this specific error.

The sample printout shows the start of a sequence until the first sample preparation starts. Reported are beginning and end of subroutines (such as **Prepare Run**), orders to set trap temperatures and turret movements as well as vial connections and vial disconnections. See Figure 5-14.

<sup>2</sup>Refer to J. Santrock, S.A. Studley and J.M. Hayes: Isotopic analyses based on the mass spectrum of carbon dioxide. *Analytical Chemistry*, 57 (1985) 1444-1448.




Figure 5-14. Evolution of Measurement Process in Info Window

**Example** Measurement started from vial 2/1 (vial 2/line 1)

VM1: vacuum gauge of trap region


VM2: vacuum gauge of rotary pump

1. As soon as the user starts the acquisition, the Kiel IV Carbonate Device will be initialized.
  - a. Vial 1/1 (vial 1/line 1) and vial 1/2 (vial 1/line 2) are connected and pumped.
  - b. Trap 1 and trap 2 are heated out.
2. Oven temperature and stability are checked.
3. Vial 1/1 is removed and magazine rotates to position 2.
4. Vial 2/1 is connected to acid valve housing.
5. Rotary pump evacuates vial 2/1 and the corresponding lines.
6. Vial 2/1 is evacuated to high vacuum, if no leak is present.
7. A leak check is performed. A leak rate is reported in the **Info** window and the result file.
8. Trap 1 is cooled down to  $-196\text{ }^{\circ}\text{C}$  by liquid nitrogen.
9. The reaction between acid and sample takes place. The produced gases  $\text{CO}_2$  and  $\text{H}_2\text{O}$  are trapped by trap 1.
10. Non-condensable gases are pumped out of trap 1.
11. Trap 1 is heated and  $\text{CO}_2$  will be released, whereas  $\text{H}_2\text{O}$  is still frozen.
12. The  $\text{CO}_2$  pressure is measured via VM1. If the pressure is too high,  $\text{CO}_2$  will be pumped out until an acceptable pressure is achieved.
13. Trap 2 is cooled down to  $-196\text{ }^{\circ}\text{C}$  by liquid nitrogen.
14.  $\text{CO}_2$  is transferred from trap 1 to trap 2.

15. Trap 1 is heated and H<sub>2</sub>O is pumped out. In the meantime, trap 2 is heated and CO<sub>2</sub> is released to the IRMS.
16. Vial 2/1 is removed and the magazine moves to the pump position.
17. A peak center is performed.
18. The reference gas pressure is high-end adjusted to the same value as the sample pressure.
19. Data acquisition starts, and the results are stored and/or printed.

## Reference Refill

### Reference Refill Process



**Figure 5-15.** Reference Refill Process\*

\*Before each run, sample and standard side capillaries are pumped.

The Reference Refill process, Figure 5-15, is used to transfer a reference gas from a reservoir bottle (usually the Reference Gas Refill unit) into the right-hand reference bellow. This gas is then used during the following measurement cycles as the reference (or standard) gas.

The Reference Refill process consists of the three following steps:

1. Dual Inlet is prepared for the transfer by pumping away all residual gas from the respective bellow and its surroundings using the high vacuum pump. See Figure 5-16.




Figure 5-16. Reference Refill Process - Step 1





Figure 5-17. Reference Refill Process - Step 2

2. The bellow is connected to the fore vacuum pump again, and the valve to the reservoir (Refill unit) is opened. Gas accumulated behind this valve is pumped away for some time (**Pump Overlay Time**). See Figure 5-17.
3. The bellow is slowly filled with gas for some time (**Refill Time**). Laminar flow must be ensured during this period of time in order to avoid fractionation. See Figure 5-18.



**Figure 5-18.** Reference Refill Process - Step 3

## Selecting Correct Amount of Reference Gas

This section treats the reference refill algorithm. It assures that the correct amount of reference gas will be selected.

Selecting the correct amount of gas to fill the reference bellows is of utmost importance for the quality of the results of the individual measurements. This is because for each sample measurement it is required to perform a proper press adjust.

Consequently, adjust the parameter **Refill Time** in Figure 3-29 and [Table 3-6](#) in such a way that after Reference Refill the amount of reference gas in the reference bellow yields a signal corresponding to the smallest sample signal that you want to measure.

Use plots like Figure 5-9 and Figure 5-10 to determine the sensitivity of your Kiel IV Carbonate Device and to decide how much signal corresponds to the smallest sample you want to analyze.

After this step, adjust the parameter **VM1 Expansion** (in Figure 3-34 and Table 3-11) according to this value and the compression of your bellow.

**Example:** If the smallest useful sample size results in a signal of 800 mV and the compression rate of your bellow is 10, then the expansion should take place above 8 V of signal. Or, if Figure 5-10 is the response curve of your IRMS, set **VM1 Expansion** to 950  $\mu$ bar.



## Chapter 6 Technical Information

**Note** This section is intended for use by trained Thermo Electron personnel only. Thermo Electron discourages use by and denies liability for the consequences of use by other than Thermo Electron personnel. ▲

This chapter contains the following topics:

- “Spare Parts and Consumables” on page 6-3
- “Valve Unit” on page 6-13
- “Valve Replacement” on page 6-14
- “IAEA Primary Standards” on page 6-17
- “Checking for Internal Leaks” on page 6-18
- “Maintenance” on page 6-26
- “Programming Information” on page 6-26
- “Compressed Air Supply” on page 6-33
- “Vacuum Schematic” on page 6-35
- “Circuit Diagrams” on page 6-36



## Spare Parts and Consumables

Spare parts and consumables of the Finnigan Kiel IV Carbonate Device have been subsumed in a kit with Part No. 119 1160. Some of these are frequently used and shown in [Table 6-1](#).


**Table 6-1.** Some Spare Parts and Consumables (Part No. 119 1160)

Position	Part No.	Designation	Quantity	Figure
1	059 8651	capillary for acid valve	2	Figure 6-2
2	109 4301	oil prevac. pump P 3 -3 ltrs	1	Figure 6-3
2	056 7830	flexible tubing (acid resistant)	0.5 m	Figure 6-4
3	017 2350	fiber reservoir TPH060/62/64	1	
3	059 8750	drop counter for acid valve*	2	Figure 4-17
4	115 3560	lamp for heating carbinet	2	
4	059 8671	gasket TEF 11.9x9x0.25 (pack of 3)	2	Figure 6-5
5	055 4440	O ring, 12.37x2.62 Viton	12	
6	054 5270	gasket gold o.d. 38xi.d. 36 5	2	Figure 6-6
7	067 1182	capillary tube, heatable	1	Figure 6-7
8	119 1170	pinch valve	1	Figure 6-8
9	065 3010	membrane, complete	2	Figure 6-9
10	065 3041	stamp Au, seal tip d=3 9	2	Figure 6-10
11	055 3140	jacket ring	2	Figure 6-11
12	115 7670	frit 1/4" o.d.x1/32", stainless steel	1	Figure 4-20
13	075 4960	sample vial, KS 19/9- 50 LG	24	Figure 6-12
14	106 9490	proximity switch, 4 mm	2	Figure 6-13

\* 6 in Figure 6-1.




**Figure 6-2.** Capillary for Acid Valve (Part No. 059 8651)



**Figure 6-3.** Oil for Fore Vacuum Pump (Part No. 109 4301)




**Figure 6-4.** Flexible Acid Resistant Tubing (Part No. 056 7830)




**Figure 6-5.** Teflon Gasket for Acid Valve (Part No. 059 8671)



**Figure 6-6.** Gold Gasket (Part No. 054 5270)



**Figure 6-7.** Heatable Capillary (Part No. 067 1182)



**Figure 6-8.** Pinch Valve (Part No. 119 1170)



**Figure 6-9.** Complete Membrane (Part No. 065 3010)



**Figure 6-10.** Gold Stamp (Part No. 065 3041)



**Figure 6-11.** Jacket Ring (Part No. 055 3140)



**Figure 6-12.** Sample Vial (Part No. 075 4960)




**Figure 6-13.** Proximity Switch (Part No. 106 9490)



**Figure 6-14.** Gold Gasket for Microvolume (Part No. 055 1010)



**Figure 6-15.** Copper Shim (Part No. 100 7730)



**Figure 6-16.** Tool to Adjust Turret Readout (Part No. 115 7390)

## Installation Kit

The Installation Kit, Part No. 115 7800, comprises some important parts which are summarized in [Table 6-2](#).

**Table 6-2.** Parts of Installation Kit (Part No. 115 7800)

Position	Part No.	Designation	Quantity
1	048 2610	distributor liquid nitrogen	1
4	114 5600	installation kit	1
5	114 7090	calcium carbonate	1
6	111 2640	ortho-phosphoric acid	500 g
7	111 3791	Finnigan Kiel IV Carbonate Device Operating Manual	1
8	041 4130	valve for liquid nitrogen	1

## Acid Valve




Figure 6-17. Parts of Acid Valve




Figure 6-18. Acid Valve (Part No. 106 9450)


The acid valve, Part No. 106 9450, is shown in Figure 6-17 and Figure 6-18. Some important parts are summarized in [Table 6-3](#).

**Table 6-3.** Parts of Acid Valve (Part No. 106 9450)


Position	Part No.	Designation	Quantity
1	106 9440	valve body for acid valve	1
2	106 0480	ring	1
3	115 7440	valve mount.	1
4	059 8651	capillary for acid valve	1
5	059 8700	contact pin for acid valve	1
6	059 8750	drop counter for acid valve	1
11	119 1170	pinch valve	1
12	059 8671	gasket, TEF, 11.9x9x0.25 (pack of 3)	1
17	045 3750	cylindrical screw, ISO 1207, M3x8, Nylon	1
18	045 0750	cylindrical screw, ISO 1207, M3x6, A4	3
19	046 0610	hexagonal nut, ISO 4032, M3, A4	2
20	115 7680	countersunk socket screw, ISO 10642, M5x16, A2	4
22	047 0210	disc, ISO 7092-3, A4	5
24	056 7830	flexible tubing, acid resistant	0.250 m
25	106 9460	distance ring, DN 6	1
26	055 4440	O ring, 12.37x2.62, Viton	1
30	052 1160	stainless nut, 1/4" o.d.	1
31	052 1330	front ferrule, Teflon, tube o.d. 1/4"	1
32	052 1340	back ferrule, Teflon, tube o.d. 1/4"	1
33	111 3460	mounting angle, acid valve	1
34	115 7670	frit, 1/4" o.d. x 1/32", stainless steel	1

## Trapping Volume

The double-walled trapping volume, Part No. 100 7740, is shown in Figure 6-19. It comprises some important parts which are depicted in Figure 6-20 and summarized in [Table 6-4](#).



**Figure 6-19.** Trapping Volume (Part No. 100 7740)



**Figure 6-20.** Important Parts of Trapping Volume\*


\*The numbers refer to [Table 6-4](#).

**Table 6-4.** Parts of Trapping Volume (Part No. 100 7740)

Position	Part No.	Designation	Quantity
1	100 7700	distribution head	1
2	100 5520	inner tube	1
3	100 5530	outer tube	1
4	100 5540	mounting screw, M14x1.5	1
5	047 0850	disc, ISO 7092-6-A2	1
6	055 1240	gasket, gold, o.d. 11.5x i.d. 10	1
7	100 5570	gasket, gold, o.d. 5x i.d. 3.5	1
8	100 7730	copper shim, 0.5 mm	1

## Cooling Unit

The cooling unit, Part No. 079 2400, is shown in Figure 6-21. It comprises some important parts which are summarized in [Table 6-5](#).



**Figure 6-21.** Cooling Unit (Part No. 079 2400)\*

\* Positions 7, 9, 11 and 37 are enclosed separately. See [Table 6-5](#).

**Table 6-5.** Important Parts of Cooling Unit (Part No. 079 2400)\*

Position	Part No.	Designation	Quantity
1	067 4300	cooling cascade	1
3	067 4390	tube	1
4	067 4400	funnel cooling unit	1
5	067 4410	lid	1
6	037 0650	clamp for cable, 6.4	1
7	050 4710	clamp for lid in cover	2
8	045 0760	cylindrical screw, ISO 1207, M3x8-A4	4
9	045 0790	cylindrical screw, ISO 1207, M4x8-A4	2
10	047 0210	disc, ISO 7092-3-A4	4

## Technical Information

### Spare Parts and Consumables

**Table 6-5.** Important Parts of Cooling Unit (Part No. 079 2400) \*, continued

Position	Part No.	Designation	Quantity
11	047 0040	disc, ISO 7089-4-A4	2
14	033 0060	eyelet, 3A2, DIN 41496	1
16	034 1630	shrink tubing, 3.2 NF	0.1 m
17	034 2030	shrink tubing, 6.4 SW	2 cm
18	034 2040	shrink tubing, 9.5 SW	0.1 m
19	035 0610	wire, 1.0 E-CU	0.25 m
20	032 2230	plug, T 3475/2	1
21	038 5820	heating cartridge, 80 W, 24 V, 40 mm	1
22	006 5160	measuring resistor, Pt 100	1
24	046 7500	feedthrough	2
25	033 0740	socket contact	2
26	035 1410	Teflon cord, 0.6 NF	0.8 m
27	047 2950	disc, 2.2x5.5x0.5, VKF	2
30	075 1960	DR WID, 220 0 5 W	1
31	075 1970	tubing, 4.0x2.5, TEF	10 cm
37	102 8300	disc	1


\*See Figure 6-21.

## Microvolume

The Microvolume, Part No. 116 3150, is shown in Figure 6-22 and Figure 6-23 together with its important parts. These are summarized in [Table 6-6](#). Refer to “[Microvolume \(Trap 2\)](#)” on [page 2-23](#) as well.



**Table 6-6.** Important Parts of Microvolume (Part No. 116 3150)

Position	Part No.	Designation	Quantity
1	116 3160	valve block	1
2	065 3001	valve unit	2
3	079 2800	front ferrule, 1/4", Swagelok	3
4	079 2810	back ferrule, 1/4", Swagelok	3
5	052 1160	stainless nut, 1/4" o.d.	3
7	055 1010	gasket, gold, o.d. 6 25, i.d. 4 75	2
8	075 4460	pressurizing screw	1
10	041 2300	heat transfer tube	1
11	078 3330	trapping volume, CO <sub>2</sub>	1
12	058 2330	washer	1
13	117 7150	pressurized air fitting, KQ2S	2



**Figure 6-22.** Microvolume (Part No. 116 3150)\*

\*The numbers refer to [Table 6-6](#).



- 1 valve block
- 8 pressurizing screw
- 10 heat transfer tube
- 11 trapping volume CO<sub>2</sub>,  
50 µl  
(with gold gasket, 7)
- 12 washer


**Figure 6-23.** Microvolume and its Parts\*

\*The numbers refer to [Table 6-6](#).

## Valve Unit

The valve unit, Part No. 065 3001, is shown in Figure 6-24. Important parts of it have been summarized in [Table 6-7](#).

**Technical Information**  
Valve Replacement



**Figure 6-24.** Valve Unit (Part No. 065 3001)\*

\* 1) indicates parts that have been slightly greased with Li grease. Parts 2, 4, 11 and 12 are mounted.

**Table 6-7.** Important Parts of Valve Unit (Part No. 065 3001)

Position	Part No.	Designation	Quantity
2	065 3010	membrane, complete	1
4	065 3050	pressure unit	1
5	065 3030	piston	1
6	065 3041	stamp, gold, seal tip, d=39	1
7	065 3060	cover plate	1
10	045 3420	cylindrical screw, ISO4762, M6x20-A4	4
11	047 3430	guard ring, 10x1, DIN 471	1
12	043 1570	spring plate	8
13	055 3140	set of sealing rings	1
14	054 5270	gasket, gold, o.d. = 38xi.d. 36 5	1

## Valve Replacement

**Table 6-8.** Parts Needed for Valve Replacement

Description	Part No.
membrane (complete)	065 3010
gold stamp	065 3041
gold gasket	054 5270

In order to rebuild one of the stainless steel valves inside the Kiel IV Carbonate Device oven, you need the parts listed in [Table 6-8](#). Follow the steps below:

1. Detach the vial on the side where the valve is leaking. As this sets the rear side of the valve to atmospheric pressure, you can take it apart.
2. Open the valve. This closes the compressed air to the normally open valves.
3. Consider removing the acid reservoir (to avoid breakage)
4. Remove the upper oven shelf to get better access to the work area.
5. Remove the four Allen screws from the faceplate. We use a special, long T-handled 5 mm Allen wrench which has a ball on the end so the wrench can be used at many angles. See Figure 6-25.



**Figure 6-25.** Tools for Valve Replacement

6. Pull the square faceplate away from the round piston assembly by hand. There is no need to disconnect the compressed air lines.
7. Rotate the piston assembly and find the air hole. Use the long tube which comes with “canned air” to blow the membrane assembly apart from the valve housing.
8. Place the membrane assembly with the set of compressed spring washers face up on a counter. Use the special tool, which looks like ‘needle nose pliers’, to remove the guard ring which frees the valve membrane from the membrane assembly. These pliers work opposite normal ones: when you squeeze the pliers they open the ring. The ends of the pliers fit into small holes in the guard ring.

**Note** Do not open the ring too far as it will stretch and have to be replaced. Remember that when you remove the ring, the compressed washers below will expand upwards. ▲

9. Remove the series of curved spring washers by simply lifting the valve housing off of the valve membrane, keeping the 'spring washer' unit intact. The valve membrane assembly will be left on the counter to be disposed of.
10. Gently place the new gold stamp into the new valve membrane. Push the gold stamp into the membrane using finger pressure only. Place the membrane assembly into a small vise and exert equal pressure on the gold stamp and the back of the valve membrane.

**Caution** Do not apply pressure to the edges of the stamp! Instead, apply firm even pressure to the middle of the stamp to avoid damaging the stamp edges. ▲

11. Insert the membrane assembly up through the valve housing, ensuring that the membrane post fits through the spring washers.
12. Using the plier tool, slightly open the guard ring and press down on the washers so the ring can fit into the groove of the membrane post.

**Caution** Be very careful while handling the attached gold stamp on the bottom! Do not scuff the edges! ▲

13. Clean the piston assembly housing and lightly grease the outer edge of the O ring seal inside. The valve will not work without some lubrication. Slide the piston assembly back into the housing
14. Dust off the valve block assembly and place the new gold O ring seal against the block body.
15. Place the piston assembly against the gold O ring seal and the valve block body.
16. Place the valve faceplate with the still attached compressed air line against the valve membrane assembly, making sure that the gold O ring seal does not slip. Reattach the four Allen screws as you would a car tire, rotating the tightening around each screw. Tighten

rather hard. The four Allen screws may need to be re-tightened after the assembly has been at vacuum.

17. Seat the new gold stamp by opening/closing the valve 20-30 times and check for leaks. Heat up the Kiel IV Carbonate Device oven to 70 °C to remove water before running samples.

## Definitions

**Table 6-9.** Definitions

Term	Definition
faceplate	Square stainless steel plate to which compressed air is attached. It is attached to the valve block assembly by four Allen screws
valve block assembly	Part of valve left in oven when faceplate and piston assembly are removed
piston assembly	membrane assembly + piston housing
membrane assembly	valve membrane + gold stamp + spring washers + guard ring

## IAEA Primary Standards

Name	Nature		Isotopic ratio	δ‰	Reference standard
V-SMOW	water	2H/1H	(155.761 ± 0.05) × 10e-6 (1) (155.751 ± 0.08) × 10e-6 (2) (155.601 ± 0.12) × 10e-6 (3)	0	V-SMOW
		18O/16O	(2005.20 ± 0.45) × 10e-6 (4)	0	V-SMOW
		17O/16O	(379.91 ± 0.8) × 10e-6 (5)	0	V-SMOW
SLAP	water	2H/1H	(89.021 ± 0.05) × 10e-6 (1) (89.12 ± 0.07) × 10e-6 (2) (88.88 ± 0.18) × 10e-6 (3)	-428.0 (6)	V-SMOW
		18O/16O	(1893.91 ± 0.45) × 10e-6 (7)	-55.50 (6)	V-SMOW
NBS-19	calcite	13C/12C		1.95 (8)	V-PDB
		18O/16O		-2.20 (8)	V-PDB
				28.6 (9)	V-SMOW
<b>intercomparison materials</b>					
GISP	water	2H		-189.73 ± 0.87	V-SMOW
		18O		-24.784 ± 0.075	V-SMOW
NBS-18	calcite	13C		-5.029 ± 0.049	V-PDB
		18O		-23.035 ± 0.172	V-PDB
IAEA-CO-1	calcite	13C		2.48 ± 0.025	V-PDB
		18O		-2.437 ± 0.073	V-PDB
IAEA-CO-8	calcite	13C		-5.749 ± 0.063	V-PDB
		18O		-22.667 ± 0.187	V-PDB
IAEA-CO-9	BaCO <sub>3</sub>	13C		-47.119 ± 0.149	V-PDB
		18O		-15.282 ± 0.093	V-PDB

**Figure 6-26.** IAEA Primary Standards

Calibrating vs. international standards requires users to have their own specimens of Primary Standards. Primary Standards are exclusively distributed by the IAEA via agencies in Europe and the US. The reference list shown in Figure 6-26 is taken from IAEA TECDOC 825, and the IAEA Analytical Quality Control Services Reference Materials Catalogue 2002-2003. It is very difficult to obtain homogeneity in both isotopes,  $^{13}\text{C}$  and  $^{18}\text{O}$ .

**Note** Refer to IAEA-TECDOC-825: Reference and intercomparison materials for stable isotopes of light elements. Proceedings of a consultants meeting held in Vienna, 1-3 December 1993. International Atomic Energy Agency (IAEA). ▲

**Note** Also refer to Chapter 5.2 Environmental Level, pp. 55 in: IAEA Analytical Quality Control Services Reference Materials Catalogue 2002-2003. First edition, January 2002. Edited by Analytical Quality Control Services, International Atomic Energy Agency, P.O. Box 100, A-1400 Vienna. ▲

## Checking for Internal Leaks




Figure 6-27. Checking VM1 Pressure Increase of Line 1

**Note** Pressure and temperature units are off-scale in Figure 6-27. ▲

Sources of a leak in order of probability:

- The vial rim is broken.
- The electrical feedthrough of the drop counter is not leaktight.
- The amount of phosphoric acid is too high in either of the acid drop vials (position 2).
- The Teflon washer of the drop counter is missing/not leaktight.
- The acid valve unit is dirty due to phosphoric acid accumulation. To clean it, see “[Cleaning Acid Valve](#)” on [page 4-14](#).
- The Swagelok connections to the pneumatic valve unit (V12, V13 and V22, V23) are leaking.
- The pneumatic valve unit is not leaktight.

**Note** If the pneumatic valve unit is not leaktight, perform a leak test. ▲

In case of a leak, use the following procedure to resolve the problem.

Depending on the line that shows the problem, set the Kiel IV Carbonate Device valves to the positions shown in Figure 6-27. For a problem in line 1, open V12, V1 and V4. For a problem in line 2, open V22, V1 and V4.

If all parts operate properly, close all other valves. The pressure of VM1 should not exceed 50-70 mbar. See upper right corner in Figure 6-27.

## High Vacuum Side

### Leak Check of Pneumatic Valves


To check the pneumatic valves on the high vacuum side of the Kiel IV Carbonate Device for small leaks, proceed as follows (see Figure 6-27):

**Note** First, wait for 2 min. The VM1 pressure must not increase! ▲

**Note** After each step, check VM1 pressure. It must never increase! ▲


1. Pump the high vacuum side with V1, V2, V3, V4, V5 open. Disconnect one of the vials. Pump until VM1 pressure is at background level, that is below 50-70 mbar. See Figure 6-28.

**Technical Information**  
Checking for Internal Leaks



**Figure 6-28.** Leak Check of Pneumatic Valves - Step 1

2. Close all valves on the high vacuum side. See Figure 6-29.



**Figure 6-29.** Leak Check of Pneumatic Valves - Step 2

3. To be able to open V12 or V22, lift the spring plate with one finger towards the proximity switch (Figure 4-18). Fill until V1 with air.
4. Close V12 or V22 (dependent on where the vial is not connected).
5. Open V1 and V2. Leave V4 and V5 closed. See Figure 6-30.





Figure 6-30. Leak Check of Pneumatic Valves - Step 5

6. Perform a **Tune Scan** on m/z 28 with COV sample side to the Kiel IV Carbonate Device open (V15, V16 closed). See Figure 6-31.



**Figure 6-31.** Leak Check of Pneumatic Valves - Step 6

To record m/z 28 from the Kiel IV Carbonate Device, the Dual Inlet system must be adjusted as shown in Figure 6-31.

A leak is possible between V3 pneumatic valve and the IRMS (that is, COV valve block, sample capillary, tubing from Changeover Valve to ion source inlet, needle valve).

### Carbonate Side Leaks

Carbonate side leaks may occur at V3 pneumatic valve, at Microvolume or at the capillary connection to the Kiel IV Carbonate Device and the Changeover Valve (COV).

#### Checking V2

To check V2, proceed as follows:

**Note** First, wait for 2 min. The VM1 pressure must not increase! ▲

**Note** After each step, check VM1 pressure. It must never increase! ▲

1. Pump out the air in the lines from V2 to the IRMS via V5. See Figure 6-32.




Figure 6-32. Checking V2 - Step 1




Figure 6-33. Checking V2 - Step 2

2. Close V5. See Figure 6-33. Observe VM1 pressure. If it does not rise, no leak is present. If it rises however, a leak is possible behind V2 at pneumatic valve unit or at the transfer tubing.

### Checking V1




Figure 6-34. Checking V1 - Step 1





Figure 6-35. Checking V1 - Step 2

**Note** First, wait for 2 min. The VM1 pressure must not increase! ▲

**Note** After each step, check VM1 pressure. It must never increase! ▲

1. Close V1. See Figure 6-34.
2. Open V2, V3, V4 and V5. See Figure 6-35.
3. Close V4 and V5. See Figure 6-36.  
Observe VM1 pressure. If it does not rise, no leak is present. If it rises however, a leak is possible behind V1 at the four locations:
  - pneumatic valve unit or
  - transfer tubing or
  - VM1 connection or
  - water trap



**Figure 6-36.** Checking V1 - Step 3

## Maintenance

### Turbo Pumps and Fore Pumps

Change the oil reservoir of the turbo pumps once a year. Check the oil level of the fore pumps frequently and exchange the complete oil once a year. Refer to manufacturer's manual.

### O Ring Seals

Remove acid and any residuals from the O ring seals of the acid valve prior to starting a run. The O ring seals may be slightly greased using a high-vacuum grease, that is, Apiezon H.

### Autocool Unit

Between runs, allow the Autocool Unit to heat up to room temperature by removing the dewar. This serves to completely remove frozen water. Discard the water that accumulates in the dewar.

**Note** If no liquid nitrogen is transported from the funnel into the cascade, the freezing temperature of  $-190\text{ }^{\circ}\text{C}$  cannot be reached. Isodat 2.5 will then show a temperature error. As it can no longer be frozen,  $\text{CO}_2$  is lost. During cooling, isotope fractionation occurs. ▲

**Note** If improper transport of liquid nitrogen between funnel and cascade occurs, ice droplets are clogging inside the funnel tubing. Release the liquid nitrogen dewar in order to let the ice melt and to release most of the frozen water. ▲

**Note** PTFE tubing connects the funnel to the cooling cascade. Due to PTFE chips sticking inside the tubing, transport of liquid nitrogen into the cascade is interrupted. File off the PTFE tubing in order to let the transport of liquid nitrogen work properly again. ▲

**Note** If temperatures above  $120\text{ }^{\circ}\text{C}$  cannot be obtained by the heating cartridge<sup>1</sup>, it must be exchanged (it tends to age). As heating cartridges sometimes cannot be released, the Autocool Unit must be exchanged. At temperatures below  $120\text{ }^{\circ}\text{C}$ , it is impossible to heat out water from the trapping system. Isotope fractionation will lead to erroneous isotope ratios due to  $\text{H}_2\text{O-CO}_2$  collisions within the ion source. ▲

## Programming Information

### Board Base Addresses


The board base addresses are summarized in [Table 6-10](#).

---

<sup>1</sup>the 30 W heater, 2 in Figure 2-32

**Table 6-10.** Board Base Addresses

Usage	Code
Inlet Controller Internal	1
Inlet Controller External Options	2
Inlet Controller Kiel IV Carbonate Device primary	B
Data Logger MS	4
Inlet Controller Kiel IV Carbonate Device secondary	5
Source Control Gnd	6
Source Control HV	7
Power Distributor MS	8
Power Distributor Kiel IV Carbonate Device	9
PeriCon	C



**Figure 6-37.** Adjusting Base Address of Power Distribution Board

Figure 6-37 depicts how to adjust the base address of the Power Distribution board. **1** denotes the address selection jumper.

## Device Addresses

**Table 6-11** summarizes the hardware addresses used within Kiel IV Carbonate Device (e.g. for plugs, contacts and pins).

**Table 6-11.** Device Addresses - Complete

Usage	Board	Jack	Pins	Logical Name in Isodat	Type of IO	Hex Address
Trap 1 temperature true value	IC 1			Get Trap 1 Temperature	ADC	0xB37
Acid true temperature	IC 1	J210		Get Acid Temperature	ADC	0xB37
Trap 2 temperature true value	IC 2	J214		Get Trap 2 Temperature	ADC	0x537
Position Sensing Array	IC 1	J204		Turrect Raw	DAC	0xB07
Motor torque setting				Turret torque	DAC	0xBFB
Turret Position				Turret position	DAC	0xBFF
Trap 1 temperature Set value	IC 1	J214		Set Trap 1 Temperature	DAC	0xB30
Trap 2 temperature Set value	IC 2	J214		Set Trap 2 Temperature	DAC	0x530
Proximity Switch Line 1	IC 2	J209		Bottle Line 1	DIO	0x507
Proximity Switch Line 2	IC 2	J208		Bottle Line 2	DIO	0x507
Drop Counter Line 1	IC 1	J215		Drop Line 1	DIO	0xB07
Drop Counter Line 2	IC 1	J215		Drop Line 2	DIO	0xB07
Piston Sensor Line 1	IC 1	J208		Piston Sensor 1	DIO	0xB07
Piston Sensor Line 2	IC 1	J207		Piston Sensor 2	DIO	0xB07
Motor control direction	IC 1				DIO	0xB00
Motor control enable	IC 1				DIO	0xB00
Motor control direction	IC 1				DIO	0xB00
Trap 1 Heater Disable	IC 1			Trap1 Heater Off	DIO	0xB00
Trap 2 Heater Disable	IC 2			Trap2 Heater Off	DIO	0x500
Valve 22	IC 1	J225	4	Valve 22	DIO	0xB28
Valve 23	IC 1	J225	5	Valve 23	DIO	0xB28
Valve 12	IC 1	J225	6	Valve 12	DIO	0xB28
Valve 13	IC 1	J225	7	Valve 13	DIO	0xB28
Piston Line 1 up	IC 1	J225	20	Piston Line 1 up	DIO	0xB28
Piston Line 2 up	IC 1	J225	21	Piston Line 2 up	DIO	0xB28
Piston Line 1 down	IC 1	J225	22	Piston Line 1 down	DIO	0xB28
Piston Line 2 down	IC 1	J225	23	Piston Line 2 down	DIO	0xB28
Pinch valve Line 1	IC 1	J217	4	Acid valve Line 1	DIO	0xB29
Pinch valve Line 2	IC 1	J217	5	Acid valve Line 2	DIO	0xB29
Valve 7	IC 1	J217	6	Valve 7	DIO	0xB29
Valve 8	IC 1	J217	7	Valve 8	DIO	0xB29
Valve 3	IC 1	J227	6	Valve 3	DIO	0xB2A
Valve 5	IC 1	J227	7	Valve 5	DIO	0xB2A
Valve 1	IC 1	J227	20	Valve 1	DIO	0xB2A
Valve 2	IC 1	J227	21	Valve 2	DIO	0xB2A
Valve 4	IC 1	J227	22	Valve 4	DIO	0xB2A

**Table 6-11.** Device Addresses - Complete, continued

Usage	Board	Jack	Pins	Logical Name in Isodat	Type of IO	Hex Address
Trap 2 cooling resistor	IC 2	J214			DIO	0x52E
Trap 1 cooling resistor	IC 1	J214			DIO	0xB02
Turbo pump Error	PD				DIO	0x902
Turbo pump 50%	PD				DIO	0x903
Control Refill Valve	PD			Refill Direct	DIO	0x901
LN2 Refill Unit Enable (LED)	PD			LN2 Refill Direct	DIO	0x901
Host Connection	PD			Host Connection	DIO	0x901
VM 1 pressure true value	IC 1	J215		VM 1	PM	0xB39
VM 2 pressure true value	IC 1	J215		VM 2	PM	0xB39
LN2 Level Sensor				LN 2 Pressure	PM	0x539

Table 6-12 shows DIO parameters as subset a of the device addresses.

**Table 6-12.** DIO Parameters

Usage	Opcode	Control Code	Address Code
Proximity Switch Line 1	1	10	0
Proximity Switch Line 2	1	11	0
Drop Counter Line 1	1	8	0
Drop Counter Line 2	1	9	0
Piston Sensor Line 1		11	
Piston Sensor Line 2		12	
Motor control direction		0	
Motor control enable		1	
Motor control direction		2	
Trap 1 Heater Disable	1	9	0
Trap 2 Heater Disable	1	9	0
Valve 22	1	0	
Valve 23	1	1	0
Valve 12	1	2	0
Valve 13	1	3	0
Piston Line 1 up	1	4	0
Piston Line 2 up	1	5	0
Piston Line 1 down	1	6	0
Piston Line 2 down	1	7	0
Pinch valve Line 1	1	0	0
Pinch valve Line 2	1	1	0
Valve 7	1	2	0

**Table 6-12.** DIO Parameters, continued

Usage	Opcode	Control Code	Address Code
Valve 8	1	3	0
Valve 3	1	2	0
Valve 5	1	3	0
Valve 1	1	4	0
Valve 2	1	5	0
Valve 4	1	6	0
Trap 2 cooling resistor		7	
Trap 1 cooling resistor		7	
Turbo pump Error		7	
Turbo pump 50%		9	
Control Refill Valve	1	4	0
LN2 Refill Unit Enable (LED)	1	12	0
Host Connection	1	0	0

Table 6-13 shows ADC parameters as subset a of the device addresses.


**Table 6-13.** ADC Parameters

Usage	Min.	Max.	Offset	Gradient 1	Gradient 2	Gradient 3	Log Calc.
Trap 1 temperature true value	0	4096	-200	0.0866699	0	0	
Acid true temperature	0	4096	-100	0.051	0	0	
Trap 2 temperature true value	0	4095	-200	0.0866699	0	0	
Trap 1 temperature Set value			-200	0.0866699	0	0	
Trap 2 temperature Set value			-200	0.0866699	0	0	
VM 1 pressure true value	0	4096	-532	1000	0	0	
VM 2 pressure true value	0	4096	-12.888	3.7087	-0.3733	0.0169	y
LN2 Level Sensor	0	4095	175.47	-0.188679	0	0	

Table 6-14 shows DAC parameters as subset a of the device addresses.

**Table 6-14.** DAC Parameters

Usage	Min.	Max.	Address Lo	Address Hi	Cut DAC Steps
Position Sensing Array	0	65535	0	0	0
Motor torque setting	0	65535	0	0	0
Turret Position	1	24	7	0	0
Trap 1 temperature Set value	0	4095	0	0	104
Trap 2 temperature Set value	0	4095	0	0	104



- 1 motor
- 2 trap
- 3 data in
- 4 power
- 5 drop counter and VM1
- 6 position sensor array
- 7 valves

**Figure 6-38.** Connections on Inlet Control Board

Figure 6-38 shows the connections on the Inlet Control board. More precisely, two identical boards are arranged in stacked order on top of each

other. The upper one has the board base address B, whereas the board base address of the lower one is 5 (see Table 6-10).

As a detail of Figure 6-38, Figure 6-39 shows where the sensitivity of the drop counters is adjusted.




Figure 6-39. Adjusting Sensitivity of Drop Counters

## Registry Use

### LineX Sample Position

- Function **Disconnect** resets variable according to **Line** status to 0 when completed successfully.
- Function **Connect** sets respective variable to **Position** when completed successfully.


**Note** Manual connection/disconnection is not detected. ▲

### Carousel Position


- controlled by function **MoveToPosition**
- not actively used with Kiel IV Carbonate Device since readout of position is always possible

## Compressed Air Supply

The compressed air supply is shown in Figure 6-40 and Figure 6-41.




**Figure 6-40.** Compressed Air Supply



**Figure 6-41.** Compressed Air Supply - Schematic

# Vacuum Schematic



**Figure 6-42.** Kiel IV Carbonate Device - Vacuum Schematic







# Glossary

**A** ampere

**ac** alternating current

**ADC** analog-to-digital converter

**AP** acquisition processor

**APCI** atmospheric pressure chemical ionization

**API** atmospheric pressure ionization

**ASCII** American Standard Code for Information Interchange

**b** bit

**B** byte (8 b)

**baud rate** data transmission speed in events per second

**°C** degrees Celsius

**cfm** cubic feet per minute

**CI** chemical ionization

**CIP** carriage and insurance paid to

**cm** centimeter

**cm<sup>3</sup>** cubic centimeter

**CPU** central processing unit (of a computer)

**CRC** cyclic redundancy check

**CRM** consecutive reaction monitoring

**<Ctrl>** control key on the terminal keyboard

**d** depth

**Da** dalton

**DAC** digital-to-analog converter

**dc** direct current

**DDS** direct digital synthesizer

**DEP™** direct exposure probe

**DS** data system

**DSP** digital signal processor

**EI** electron ionization

**EMBL** European Molecular Biology Laboratory

**<Enter>** Enter key on the terminal keyboard

**ESD** electrostatic discharge

**ESI** electrospray ionization

**eV** electron volt

**f** femto (10<sup>-15</sup>)

**°F** degrees Fahrenheit

**.fasta file** extension of a SEQUEST search database file

**FOB** free on board

**ft** foot

**FTP** file transfer protocol

**g** gram

**G** Gauss; giga (10<sup>9</sup>)

**Glossary:** GC

**GC** gas chromatograph; gas chromatography

**GC/MS** gas chromatograph / mass spectrometer

**GND** electrical ground

**GPIB** general-purpose interface bus

**GUI** graphical user interface

**h** hour

**h** height

**HPLC** high-performance liquid chromatograph

**HV** high voltage

**Hz** hertz (cycles per second)

**ICIS™** Interactive Chemical Information System

**ICL™** Instrument Control Language™

**ICP** inductively coupled plasma

**ICP-OES** inductively coupled plasma optical emission spectroscopy

**ID** inside diameter

**IEC** International Electrotechnical Commission

**IEEE** Institute of Electrical and Electronics Engineers

**in.** inch

**I/O** input/output

**IRMS** isotope ratio mass spectrometer

**k** kilo ( $10^3$ , 1000)

**K** kilo ( $2^{10}$ , 1024)

**KEGG** Kyoto Encyclopedia of Genes and Genomes

**kg** kilogram

**l** length

**L** liter

**LAN** local area network

**lb** pound

**LC** liquid chromatograph; liquid chromatography

**LC/MS** liquid chromatograph / mass spectrometer

**LED** light-emitting diode

**LHe** liquid helium

**LN2** liquid nitrogen

**μ** micro ( $10^{-6}$ )

**m** meter

**m** milli ( $10^{-3}$ )

**M** mega ( $10^6$ )

**M+** molecular ion

**MB** Megabyte (1 048 576 bytes)

**MH+** protonated molecular ion

**min** minute

**ml** milliliter

**mm** millimeter

**MS** mass spectrometer; mass spectrometry

**MS**  $MS^n$  power: where  $n = 1$

**MS/MS**  $MS^n$  power: where  $n = 2$

**MS<sup>n</sup>**  $MS^n$  power: where  $n = 1$  through 10

**m/z** mass-to-charge ratio

<b>n</b> nano ( $10^{-9}$ )	<b>RS-232</b> industry standard for serial communications
<b>NCBI</b> National Center for Biotechnology Information (USA)	<b>s</b> second
<b>NIST</b> National Institute of Standards and Technology (USA)	<b>SIM</b> selected ion monitoring
<b>OD</b> outside diameter	<b>solids probe</b> direct insertion probe
$\Omega$ ohm	<b>SRM</b> selected reaction monitoring
<b>p</b> pico ( $10^{-12}$ )	<b>SSQ</b> ™ single stage quadrupole
<b>Pa</b> pascal	<b>TCP/IP</b> transmission control protocol / Internet protocol
<b>PCB</b> printed circuit board	<b>TIC</b> total ion current
<b>PID</b> proportional / integral / differential	<b>Torr</b> torr
<b>P/N</b> part number	<b>TSQ</b> ™ triple stage quadrupole
<b>P/P</b> peak-to-peak voltage	<b>u</b> atomic mass unit
<b>ppm</b> parts per million	<b>V</b> volt
<b>psig</b> pounds per square inch, gauge	<b>V ac</b> volts alternating current
<b>RAM</b> random access memory	<b>V dc</b> volts direct current
<b>RF</b> radio frequency	<b>vol</b> volume
<b>RMS</b> root mean square	<b>w</b> width
<b>ROM</b> read-only memory	<b>W</b> watt



# Index

## Numerics

- 0 V input 2-15
- 006 5160 6-12
- 017 2350 6-3
- 032 2230 6-12
- 033 0060 6-12
- 033 0740 6-12
- 034 1630 6-12
- 034 2030 6-12
- 034 2040 6-12
- 035 0610 6-12
- 035 1410 6-12
- 037 0650 6-11
- 038 5820 6-12
- 041 2300 6-12
- 041 4130 6-7
- 043 1570 6-14
- 045 0750 6-9
- 045 0760 6-11
- 045 0790 6-11
- 045 3420 6-14
- 045 3750 6-9
- 046 0610 6-9
- 046 7500 6-12
- 047 0040 6-12
- 047 0210 6-9, 6-11
- 047 0850 6-10
- 047 2950 6-12
- 047 3430 6-14
- 048 2610 2-38, 6-7
- 050 4710 6-11
- 052 1160 6-9, 6-12
- 052 1330 6-9
- 052 1340 6-9
- 054 5270 6-3–6-4, 6-14
- 055 1010 6-7, 6-12
- 055 1240 6-10
- 055 3140 6-3, 6-6, 6-14
- 055 4440 6-3, 6-9
- 056 7830 6-3–6-4, 6-9
- 058 2330 6-12
- 059 8651 6-3, 6-9
- 059 8671 6-3–6-4, 6-9
- 059 8700 6-9
- 059 8750 6-3, 6-9
- 065 3001 6-12–6-14
- 065 3010 6-3, 6-5, 6-14
- 065 3030 6-14
- 065 3041 6-3, 6-6, 6-14
- 065 3050 6-14
- 065 3060 6-14
- 067 1182 6-3, 6-5
- 067 4300 6-11
- 067 4390 6-11
- 067 4400 6-11
- 067 4410 6-11
- 075 1960 6-12
- 075 1970 6-12
- 075 4460 6-12
- 075 4960 5-5, 6-3, 6-6
- 078 3330 6-12
- 079 2400 6-11
- 079 2800 6-12
- 079 2810 6-12
- 1/16" tubing 2-13
- 100 5520 6-10
- 100 5530 6-10
- 100 5540 6-10
- 100 5570 6-10
- 100 7700 6-10
- 100 7730 6-7, 6-10
- 100 7740 6-9–6-10
- 102 8300 6-12
- 106 0480 6-9
- 106 9440 6-9
- 106 9450 6-8–6-9
- 106 9460 6-9
- 106 9490 6-3, 6-6
- 108 2840 2-21
- 109 2481 3-4–3-6, 3-14, 3-20, 5-11
- 109 4301 6-3–6-4
- 111 2640 4-28, 6-7
- 111 3460 6-9
- 111 3791 6-7
- 114 5600 6-7
- 114 7090 6-7
- 115 3560 6-3
- 115 4991 3-4, 3-14–3-15, 3-42
- 115 7390 2-34, 6-7
- 115 7440 6-9
- 115 7670 4-17, 6-3, 6-9
- 115 7680 6-9
- 115 7800 6-7
- 116 3150 6-12–6-13
- 116 3160 6-12

## Index: A

117 7150 6-12  
119 1160 6-3  
119 1170 6-3, 6-5, 6-9  
13C 5-3, 5-13–5-14, 5-18, 6-18  
18O 5-3, 5-13–5-14, 6-18  
7.7 V input 2-14–2-15

## A

About Acid tab 3-47–3-48  
absolute standard 5-14, 5-17  
Accessories bar 3-5–3-8, 3-16  
accuracy 5-4, 5-14  
acetic acid 4-29  
acetone 1-3, 4-18–4-19, 4-29–4-30  
acid capillary 4-14, 4-16, 4-19, 6-3  
acid contamination 3-48  
acid drop 2-16, 2-22, 3-27, 3-48, 3-52, 4-19, 4-23  
Acid Drop Test 3-48–3-49, 3-52, 4-21  
acid dropper capillary 2-29, 4-16, 4-19  
acid dropping timeout 3-48  
acid flow 2-3, 2-13, 2-43–2-44, 4-21  
acid flow limitation 2-8  
acid glass ports 2-13  
acid reservoir 2-3, 2-13, 2-16, 2-43, 4-21, 5-7, 6-15  
acid temperature 3-26, 3-59  
acid temperature tolerance 3-48  
acid tubing 2-13, 2-44–2-45, 3-53, 4-15, 4-21, 6-3–6-4, 6-9  
acid valve 2-2–2-3, 2-8, 2-13–2-14, 2-32, 2-43–2-45, 3-48, 3-50, 3-52, 3-55, 4-7–4-8, 4-14–4-22, 4-25, 5-21, 6-3–6-4, 6-8–6-9, 6-19, 6-26  
acid valve holder 4-15  
acid-soluble minerals 4-29  
acquisition cycle 3-22  
acquisition file 3-64  
Acquisition mode 3-4, 3-16, 3-20, 3-40  
acquisition script 3-15, 3-22, 3-64, 3-66  
acquisition time 4-4  
Acquisition window 3-5  
acquisition.isl 3-15, 3-22  
Active Pirani Gauge (APG-M) 2-10, 2-16  
actuator for compressed air 2-19  
additional valve 3-9  
address selection jumper 6-27  
Administrative Panels button 3-8  
air compressor 1-3  
air leak 4-4  
air release 2-13, 2-16  
amplification 4-4  
amplitude 3-11  
Analyzer Pumps switch 2-6  
Apiezon H 4-18, 6-26  
aragonite 5-3

argon background 4-4  
argon intensity 4-5–4-6  
argon tank 1-3, 4-4  
Arrow Down key 2-27–2-28  
Arrow Up key 2-27–2-28  
ASCII format 3-19, 3-37  
Autocool Unit 2-10, 2-21–2-26, 2-38, 4-26, 6-26  
autosampler 2-2, 2-29, 2-34, 3-33, 4-7, 4-23, 5-2, 5-7  
Autotune 2-28  
auxiliary gas 2-12, 3-45  
auxiliary gas pressure 2-12  
average  $\delta$  value 3-13

## B

background determination 3-25  
background gas 4-3  
background gas composition 4-2–4-3  
background level 4-3, 6-19  
background measurement 3-59  
background predelay 3-33  
background scan 3-33  
background vacuum 3-24  
bakeout 2-40, 4-3, 4-6–4-7, 4-12, 4-22, 5-10  
Bars tab 3-5, 3-16  
Basic bar 3-5  
Basic Functionality tab 3-44–3-45  
basic measurement 4-1  
basic test 4-21  
beam intensity 5-12  
bellow adjustment 3-45  
board base address 6-32  
boat 5-4–5-5  
Browser tab 3-16  
brush 5-5  
Burman, J. 4-26

## C

C111 2-27  
C112 2-27  
C113 2-27  
C115 2-27–2-28  
C116 2-27  
calcite 5-3  
capillary crimp 2-23  
capillary heating transformer 2-14–2-15, 4-6, 6-5  
capillary matching 4-9, 4-12  
capillary port 3-2–3-3  
capillary pump out time 3-23  
carbonate analysis 2-22, 3-25, 3-44, 3-46, 4-23, 4-28, 5-4, 5-7, 5-13  
carbonate ion 5-3

carbonate sample 3-20, 3-58, 5-3, 5-9, 5-11  
carbonate sequence 4-21  
carbonate side leak 6-22  
carousel position 6-32  
cascade 2-24–2-26, 4-20, 6-11, 6-26  
center cup 3-12, 3-23, 4-4  
Changeover Closed 3-10  
Changeover Left 3-10  
Changeover Right 3-10  
Changeover Valve 2-12–2-13, 2-23, 3-9–3-10, 3-24, 3-33, 3-59, 3-65, 4-12, 5-10, 6-21–6-22  
checking V1 6-24–6-25  
checking V2 6-23  
circuit diagram 6-36–6-37  
citric acid 4-29  
clear buffer 3-17  
Close V4 Temperature 3-46  
CO<sub>2</sub> 4-3  
CO<sub>2</sub> after Exp. 3-59  
CO<sub>2</sub> Freeze Temperature 3-46  
CO<sub>2</sub> Release Temperature 3-46  
CO2\_Craig 3-28  
CO2\_SSH 3-28  
code example 3-65  
compressed air 1-3, 2-4, 2-6, 2-12, 2-18, 2-20–2-21, 6-17  
compressed air cable 2-12  
compressed air connection 2-21, 2-31  
compressed air connector 2-12  
compressed air distributor 1-3, 2-21  
compressed air inlet 2-6  
compressed air lever 2-31–2-32  
compressed air line 6-15–6-16  
compressed air plunger 2-18  
compressed air reservoir 2-7  
compressed air signal 2-21  
compressed air supply 1-3, 6-33–6-34  
compressed air tubing 1-3  
compressed air valve 2-21, 6-15  
compression rate 5-25  
compressor 2-20  
computer failure 2-8  
configuration 2-13, 3-2, 3-5–3-6, 3-20–3-21, 3-24, 4-4  
Configurator 3-2–3-3, 3-6  
connect algorithm 4-8, 4-23  
connect/disconnect procedure 3-44, 4-23, 4-25–4-26  
consumables 6-3  
contact pipe 2-26  
contamination 2-21, 2-42, 4-22, 4-29–4-30  
Continuous Flow mode 3-60  
Control panel 2-2, 2-5–2-6, 5-2  
Control Refill Valve 6-29–6-30  
cooling finger 2-22  
cooling unit 6-11

copper shim 6-7, 6-10  
cover plate 2-29–2-30, 6-14  
covering cap 2-19  
Craig, H. 5-18  
crimp 3-25, 4-11, 4-21  
crimp position 3-9, 6-5  
cross-contamination 5-10  
CY.1 2-28  
CY.2 2-28

## D

$\delta$  value 5-14, 5-16–5-17  
 $\delta$  definition 5-16–5-17  
d.F 2-28  
d.t 2-28  
 $\delta^{13}\text{C}$  3-12, 5-14  
 $\delta^{18}\text{O}$  3-12, 5-14  
DAC parameter 6-30  
DAC step 3-59, 6-30  
data bus 2-20  
data logger MS 6-27  
data point triplet 3-22  
db 2-28  
debugging 3-22  
deionized water 4-18–4-19, 4-29  
deionized water of high quality 4-29  
delay time 3-25, 3-45, 3-48  
Delta Plus 3-45  
Delta Plus XL 3-45  
DELTA V 1-2, 3-59, 4-2–4-3  
DELTA V Operating Manual 2-12  
Delta V Plus 1-2, 3-2  
density 3-59, 4-27–4-28  
detergent 4-29  
dewar 2-3–2-5, 2-26, 2-35, 2-38–2-39, 4-26–4-27, 5-12, 6-26  
disconnect algorithm 4-23  
dishwasher 4-29  
dissolved inorganic carbon 4-29  
distilled water 4-29  
distilled water of high quality 4-29  
dolomite 5-3  
drop counter 2-44, 4-16, 4-19–4-20, 5-6, 6-3, 6-9, 6-19, 6-28–6-29, 6-31  
drop counter connection 2-44  
drop counter spring 2-29  
droplet generation 5-8  
dropping 2-13, 2-16, 3-48, 4-19  
dropping capillary 4-20  
drying chamber 4-29–4-30  
Dual Inlet argon signal 4-4  
Dual Inlet board 2-20–2-21  
Dual Inlet measurement 3-10–3-11, 3-25, 3-27, 5-11

## Index: E

Dual Inlet method 3-10  
Dual Inlet mode 2-23  
Dual Inlet Sets 3-3  
Dual Inlet Sets tab 3-3  
Dual Inlet system 1-2, 2-18–2-19, 2-21, 2-40, 3-9, 3-23–3-24, 3-47, 4-10–4-11, 5-9–5-10, 5-23, 6-2, 6-22  
Dual Inlet system turbo pump 3-9–3-10  
Dual Inlet system valve 2-19  
Dual Inlet valve failure 3-62  
Dual Inlet window 3-7–3-9  
DUO 2.5 2-16  
δvalue 3-13, 5-12–5-17  
Dyn Externals 3-44, 3-49, 4-7

## E

Editors icon 3-40  
electrical feedthrough 4-16, 6-19  
electronics cabinet 2-7–2-8  
elemental isotopic composition 3-57  
equilibration time 3-46  
error message 3-26, 3-33, 3-51, 3-59  
error of mean determination 3-57  
error status 2-5  
error tracking 3-22  
European Union Waste Electrical & Electronic Equipment (WEEE) Directive 1-1  
Evaluation tab 3-21, 3-27–3-29  
evaluation type 3-28, 3-57  
example file 3-14, 3-20, 3-29, 3-31, 3-34  
example folder 3-13, 3-20, 3-30–3-31, 3-33–3-34  
Excel Export 3-14, 3-42, 5-11  
Excel Export Editor 3-42  
Excel Export template 3-42–3-43  
Excel format 3-36–3-37  
excess pressure 2-37, 4-4  
Expansion Equilibration Delay 3-47  
Expansion Pump Time 3-47  
export file 3-36  
Export tab 3-14  
export template 3-14, 3-42, 5-11

## F

face mask 4-26  
face shield 2-26, 2-34  
fatal error 3-48  
fatal script error 3-24  
fiberline 2-8–2-9  
fiberline cable 2-12  
File Browser 3-7, 3-13, 3-36  
File Search button 3-15  
Fill Grid with Data 3-33

fill level 2-35, 3-33, 4-23  
fill level sensor 2-3–2-4  
flow limiter 2-31  
Focus Delta 3-7  
foraminifera 5-19  
fore vacuum 2-16  
fore vacuum connection 2-10  
fore vacuum exhaust 2-12  
fore vacuum gauge 2-11, 3-9  
fore vacuum pressure 3-9  
fore vacuum pump 2-4–2-5, 2-7, 2-10–2-11, 2-16, 2-42, 3-23, 5-24, 6-4  
fore vacuum section 2-7, 3-45, 4-8  
fractionation 2-40, 5-24  
frostbite 2-26, 2-34  
fume hood 1-3, 4-26  
funnel 2-24–2-26, 4-20, 6-11, 6-26  
fuse 2-15  
FV Threshold 3-23–3-24  
FV Pump Time 3-24

## G

Gas Configuration 3-6, 3-20–3-21, 3-57  
gas exhaust 2-19  
gas flow time 3-23  
gas inlet 2-19  
gas line 2-36  
gas loss 5-12  
gas outlet 4-5  
gas requirements 1-2–1-3  
gas supply 2-12  
gas tank 2-40, 5-10  
Generic Editor 3-38  
glass container 2-2  
gloves 4-26  
goggles 4-26  
gold gasket 2-23, 6-4, 6-7, 6-13–6-14  
gold seal 2-23  
gold stamp 6-6, 6-14, 6-16–6-17  
grain size 5-5  
Grid Errors tab 3-58  
Grid Infos tab 3-58–3-59, 5-11  
grounding 4-6, 4-20–4-21  
guard ring 6-14–6-17  
guide sleeve 2-18  
Gustafsson, O. 4-26

## H

hardware failure 2-32  
Hayes, J.M. 4-27, 5-19  
heat gun 4-6

heat transfer tube 2-23–2-24, 6-12–6-13  
 heatable capillary 6-5  
 heater 2-24, 2-26  
 heater element 2-26  
 heating band 4-6–4-7  
 heating cabinet 2-2–2-3, 3-26, 3-48, 4-19, 6-3  
 heating cartridge 2-25, 4-6, 6-12, 6-26  
 heating wire 2-14  
 Heatout Temperature 3-46  
 high vacuum side 2-18, 6-19–6-20  
 high voltage 3-7, 3-59  
 highest amplified cup 4-4  
 host connection 2-5, 6-29–6-30  
 Host Connection LED 2-12, 6-29  
 HV pump time 3-23–3-24  
 HyS.1 2-28  
 HyS.2 2-28

**I**

IAEA primary standard 6-17  
 idle time 3-10–3-11, 3-24, 4-21  
 Info Properties window 3-18  
 Info window 3-16–3-17, 3-19, 3-25, 5-19–5-21  
 Information bar 3-5  
 Inlet Control board 2-8–2-9, 6-31  
 Inlet Controller External Options 6-27  
 Inlet Controller Internal 6-27  
 inlet port 3-23  
 inlet pressure 5-10  
 inlet system 2-40, 4-12, 5-10  
 inlet valve 4-3–4-4  
 Installation Kit 6-7  
 Instrument Control 2-12, 2-26, 2-30, 4-4–4-6  
 Instrument tab 3-21–3-23  
 integration time 3-10–3-12, 3-22, 3-25–3-26, 3-57, 3-60  
 intensity ratio 4-3  
 interface standby 3-38  
 interfering masses 3-64–3-65  
 interfering masses function 3-64, 3-66  
 Interfering Masses tab 3-64  
 International Atomic Energy Agency (IAEA) 1-3, 6-18  
 international standard 5-14, 5-19  
 ion correction 3-28, 3-57, 5-17–5-18  
 ion correction location 3-28  
 ion intensity 3-24  
 ion source 3-7, 3-10, 3-59, 4-3–4-4, 6-22, 6-26  
 ion source pressure 3-7  
 IRMS 1-2–1-3, 2-12–2-13, 2-20, 2-23, 3-2–3-3, 3-15, 3-44–3-45, 3-59–3-60, 4-3–4-4, 4-6, 4-9–4-10, 4-12, 4-21–4-22, 5-2, 5-6, 5-10, 5-12, 5-18, 5-22, 5-25, 6-22  
 IRMS connector 2-12  
 IRMS distributor 1-2, 2-12

IRMS method 3-33  
 ISL script 3-38, 4-24  
 ISL script language 3-38  
 ISL tab 3-15  
 ISODAT NT Operating Manual 3-4  
 ISODAT NT Operating Manual - Upgrade to Version 2.0 3-4  
 isotope fractionation 2-22, 3-46, 5-10, 6-26  
 Isotope MS 3-22  
 isotope ratio 5-17  
 isotope ratio determination 3-20, 5-3  
 isotopic composition 3-57, 5-17  
 isotopic standard 3-40, 5-18  
 isotopic value 3-29, 5-3

**J**

jacket ring 6-3, 6-6  
 Jumo itron 16 temperature controller 2-6, 2-27

**K**

Kiel IV Carbonate window 3-7–3-8, 3-44, 4-22–4-23  
 Kimwipes® 4-22  
 knife edge 2-19

**L**

lab boy 2-4, 2-26  
 laminar flow 5-24  
 largest signal achievable 3-45, 5-12  
 layout 2-2, 2-18, 5-3  
 leak rate 3-47, 3-59, 4-14, 4-23, 5-21  
 leak test 2-13, 3-26, 3-46, 3-59, 4-2, 4-4–4-5, 5-8, 5-21, 6-19–6-22  
 Leak Test Rise Time 3-47  
 leak-tight 2-13, 2-16, 2-39, 4-8, 6-19  
 left side capillary 3-10  
 level indicator 2-36  
 level sensor 2-38  
 LIMS 3-14  
 line 1 2-4, 2-29, 2-31, 3-26, 3-33, 3-52, 5-8, 5-20–5-21, 6-18–6-19  
 line 2 2-29, 2-31, 3-26, 3-33, 3-52, 3-54, 5-8, 5-21, 6-19  
 linearity error 5-14  
 linearity region 5-13–5-14  
 LineX Sample Position 6-32  
 lint-free cloth 4-22  
 lint-free wiper 4-19  
 liquid fill and decant valve 2-35–2-36, 2-38  
 liquid nitrogen 2-4, 2-23–2-24, 2-26, 2-36, 2-38–2-39, 3-58, 4-20, 4-26–4-27, 5-6, 5-9, 5-21, 6-26  
 liquid nitrogen distributor 2-36–2-38, 6-7  
 liquid nitrogen filling tube 5-6

liquid nitrogen flow 2-24, 2-36  
 liquid nitrogen level 2-38  
 liquid nitrogen refill 2-38–2-39, 4-20, 4-23  
 liquid nitrogen refill device 2-34–2-35, 5-6  
 liquid nitrogen refill valve 4-20  
 liquid nitrogen safety unit 2-35  
 liquid nitrogen storage 2-36  
 liquid nitrogen supply 2-36  
 liquid nitrogen tank 1-3, 2-36, 4-20, 5-6  
 liquid nitrogen temperature 2-39  
 liquid nitrogen transport 6-26  
 LN2 level sensor 6-29–6-30  
 LN2 overpressure gauge 2-35  
 LN2 Refill Unit Enable (LED) 6-29–6-30  
 Load Magazine 2-14, 2-30, 3-49, 3-51  
 log file 3-17–3-19  
 Lotus format 3-37  
 Luer-type syringe 2-13

## M

m/z 16 4-3  
 m/z 17 4-3  
 m/z 18 4-3  
 m/z 28 4-3, 6-21–6-22  
 m/z 32 4-3  
 m/z 40 4-3–4-4, 4-6  
 m/z 44 3-24–3-25, 3-45, 3-59, 3-62–3-63, 4-3, 5-12  
 m/z 45 3-24, 3-59, 3-62–3-63, 5-18  
 m/z 46 3-24, 3-59, 3-62–3-63, 4-4  
 magazine 2-14, 2-29–2-30, 3-50–3-51, 3-53–3-54, 4-19, 4-23, 5-6–5-8, 5-21–5-22  
 magnet current 3-7, 4-3  
 magnet step 4-3  
 magnet valve 5-6  
 main switch 2-6, 2-15  
 mains cable 2-12  
 manifold block 2-20  
 manometer 2-42  
 manual cleaning 4-28–4-29  
 manual cutoff valve 2-36–2-37  
 manual valve 2-3, 2-16, 2-40–2-41, 2-43, 4-7, 4-20, 5-6  
 marine microfossils 5-3  
 mass scan 4-3–4-4  
 mass spectrometer background 3-24  
 master 3-25  
 matching sample capillary to standard capillary 2-16  
 maximum acceptable ion intensity difference 3-25  
 maximum intensity 4-3  
 maximum ion intensity achievable 3-45  
 maximum number of drops 3-48  
 maximum pressure rise accepted 3-26  
 measure only selection 3-37

measurement temperature 3-46  
 membrane assembly 6-15–6-17  
 Method icon 3-21  
 method, new 3-14, 3-20–3-21, 3-29, 4-25  
 Methods tab 3-13  
 methylchlorine 4-29  
 micro balance 1-3  
 Microvolume 2-13, 2-17, 2-23–2-25, 3-59, 4-6, 4-9, 4-11, 5-9–5-10, 6-7, 6-12–6-13, 6-22  
 minimum ion intensity acceptable 3-45  
 minimum pressure 3-23  
 molecular  $\delta$  value 3-57  
 monoblock 2-20  
 monoflop 2-8  
 morphotype 5-19  
 Motor control direction 6-28–6-29  
 Motor control enable 6-28–6-29  
 Motor torque setting 6-28, 6-30  
 MoveToPosition 6-32

## N

National Bureau of Standards (NBS) 1-3  
 NBS-19 1-3, 5-5, 5-11  
 needle valve 6-22  
 neogloboquadrina pachyderma 5-19  
 New button 3-20, 3-31  
 NIST 5-14  
 nominal closing pressure 2-18  
 non-condensable gases 2-22, 3-27, 5-3, 5-21  
 Nr. of Exp. 3-59  
 number of acid drops 3-27, 3-48  
 number of cycles 3-11–3-12, 3-24  
 number of repetitions 3-12, 3-25, 3-57–3-58  
 number of samples 3-32

## O

O ring seal 2-18, 2-35–2-36, 4-6, 5-6, 6-16  
 Object Properties 3-7  
 offline view 3-19  
 OFFS 2-27  
 oil level 6-26  
 oil reservoir 6-26  
 oil trap 1-3  
 optical fiber 2-12  
 Options button 3-16  
 orthophosphoric acid 4-26, 6-7  
 outlet pressure range 1-2  
 outlier 3-27, 3-57  
 outlier test 3-26–3-27, 3-57  
 output grid 3-12  
 oven 2-3, 2-11, 2-13–2-14, 2-27, 2-29–2-30, 2-34, 3-50–3-51,

3-59, 4-18, 4-23, 5-6–5-7, 6-15, 6-17  
 oven control 2-6, 2-27, 4-23  
 oven rack 4-29–4-30  
 oven section 2-8, 2-11  
 oven temperature 2-6, 2-13, 2-27, 3-26, 3-48, 3-50, 5-6–5-7, 5-21

## P

P key 2-27–2-28  
 P VM1 Pre-Reaction 3-59  
 paleoclimatological research 5-5  
 Parafilm 4-28  
 Pb.1 2-28  
 Pb.2 2-28  
 PDB 5-14  
 peak center 3-23, 3-32, 3-59, 5-22  
 peak center cup 3-23  
 peak center cycle 3-23  
 Peak Detection tab 3-29  
 PeriCon 6-27  
 Peripherals tab 3-21, 3-23–3-26  
 personal injury 2-26, 2-34  
 phosphoric acid 2-2, 2-32, 3-48, 4-14, 4-21, 4-26–4-28, 5-3, 5-7, 6-19  
 phosphoric acid preparation 1-3, 2-13, 4-26  
 phosphorous pentoxide 1-3, 4-26–4-28  
 pinch valve 2-3, 2-8, 2-44–2-45, 4-14–4-16, 4-21, 6-3, 6-5, 6-9, 6-28  
 Pinch valve Line 1 6-28–6-29  
 Pinch valve Line 2 6-28–6-29  
 Pirani gauge 3-9  
 piston assembly 6-15–6-17  
 piston height 2-32  
 Piston Line 1 down 6-28–6-29  
 Piston Line 1 up 6-28–6-29  
 Piston Line 2 down 6-28–6-29  
 Piston Line 2 up 6-28–6-29  
 Piston Sensor Line 1 6-28–6-29  
 Piston Sensor Line 2 6-28–6-29  
 piston speed 2-31  
 pliers 2-18, 4-21, 6-15  
 plunger 2-18–2-19  
 pneumatic lever 2-32  
 pneumatic valve 2-18, 6-19, 6-22, 6-24–6-25  
 position sensor 2-3, 2-31, 2-33–2-34  
 position sensor array 2-33, 6-28, 6-30–6-31  
 postdelay 3-23  
 power consumption 1-2  
 power corruption 1-2  
 power disruption 1-2  
 Power Distribution board 2-8, 2-17, 6-27  
 Power Distributor MS 6-27  
 power failure 2-21

power indicator 2-5  
 power requirements 1-2  
 Pre process drop interval 3-48  
 preamplifier dynamic range 3-45  
 precision 1-3, 3-11–3-12, 4-6, 5-4, 5-6, 5-10, 5-14  
 predelay 3-10, 3-23–3-25  
 press adjust 3-25, 3-33, 3-59, 5-11–5-13, 5-24  
 pressure control valve 2-37  
 Pressure Fast Adjust tab 3-45  
 pressure gauge 2-35–2-36, 2-39  
 pressure limiter 2-36  
 pressure meter VM2 3-51  
 pressure raising valve 2-36  
 pressure regulator 2-36, 2-42  
 pressure threshold 3-45, 4-8  
 pressure transducer 3-9, 5-9  
 pressurizing screw 2-23, 6-12–6-13  
 primary ratio 3-39  
 primary standard 3-38–3-39, 3-42, 5-14–5-15, 6-17–6-18  
 primary standard ratio 3-41–3-42  
 primary standards database 3-39  
 Printout tab 3-21, 3-29  
 printout template 3-29  
 Process Timing tab 3-47  
 protective clothing 2-26, 2-34  
 protective gloves 2-26, 2-34  
 proximity switch 2-32–2-33, 3-51, 4-8, 4-20, 4-25–4-26, 6-3, 6-6, 6-21, 6-28–6-29  
 PTFE tubing 6-26  
 pump controller 2-10, 2-13  
 pump controller status 2-17  
 pump mode 3-48  
 Pump Overlay Time 3-23, 5-24  
 pump temperature 3-47  
 pump time 3-23, 3-27, 3-47  
 Pyrex tubing 4-26

## Q

Quick Access buttons 3-7  
 quick release connection 2-12

## R

r.t. 2-28  
 Ratio Editor 3-57  
 Ratio Group 3-57  
 Ratios tab 3-63–3-64  
 Raw Complete tab 3-61  
 raw ratio 5-16  
 Raw Ratios Complete tab 3-63  
 Raw Ratios Reference tab 3-62–3-63

## Index: S

- Raw Ratios Sample tab 3-62
- Raw Reference tab 3-60–3-61
- Raw Sample tab 3-60
- Raw tab 3-56–3-57
- reaction time 1 3-27
- reaction time 2 3-27
- Ref. Name 3-29
- reference bellows 5-11, 5-22, 5-24
- reference gas 1-3, 2-39–2-40, 3-23–3-24, 3-29, 3-57, 3-61, 4-3, 5-12, 5-14, 5-16, 5-22, 5-24
- reference gas amount 3-33, 5-24
- reference gas capillary 3-61
- reference gas pressure 5-22
- reference gas refill 2-39–2-40, 5-22
- reference gas refill unit 2-40
- reference gas supply 2-40
- reference gas tank 2-40, 4-3
- reference matching 4-12
- Reference Refill 2-40, 3-3, 3-21, 3-23, 3-33, 3-45, 5-13, 5-22–5-24
- Reference Refill algorithm 5-24
- Reference Refill device 2-41–2-42
- Reference Refill parameter 2-41, 5-13
- Reference Refill tank 3-23
- Reference Refill tank capillary 2-12
- reference sample 1-3, 5-13
- reference side 3-60, 4-12
- reference tank 4-7
- refill device 2-34–2-36, 2-38
- refill sensor 2-38–2-39, 4-20
- Refill Time 3-23, 3-45, 5-24
- refill valve 2-38
- registry use 6-32
- release temperature 2-22
- reservoir pressure 2-36
- resistor 2-24, 3-12, 3-26, 4-20, 6-12
- response curve 5-11, 5-25
- result file 3-57, 3-64, 5-21
- result path 3-14–3-15
- Result Workshop 3-14, 3-16, 3-29, 3-31, 3-35
- results export 3-35–3-36
- results storage 3-36
- Results tab 3-14, 3-36
- reverse ion correction 3-57
- right side capillary 3-10
- rotary pump 2-16, 5-20–5-21
- safety valve 2-37–2-38
- sample amount 3-33, 5-3
- sample flow 4-9
- sample inlet 4-7
- sample ion intensity 3-25
- sample measurement 3-58, 3-63, 3-65, 5-8, 5-19, 5-24
- sample port 4-10
- sample preparation 3-33, 3-59, 5-2, 5-6, 5-8–5-9, 5-19
- sample pressure 5-22
- sample side 3-60, 4-11, 6-21
- sample signal 3-59, 5-11–5-12, 5-24
- sample size 3-59, 5-12, 5-25
- sample vial 3-26, 3-59, 5-5, 6-3, 6-6
- sample vial vent gas 1-3, 2-10
- sample weight 3-33, 5-11
- Santrock, J. 5-19
- Scans tab 3-15
- Schmitz, B. 4-26
- script editing 3-22
- Search tab 3-15
- Segl, M. 4-26
- sensitivity 4-3, 5-24, 6-32
- sensor status 2-33
- sequence grid 3-32–3-33, 3-59, 5-8
- Sequence icon 3-31
- sequence line 3-60
- Sequence Line tab 3-60
- sequence, new 3-14, 3-29, 3-31, 3-33, 4-25
- Sequences tab 3-13
- serial bus interface 2-8–2-9, 2-12
- service script 3-38, 3-49, 3-51–3-53, 3-55
- Set Path 3-14
- setpoint 2-27–2-28
- sharing violation 3-18
- short circuit 4-20
- shortcut link tube 4-10
- shot noise 3-11–3-12
- signal 3-25, 3-45, 3-57, 5-11–5-13, 5-24–5-25
- signal decrease 3-61
- signal height 3-12, 4-9, 4-12, 5-13
- signal settling time 4-12
- signal up 3-25
- silicone tube 2-13
- single measurement 3-13, 3-57
- site requirements 1-2
- sleeve 2-18–2-19
- slope threshold 3-26
- smallest useful signal to set 3-45
- SNOOP 4-4
- software version 3-14, 3-20, 3-29, 3-31, 3-34
- solenoid valve 2-36–2-37
- Solnhofen limestone 5-15
- Source Control Gnd 6-27

## S

- safety 2-35
- safety device 2-38
- safety nut 2-32
- safety pressure relieve valve 2-36
- safety switch-off 2-8

Source Control HV 6-27  
 SP.H 2-27  
 SP.L 2-27  
 spare parts 2-1, 6-3  
 spring 6-15–6-17  
 spring assembly 2-18, 2-32  
 spring plate 4-16–4-17, 4-26, 6-14, 6-21  
 stainless steel capillary 2-13, 2-23  
 stainless steel frit 4-17  
 stainless steel membrane 2-19  
 stand 2-29  
 standard bellow 3-23, 3-25  
 standard database 3-29  
 standard deviation 3-11–3-13, 3-27, 3-57–3-58  
 Standard Editor 3-29, 3-38, 3-40–3-42  
 standard error 3-12–3-13, 3-57–3-58  
 standard flow 4-9  
 standard gas 3-25  
 standard gas refill 2-39  
 standard inlet 4-7  
 standard parameter 3-28–3-29  
 standard sample 5-17  
 standard side 2-12, 2-40, 3-24, 4-9, 5-22  
 standard soap solution 4-4  
 Standby and Drop 3-48–3-49, 3-53  
 Standby and Pump 3-49, 3-55  
 Standby Drop Interval 3-48, 3-54  
 Standby mode 2-13, 3-55  
 Standby Temperature 3-46  
 StandbyAndDrop.isl 3-54  
 Start button 3-35  
 Start Transfer Temperature 3-46  
 statistical noise 3-12  
 Status bar 3-5–3-6, 3-22  
 stirrer 1-3, 4-26–4-28  
 Studley, S.A. 5-19  
 suffocation 2-35  
 support rod 2-30  
 Swagelok® 6-12  
 Swagelok® connection 4-16–4-17, 6-19  
 Swagelok® connector 1-2, 2-19–2-20, 4-7  
 Swagelok socket 2-23  
 synthetic wiper 4-19

## T

T piece 2-11  
 Take Magazine 2-30, 3-49–3-50  
 tank head pressure 4-20  
 Teflon®gasket 4-16, 6-4  
 Teflon®tube 2-16  
 Teflon washer 4-19, 6-19

Temp Settings tab 3-46  
 temperature controller 2-6, 2-27–2-28, 3-26, 4-23  
 temperature error 6-26  
 temperature sensor 2-24–2-27  
 temperature set point 3-26  
 temperature stability 5-10  
 Terminate Script 3-54  
 thermal contact 2-26  
 Thermo Electron 2-36–2-37, 3-14, 3-20, 3-29, 3-31  
 Time allowed for pumping trap(s) 3-47  
 time constant 2-28  
 time events 3-22  
 time slicing 3-26, 3-60  
 Time Slicing mode 3-60  
 TMH 071 P 2-16  
 tolerance 3-25, 4-11  
 Total CO<sub>2</sub> 3-59  
 total volume 2-23  
 trace element content 5-4, 5-9  
 transfer time 3-27  
 transfer tubing 6-24–6-25  
 Trap 1 cooling resistor 6-29–6-30  
 Trap 1 Heater Disable 6-28–6-29  
 Trap 1 temperature set value 6-28, 6-30  
 Trap 1 temperature true value 6-28, 6-30  
 Trap 2 cooling resistor 6-29–6-30  
 Trap 2 Heater Disable 6-28–6-29  
 Trap 2 temperature set value 6-28, 6-30  
 Trap 2 temperature true value 6-28, 6-30  
 trap arrangement 2-17  
 trap region 5-20  
 trap temperature 2-22, 3-46, 4-23  
 trap 1 2-3–2-4, 2-10, 2-13, 2-16–2-17, 2-22–2-23, 2-26, 3-27, 3-46, 3-59, 4-6, 4-10, 4-14, 4-23, 4-26, 5-7, 5-9, 5-11–5-12, 5-21–5-22  
 trap 2 2-10, 2-13, 2-16–2-17, 2-21–2-23, 2-26, 3-27, 3-46, 3-59, 4-6, 4-10–4-11, 4-14, 4-23, 4-26, 5-7, 5-9, 5-21–5-22  
 trapping efficiency 2-22  
 trapping section 2-16, 5-10, 6-26  
 trapping volume 2-21–2-23, 6-9–6-10, 6-12–6-13  
 trapping volume valve block 2-21  
 TubeCracker 2-21  
 Tune Scan 2-28, 4-4, 6-21  
 turbo pump 2-7, 2-10, 2-16–2-17, 2-42, 3-23, 4-9, 5-23  
 Turbo pump 50% 6-29–6-30  
 turbo pump controller 2-17  
 turbo pump error 6-29–6-30  
 turbo pump status 2-5  
 turret 2-2–2-3, 2-29–2-30, 3-33, 3-50, 4-23, 5-5, 6-28  
 turret control 4-23, 5-19  
 turret motor 2-3, 2-33–2-34  
 turret position 2-34, 6-7, 6-28, 6-30

## U

umbilical view 5-19

uninterruptible power supply (UPS) 1-2

## V

vacuum connection 2-23

vacuum gauge VM2 3-26

vacuum line 2-42

vacuum scheme 4-13, 6-2, 6-35

vacuum system 2-2, 2-6, 2-13, 2-16, 4-6, 5-9

valve 1 6-28, 6-30

valve 2 6-28, 6-30

valve 3 6-28, 6-30

valve 4 6-28, 6-30

valve 5 6-28, 6-30

valve 7 6-28–6-29

valve 8 6-28, 6-30

valve block 2-19–2-20, 2-23, 4-6, 4-17, 6-12–6-13, 6-16

valve block assembly 6-17

valve flange 4-18

valve plug 2-19

valve system 2-10

valve unit 6-12–6-14, 6-19

vent gas 2-6, 4-8

vent gas connection 2-6

vent gas pressure 4-8

vent valve 2-10, 2-17, 2-35

vial connection 4-21, 4-23

vial disconnection 4-21, 4-23

vial section 2-5, 4-8

vial test 2-16, 4-23–4-25

vial test method 4-24–4-25

vial vent gas 2-11

vial 1/1 3-51, 5-8, 5-21

vial 2/1 3-51, 3-54, 5-20–5-22

vialtest.isl 4-24

viscous flow 2-23

Viton®ring 4-18–4-19

Viton tubing 2-3, 2-13–2-14, 2-16, 2-44, 4-21

VM 1 pressure true value 6-29–6-30

VM 2 pressure true value 6-29–6-30

VM1 connection 6-25

VM1 Leak Threshold 4-25

VM1 expansion 3-27, 3-59, 5-25

VM1 pressure 3-59, 5-11, 6-18–6-19, 6-22, 6-24–6-25

VM2 Leak Threshold 3-26

VPDB 5-14–5-15, 5-17

## W

Wachter, E.A. 4-27

waiting time 3-23–3-25, 3-47

washer 2-23, 2-25, 6-12–6-13

water level 4-4

WEEE compliance 1-1

weighing 5-4–5-5

weighing instruments 1-3

Winchester 4-26

working standard 3-40–3-41, 5-14, 5-17, 5-19

## Y

Y.0 2-27–2-28

Y.1 2-28

Y.2 2-27–2-28

Thermo Electron Corporation  
81 Wyman Street  
P.O. Box 9046  
Waltham, Massachusetts 02454-9046  
United States

[www.thermo.com](http://www.thermo.com)

Analyze • Detect • Measure • Control™

**Thermo**  
ELECTRON CORPORATION